

Leading Gastric Pressure Relief Device

The FARRELL® Valve System is a closed reservoir overflow system which is designed to help patients who suffer from poor gastric motility, pain, and bloating. The FARRELL Valve maximizes enteral feeding by providing a channel to constantly decompress the stomach, allowing the stomach to fill at its own pace. By providing a reservoir channel, the clinician can maximize enteral feeding and minimize the need for TPN.

ENFit™ connectors are designed to help prevent feeding tube misconnections and enhance safety.*

The FARRELL Valve System seamlessly transitions patients across the continuum of care.

BENEFITS OF REDUCING GASTRIC DISTENSION AND EXCESS FLUID:

- Reduces risk of pulmonary aspiration¹
- Relieves pressure on recent gastric surgeries (e.g. fundoplication^{2,3})
- Improves patient comfort and reduces pain^{5,6,7}
- Facilitates enteral feeding tolerance and helps patients reach caloric goals¹
- Patients able to tolerate enteral feeding may not need alternative nutrition such as TPN⁸

THE FARRELL VALVE SYSTEM:

Closed Overflow Reservoir

- First closed system for gastric fluid containment and simultaneous venting
- Contains and retains fluid until patient can tolerate fluid volume
- Prevents loss of medications and formula
- Prevents caregiver exposure to gastric contents

Equilibrium Technology

- Equalizes pressure between patient's gastric system and the external environment
- Reduction of gas and fluid containment may allow patients to achieve caloric goals


Precise Graduation Markings

- Increased resolution on bag volume readings
- Improved monitoring of reflux levels



Self Centering Hanger Tab

- Improved ambulatory use



Adjustable Roller Clamp

- Improved fluid flow and control, upon refeeding



Bi-directional Flow Indicators

- Inform users on fluid path




"My son had retching, gagging and bloating post Fundoplication. We were given medication to help with this. It didn't totally eliminate his symptoms and our community pediatric nurse suggested "FARRELL Valve bags." They are amazing and help him lots with his symptoms making him more comfortable as it decompresses his stomach. I would highly recommend these to anyone who is enteral fed. We just hook the Farrell valve up at the same time as his feed and it vents him at the same time. When he is particularly symptomatic, we attach the FARRELL directly to his g-tube. Who would think a plastic bag would be an amazing, can't-do-without piece of equipment!"

Travis' Mum

1. Kazi N et al. Enteral feeding associated gastroesophageal reflux and aspiration pneumonia: a review. Nutrition Review. 1996;54(10):324-328. 2. Dunn et al. Long-term quantitative results following fundoplication and antroplasty for gastroesophageal reflux and delayed gastric emptying in children. Am J Surg. 1998;175:27-29. 3. Fonkalsrud E et al. Surgical treatment of gastroesophageal reflux in children: a combined hospital study of 7467 patients. Pediatrics. 1998;101(3):419-422. 4. Verkalon T, Walden M, eds. Core Curriculum for Neonatal Intensive Care Nursing, 4th ed. National Association of Neonatal Nurses. 5. Salet G et al. Responses to gastric distension in functional dyspepsia. Gut. 1998; 42(6):823-829. 6. Ladabaum U et al. Gastric distension correlates with activation of multiple cortical and subcortical regions. Gastroenterology. Feb 2001; 2:369-376. 7. Mollitt DL. Symptomatic gastroduodenal pancreatic rest in children. J Ped Surg. Aug 1984; 19(4):449-450. 8. Guidelines for the use of parenteral and enteral nutrition in adult and pediatric patients. JPEN 2002 26: 15A 9. Reese Parish C. Enteral feeding: The art and science. Nutrition in Clinical Practice 18:76-85, Feb 2003. NOTE: Citations are for reference purposes only. These studies were not about the FARRELL Valve System, however, they do support the need for relieving gastric pressure in certain patient groups.

*ISO 80369-compliant



“Intolerance of enteral feeding may result in extended hospital stays, increased costs, and potential malnutrition.”

Gerry Hennies, RN
Cincinnati Children’s Hospital

Use the FARRELL Valve System to relieve gastric distension that may develop due to the following medical therapies or conditions:

- Enteral Feeding
- Funduplications
- Gastroesophageal Reflux Disease (GERD)
- Intestinal Motility Disorder
- Delayed Gastric Emptying (DGE)
- Obstructive Sleep APNEA/CPAP
- Cystic Fibrosis
- Muscular Dystrophy
- Intestinal Pseudo-Obstructions
- Swallowing Disorders
- Aspirate Pneumonia
- Cyclic Vomiting Syndrome
- Ventilated Patients
- Surgical Repair of Intestinal Atresia/Obstruction
- Acute Respiratory Distress Syndrome (ARDS)
- Neurologically Impaired
- Hirschsprung Disease
- Unable to Vomit
- Comatose
- Congestive Heart Failure
- Intubated Patients
- Transitioning TPN to Enteral Feeding

PRODUCT CODE	DESCRIPTION	UNIT OF SUPPLY
44-4100	FARRELL Valve System with ENFit [™] Connectors [†]	30/case

[†]Product not made with: Natural rubber latex
Phthalates (DEHP)
Bisphenol A (BPA)



TO ORDER

IN THE UK
0800 144 4480
Fax: 0800 151 3191
orders@corpakuk.com
CORPAK MedSystems UK
The Beehive, City Place
Gatwick RH6 0PA

NORTHERN IRELAND HOSPITAL ACCOUNTS:
028 9056 7135
Fax: 028 9070 5623
patricia.cartlidge@sangers.co.uk

NORTHERN IRELAND OTHER ACCOUNTS:
028 9040 3409
Fax: 028 9040 1958
bradburycustomerservices@sangers.co.uk

THE REPUBLIC OF IRELAND
01 468 8456
Fax: 01 466 0782
allpharservices@uniphar.ie
Allphar Services Limited
Orders Department
4045 Kingswood Road Citywest
Business Park
Co. Dublin

For more product information, please visit
www.corpak.com

CORPAK University Training
& Education Portal:
www.corpakU.com



FARRELL, FARRELL logo, EQUILIBRIUM TECHNOLOGY, CORPAK, CORPAK logo, CORPAK University, CORPAK University logo, and the color YELLOW for feeding tubes are either registered trademarks or trademarks of CORPAK MedSystems, Inc. in the United States and/or other countries. ENFit[™] is a trademark of GEDSA.