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Referenser

1. Waring MJ and Parsons D. Physico-chemical characterisation of carboxymethylated spun cellulose fibres. *Biomaterials*. 2000;22:903-912.
2. Walker M, Hobot JA, Newman GR, Bowler PG. Scanning electron microscopic examination of bacterial immobilisation in a carboxymethylcellulose (Aquacel) and alginate dressings. *Biomaterials*. 2003;24(5):883-890.
3. World Union of Wound Healing Societies. Principles of best practice: Wound exudate and the role of dressings. A consensus document. London: MEP Ltd; 2007.
4. Bowler PG, Jones SA, Davies BJ, Coyle E. Infection control properties of some wound dressings. *J Wound Care*. 1999;8(10):499-502.
5. Lawrence JC, Lilly HA, Kidson A. Wound dressings and airborne dispersal of bacteria. *Lancet*. 1992;339(8796):807.
6. Coutts P, Sibbald RG. The effect of a silver-containing Hydrofiber dressing on superficial wound bed and bacterial balance of chronic wounds. *Int Wound J*. 2005;2(4):348-356.
7. Robinson BJ. The use of a hydrofibre dressing in wound management. *J Wound Care*. 2000;9(1):32-34.
8. Hoekstra MJ, Hermans MHE, Richters CD, Dutrieux RP. A histological comparison of acute inflammatory responses with a Hydrofiber or tulle gauze dressing. *J Wound Care*. 2002;11(2):113-119.
9. Reddy M, Kohr R, Queen D, Keast D, Sibbald G. Practical treatment of wound pain and trauma: a patient-centered approach. An overview. *Ostomy Wound Manage*. 2003;49(suppl 4A):S2-S15.
10. Parsons D, Bowler P, Myles V, Jones S. Silver antimicrobial dressings in wound management: A comparison of antibacterial, physical and chemical characteristics. *Wounds*. 2005;17: 222-232.
11. Performance analysis of AQUACEL vs SORBSAN and SORBSAN PLUS DRESSINGS. WHRI 2725 AI128. May 5, 2005. Data on file, ConvaTec.
12. Armstrong SH, Brown DA, Hill E, Ruckley CV. A randomized trial of a new Hydrofiber dressing, Aquacel, and an alginate in the treatment of exuding leg ulcers. Presented at: 5th European Conference on Advances in Wound Management; November 21-24, 1995; Harrogate, UK.
13. Barnea Y, Amir A, Leshem D, Zaretski A, Weiss J, Shafir R, Gur E. Clinical comparative study of Aquacel and paraffin gauze dressing for split-skin donor site treatment. *Ann Plast Surg*. 2004;53(2):132-136.
14. Caruso DM, Foster KN, Blome-Eberwein SA, et al. Randomized clinical study of Hydrofiber dressing with silver or silver sulfadiazine in the management of partial-thickness burns. *J Burn Care Res*. 2006;27(3):298-309.
15. Kogan L, Moldavsky M, Szvalb S, Govrin-Yehudain J. Comparative study of Aquacel and Silverol treatment in burns. *Ann Burns Fire Disasters*. 2004;17(4):201-207.
16. Harding KG, Price P, Robinson B, Thomas S, Hofman D. Cost and dressing evaluation of Hydrofiber and alginate dressings in the management of community-based patients with chronic leg ulceration. *Wounds*. 2001;13(6): 229-236.
17. Collado GF, Morán SMJ, García RMJ. Care for cutaneous lesions by combining hydrofiber and extra fine hydrocolloid dressings. A multi-centric study. *Rev ROL Enf*. 2002;25(2):130.
18. Foster L, Moore P, Clarke S. A comparison of hydrofibre and alginate dressings on open acute surgical wounds. *J Wound Care*. 2000; 9(9): 442-445.
19. Ravenscroft MJ, Harker J, Buch KA. A prospective, randomised, controlled trial comparing wound dressings used in hip and knee surgery AQUACEL® and Tegaderm versus Cutiplast. *Ann R Coll Surg Engl*. 2006;88:18-22.
20. Jones SA, Bowler PG, Walker M, Parsons D. Controlling wound bioburden with a novel silver-containing hydrofiber® dressing. *Wound Repair Regen*. 2004;12:288-294.
21. Bowler PG, Cochrane CA. Progression toward healing: Wound infection and the role of an advanced silver-containing Hydrofiber® dressing. *Ostomy Wound Manage*. 2003;49(8)(suppl):S 2-S5.
22. Bowler PG, Jones SA, Walker M, Parsons D. Microbicidal properties of a silver-containing hydrofiber dressing against a variety of burn wound pathogens. *J Burn Care Rehabil*. 2004;25:92-196.
23. Walker M, Bowler PG, Cochrane CA. In vitro studies to show sequestration of matrix metalloproteinases by silver-containing wound care products. *Ostomy Wound Manage*. 2007;53(9):18-25.
24. Jurczak F, Dugre T, Johnstone A, Offori T, Vujovic Z, Hollander D et al. Randomised clinical trial of Hydrofiber dressing with silver versus povidoneiodine gauze in the management of open surgical and traumatic wounds. *Int. Wound J*. 2007;4:66-76.
25. Jude EB, Apelqvist J, Sprault M, Martini J. Prospective randomised controlled study of Hydrofiber® dressing containing ionic silver or calcium alginate dressings in non-ischaemic diabetic foot ulcers. *Diabetic Medicine*. 2007; 24:280-288.
26. Vanscheidt W, Lazareth I, Routkovsky-Norval. Safety evaluation of a new ionic silver dressing in the management of chronic ulcers. *Wounds*. 2003;15:371-378.

Kontaktinfo

ConvaTec Sverige AB,

Box 15138,

167 15 Bromma,

Kundservice:

020-21 22 22

Fax: 08-704 89 98,

convatec.kundservice@convatec.com

AQUACEL Ag
Hydrofiber® Dressing with Silver

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© ConvaTec Inc 2010
®/™ AQUACEL, Hydrofiber och DuoDERM är varumärken tillhörande ConvaTec Inc.