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MRI Information

MR Conditional

The 3M™ Red Dot™ Cardiac Monitoring Electrodes (2244, 2268, 2269T, 2570, 2660 or 2670) or 3M™ Red Dot™ Resting Electrodes (2330 or 2360) were determined to be MR conditional. Non-clinical testing demonstrated that the 3M Red Dot cardiac monitoring electrodes (2244, 2268, 2269T, 2570, 2660 or 2670) or 3M Red Dot resting electrodes (2330 or 2360) are MR conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720 or 4000 (See Table Below)-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, the 3M Red Dot cardiac monitoring electrode or 3M Red Dot resting electrode produced the following temperature rise during MRI performed for 15-min of scanning (i.e., per pulse sequence) in the 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR system: *Highest temperature change of +1.7°C was measured during the testing of the 3M™ Red Dot™ Neonatal, Pre-Wired, Radiolucent Monitoring Electrode 2269T.*

Therefore, the MRI-related heating experiments for the 3M Red Dot cardiac monitoring electrode or 3M Red Dot resting electrode at 3-Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9 -W/kg or less (**see Table Below**) (i.e., associated with a calorimetry measured whole body averaged value of 2.7-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.7°C.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the 3M Red Dot cardiac monitoring electrode or 3M Red Dot resting electrode. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) for each electrode extends a few millimeters relative to the size and shape of the 3M Red Dot cardiac monitoring electrode or 3M Red Dot resting electrode.

3M Electrode	Electrode Type	Maximum Spacial Gradient Field (Gauss/cm)	SAR (W/Kg)	Max Heating (°C)	Artifact (mm extended)
2244	Cardiac Monitoring	720	2.9	<1.5	2
2268	Cardiac Monitoring	4000	2	1.4	3
2269T	Pre-Wired Cardiac Monitoring	4000	2	1.7	5
2570	Cardiac Monitoring	720	2.9	<1.5	2
2660	Cardiac Monitoring	720	2.9	<1.5	2
2670	Cardiac Monitoring	720	2.9	<1.5	2
2330	Resting	4000	2	1.3	1
2360	Resting	4000	2	1.3	1

Important Note: This 3M™ Red Dot™ ECG Electrode that underwent MRI testing was not connected to an electrocardiograph (ECG) lead and/or an MR conditional monitoring system. Therefore, the safety of using this 3M Red Dot ECG electrode attached to an ECG lead and/or MR conditional monitoring system is unknown. Furthermore, an MR conditional monitoring system typically requires use only with ECG electrodes specifically designed and tested with that specific equipment. It is unknown if this 3M Red Dot ECG electrode is acceptable for use with an MR conditional ECG monitoring system designed to be used to monitor a patient undergoing an MRI procedure.

[Vad är “Alarm fatigue”?
go.3M.com/AlarmFatigue](https://go.3M.com/AlarmFatigue)

[Hur du applicerar och tar bort EKG-elektroder
\(information på engelska\)
go.3M.com/RedDotHowTo](https://go.3M.com/RedDotHowTo)

3M Svenska AB
Medical Solutions Division
Herrjärva torg 4
170 67 Solna
Tel: 08-92 21 00
www.3M.se/msd

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3M Science. Applied to Life, 3M och Red Dot är varumärken inom 3M Company. Lämnas till återvinning. Tryckt i Sverige. © 3M 2020. Alla rättigheter förbehålles. 17465207_SE20200915 / ht000001062