

Product Specifications

Brand : Cenaman®
Product : ENFit Syringe – single patient use

Product description

Cenaman® ENFit syringe – single patient use are syringes used for delivering liquid medication and enteral nutrition orally or via a feeding tube into the gastro-intestinal tract. These syringes can be reused in the same patient up to 14 days or 84 usages (whichever comes first).

Intended use

The Cenaman® ENFit syringe – single patient use is intended for the delivery of liquid medication, enteral feed and hydration and for the flushing of enteral feeding devices.

Classification

Class Is under the Medical Devices Regulations (SI 2002 No. 618, as amended) and
Class Is under Regulation (EU) 2017/745 Medical Device Regulation (MDR).

Packaging & Product Reference

Product reference	Description	Size (*)	Piece(s)	Inner carton (pcs)	Outer carton (pcs)
C063001	CENAMAN MULTI-USE ENFIT SYRINGE 1ML	1 ML	1	10	1200
C063002	CENAMAN MULTI-USE ENFIT SYRINGE 2,5ML	2,5 ML	1	10	1200
C063005	CENAMAN MULTI-USE ENFIT SYRINGE 5ML	5 ML	1	10	1200
C063010	CENAMAN MULTI-USE ENFIT SYRINGE 10ML	10 ML	1	10	800
C063020	CENAMAN MULTI-USE ENFIT SYRINGE 20ML	20 ML	1	10	360
C063020-30	CENAMAN MULTI-USE ENFIT SYRINGE 20ML X30	20 ML	1	30	480
C063061	CENAMAN MULTI-USE ENFIT SYRINGE 60ML ECCENTRIC	60 ML	1	10	240

Product Specifications

Product Reference	Description	Size (*)	Piece(s)	EAN Inner Carton	EAN Outer Carton
C063001	CENAMAN MULTI-USE ENFIT SYRINGE 1ML	1 ML	1	16971564466483	26971564466480
C063002	CENAMAN MULTI-USE ENFIT SYRINGE 2,5ML	2,5 ML	1	16971564466490	26971564466497
C063005	CENAMAN MULTI-USE ENFIT SYRINGE 5ML	5 ML	1	16971564466506	26971564466503
C063010	CENAMAN MULTI-USE ENFIT SYRINGE 10ML	10 ML	1	16971564466513	26971564466510
C063020	CENAMAN MULTI-USE ENFIT SYRINGE 20ML	20 ML	1	16971564466520	26971564466527
C063020-30	CENAMAN MULTI-USE ENFIT SYRINGE 20ML X30	20 ML	1	56971564461721	66971564461728
C063061	CENAMAN MULTI-USE ENFIT SYRINGE 60ML ECCENTRIC	60 ML	1	16971564466537	26971564466534

*Sizes are approximate

Features

- Enfit connector (ISO 80369-3)
- Clear number markings
- High visibility through barrel
- Purple coloured plunger
- Low dose tip in 1ml, 2.5ml and 5ml
- Smooth plunger movement

Precautions

- Discard the syringe if:
 - any part is damaged, distorted or blocked
 - barrel markings are not clear
 - if syringe becomes too difficult to use or once the syringe has been used for 14 days or 84 usages (whichever comes first).
- When connecting or disconnecting the syringe, always insert or remove the male ENFit connector using a “straight-on/straight-off” approach. Angled insertions or removals may cause poor functioning and/or damage to the connectors. If syringe is not correctly connected, leakage may occur. The tip caps are a choking hazard. Keep out of reach of children.
- The device is not for IV use.
- Forceful squirting oral liquid medication to the back of the throat may cause choking.
- Always refer to nutrition/medication manufacturer indications and instructions for use.
- Please note that de manufacturer does not claim compatibility with specific medication as this

Product Specifications

is the responsibility of the medicine manufacturer. The end user should confirm that the medication is suitable for enteral or oral administration via plastic syringes prior to application.

Instructions for reuse

Do not autoclave or sterilize syringe as this may cause syringe to malfunction. Do not place syringe in dishwasher as this may cause syringe to malfunction.

1.

Remove plunger from syringe barrel.
2.

Wash the barrel and plunger separately in warm soapy water until visible dirt is removed. Follow local guidelines on the appropriate soap.
3.

If visible soil remains, repeat cleaning steps using a soft-bristled brush to aid in washing the plunger and the INSIDE of the syringe barrel. Do not use abrasive products on the OUTSIDE of the barrel as this may lead to premature removal of printed markings on the syringe.
4.

Rinse barrel and plunger with cold running tap water until no traces of soap remains.
5.

Shake or/and absorb with gauze any residual water of the syringe components and then allow both parts to air dry in a clean area before putting plunger back in barrel. Repeat cleaning steps if visible soil remains.
6.

Store in a clean dry container with lid.

Specifications

Component	Material
Barrel	Polypropylene
Plunger	Polypropylene
Lubricant	Silicone Oil
Gasket	Rubber

Absence of specific components

Component	Present
Latex	No
Phthalates	No
Substances of animal origin	No
Substances of human origin	No
Blood or blood derivates	No
Medicinal substances	No

Product Specifications

Heavy metals	No
PVC	No

Sterilization

EO sterilized.

Storage conditions

Store the product in a cool, dry environment away from direct sunlight, high heat, and high humidity.

Shelf life

Cenaman® ENFit syringe – single patient use have a shelf-life of 5 years from production date, when stored under specified conditions.

Picture



Total Cenaman syringes range

Approval	Name	Date	Signature
Associate Category Management Director	Susanne van den Berg	Susanne vd Berg	Digitaal ondertekend door Susanne vd Berg Datum: 2026.05.12 16:25:41 +02'00'
RA Specialist	Alexandra Duarte	Alexandra Duarte	Digitally signed by Alexandra Duarte Date: 2026.05.05 14:25:02 +01'00'