

# Suprapubic Catheter

## 100% Silicone, PVC-free

### Open Tip



#### Area of application

Suprapubic catheters are used for percutaneous bladder drainage when the urethra is damaged or blocked, or when an intermittent catheter is not applicable. The products are for single use only.

#### Product range

Total Length, mm (±10)	Balloon Capacity, ml	CH/Fr.	Tip	Packaging	Reference
<b>All Silicone Suprapubic Catheter, Open tip</b>					
280	5	12	Open	1/10/200	41191205
280	5	14	Open	1/10/200	41191405
280	10	16	Open	1/10/200	41191610

#### Sterilization

Sterile, Sterilized by ethylene oxide, Shelf life: 3 years

#### Certifications

<b>CE-marking</b>	CE 0123
<b>Product classification</b>	Class IIA (EU)
<b>Quality system certified to</b>	ISO 13485 Medical Device Single Audit Program (MDSAP)
<b>European Medical Device Nomenclature (EMDN)</b>	U040102 - KITS WITH DRAINAGE CATHETERS AND INTRODUCERS
<b>Global Medical Device Nomenclature (GMDN)</b>	34924

#### Instructions for use

Included?	Yes
-----------	-----

#### Material specifications

<b>Contains latex?</b>	No
<b>Contains phthalates?</b>	No
<b>Contains PVC?</b>	No
<b>Contains bisphenols?</b>	No
<b>Contains cotton?</b>	No
<b>Contains paper?</b>	Packaging only
<b>Catheter shaft and funnel</b>	Silicone
<b>X-ray line</b>	Silicone
<b>Balloon</b>	Silicone
<b>Hard valve</b>	Polypropylene, Polyester, Nitrile, Stainless steel
<b>Valve sleeve</b>	Polypropylene
<b>Catheter Sleeve</b>	Ethylene-vinyl acetate

# Suprapubic Catheter

## 100% Silicone, PVC-free

### Open Tip



#### Packaging information

Reference	Primary pack./Single peelpack		Secondary pack./Shelf box		Tertiary pack./Shipper			Gross weight kg	Pallet Qty. pcs
	Qty. pcs	Dimensions, mm (LxW)	Qty. pcs	Dimensions mm (LxWxD)	Qty. pcs	Dimensions mm (LxWxD)	Volume m <sup>3</sup>		
41191205	1	580 x 80	10	585 x 85 x 85	200	605 x 445 x 365	0.098	5.00	6,400
41191405	1	580 x 80	10	585 x 85 x 85	200	605 x 445 x 365	0.098	5.30	6,400
41191610	1	580 x 80	10	585 x 85 x 85	200	605 x 445 x 365	0.098	5.50	6,400

#### GTIN codes

REF	Primary package	Secondary package	Tertiary package
41191205	06925816626098	16925816626095	26925816626092
41191405	06925816626104	16925816626101	26925816626108
41191610	06925816626111	16925816626118	26925816626115

#### Packaging material










Primary package	Sleeve: Ethylene-vinyl acetate Peel pouch: Coated paper, Film
Secondary package	Single-wall corrugated fiberboard
Tertiary package	Double-wall corrugated fiberboard

# Suprapubic Catheter

## 100% Silicone, PVC-free

### Open Tip




Storage and handling	
	This way up
	Keep dry
	Fragile, handle with care
	Keep away from sunlight
	Do not resterilize
	Do not re-use
	Do not use if package is damaged and consult instructions for use
	Temperature limit (-15° to 49°C)
	Consult instructions for use or consult electronic instructions for use

**DISTRIBUTED BY**  
**UNOQUIP GMBH**  
 Schützengraben 20  
 8200 Schaffhausen  
 Switzerland  
[info@unoquip.com](mailto:info@unoquip.com)  
[www.unoquip.com](http://www.unoquip.com)

**SPONSORED IN AUSTRALIA BY**  
 UNOQUIP Pty Ltd  
 PO Box 1117  
 Windsor VIC 3181  
  
 1300 55 80 79  
[info@unoquip.au](mailto:info@unoquip.au)



 Well Lead Medical Co., Ltd.  
 C-4 Jinhui Industrial Estate  
 Hualong, Panyu Guangzhou  
 511434 P.R. China

**EC REP** Shanghai Intl. Holding  
 Corp. GmbH  
 Eiffestrasse 80  
 20537 Hamburg  
 Germany

**CH REP** QS Engineering AG  
 Rütihof 1  
 CH 8820 Wädenswil  
 Switzerland

**UK Responsible Person (UKRP)** SUNGO Certification  
 Company Limited  
 3rd Floor, 70 Gracechurch  
 Street, London, EC3V 274HR  
 England, United Kingdom