

Qufora® Leg Bags | Multi Chamber | Drainable

Q-Flow Valve or T-Valve | Adjustable or Fixed Tube Length | Single Packed

Discreet. Comfortable. Easy-to-use.

Data Sheet

Medical Device

500 ml
750 ml

Non-sterile

Sterile

The **Qufora® leg bags** are designed to be worn discreetly under clothing and can be attached to either a urinary catheter or a urisheath.

PRODUCT BENEFITS

- ▶ The **multi chamber** system ensures that the Qufora® leg bags are almost invisible and noiseless under clothing, even when full.
- ▶ The Qufora® leg bags have a soft **adjustable or fixed inlet tube**. The adjustable tube makes it possible to personalise the tube's length. Two sterile leg bags also come with a separate connector/cap.
- ▶ A **non-return valve** ensures that urine will not flow back into the tube, reducing the risk of infection.
- ▶ The front has a **scale** for quick estimation of the urine volume in the bag and the back has a skin friendly **soft backing** for increased comfort.
- ▶ The **Q-Flow** valve provides for a **free flow** drainage of urine in a **hygienic** and **safe** manner.
- ▶ The **T-Valve** is easy and safe to operate with one hand.
- ▶ The sterile Qufora® leg bags come with a **needle-free sample port** for easy urine sampling. The sample port is specially designed for easy cleaning.



QUFORA® LEG BAGS – MULTI CHAMBER - DRAINABLE NON-STERILE - Q-FLOW VALVE OR T-VALVE - ADJUSTABLE TUBE LENGTH - SINGLE PACKED

REF	Volume	Tube length ³	White "stopper"	Soft backing	Quantity retail box	Quantity shipper box
11121501 ¹	500 ml	50 cm	Yes	Yes	10	100
11161502 ²	500 ml	50 cm	No	Yes	10	100
11321501 ¹	750 ml	50 cm	Yes	Yes	10	100
11361502 ²	750 ml	50 cm	No	Yes	10	100

All bags with separate tree connector/cap.

¹T-Valve; ²Q-Flow Valve; ³Adjustable tube length

QUFORA® LEG BAGS – MULTI CHAMBER - DRAINABLE STERILE - Q-FLOW VALVE - ADJUSTABLE OR FIXED TUBE LENGTH - SINGLE PACKED

REF	Volume	Tube length	Needle-free sample port	Soft backing	Quantity retail box	Quantity shipper box
12161802*	500 ml	10 cm ¹	Yes	Yes	10	100
12161602	500 ml	25 cm ¹	Yes	Yes	10	100
12161502	500 ml	50 cm ¹	Yes	Yes	10	100
12361802	750 ml	10 cm ¹	Yes	Yes	10	100
12361902	750 ml	25 cm ¹	Yes	Yes	10	100
12361502	750 ml	50 cm ²	Yes	Yes	10	100

* Wide inlet tube

¹Fixed tube length; ²Fixed tube length and adjustable option with separate connector/cap (no needle-free sample port)

Importer

Qufora A/S
Gydevang 28-30
DK-3450 Allerød
Denmark

Tel: +45 88 30 05 92
info@qufora.com
www.qufora.com

Qufora® Leg Bags | Multi Chamber | Drainable

Q-Flow Valve or T-Valve | Adjustable or Fixed Tube Length | Single Packed

Product description

- ▶ The Qufora® leg bags are single-use medical devices. **All materials are 100% latex- and DEHP-free**

MATERIALS

Films	PVC (plasticiser: DOTP/DEHT)
Inlet tube	PVC (plasticiser: DOTP/DEHT)
Connector	PP: Non-sterile products & separate connector for sterile products: 12161502 & 12361502
Connector (sample port)	PVC; Port: TPE (sterile products)
Connector cap	PE
Non-return valve	PVC (plasticiser: DOTP/DEHT)
Soft backing	Non-woven PET
Q-Flow valve	LLDPE/POM; Inlet tube: PVC (plasticiser: DOTP/DEHT)
T-Valve/white "stopper"	PP/PE Inlet tube: PVC (plasticiser: DOTP/DEHT)
Outlet tube	Silicone
Single bag (product)	LDPE (non-sterile products)
Single bag (product)	HDPE & PE/PET (sterile products)

PACKAGING

Retail box	Cardboard
Shipper	Corrugated cardboard
IFU	One (1) per retail box
Shelf-life	Five (5) years
Storage	Dry in temperatures between 0°C and 35°C; protected from direct sunlight; humidity level max 80%; not to be exposed to heavy pressure

STERILIZATION

Gas sterilized with ethylene oxide (EtO)

PRODUCT CLASSIFICATION

Medical Devices conforming to MDR 2017/745 in Class I /Class I sterile (Is)

APPLIED STANDARDS

ISO 13485	Quality Management System
ISO 8669-2	Urine collection bags - Part 2: Requirements and test methods
GMDN	Urinebags: Non-sterile: 58923; Sterile: 58923

Data Sheet

Medical Device

500 ml
750 ml

Non-sterile

Sterile

Connector with needle-free sample-port



Q-Flow valve



T-Valve



Importer

Qufora A/S
Gydevang 28-30
DK-3450 Allerød
Denmark

Tel: +45 88 30 05 92
info@qufora.com
www.qufora.com