

## 63623 Adhesive Aperture Drape, PRIMARY

Impermeable and absorbing drapes with a variety of apertures

### Product details

---

<b>Product group name:</b>	Adhesive aperture drape
<b>Size:</b>	50cm x 60cm, Aperture 6x8cm
<b>Descriptive feature:</b>	Absorbent, Adhesive, Centered ellipse aperture
<b>Sterile:</b>	Sterile
<b>Brand:</b>	BARRIER®

### Product descriptions

---

**Short description:** Impermeable and absorbing drapes with a variety of apertures

**Long description:** BARRIER® adhesive aperture drapes are available in a range of sizes and aperture shapes designed for different surgical needs. You can choose centred or off-centre aperture positions. We also offer variable aperture drapes which can be adjusted to different needs.


- Variety of adhesive aperture shapes and sizes to suit different procedures.
- Absorbent top layer
- Impermeable materials help prevent microbial transfer.
- Folded and packaged for easy, aseptic application by one person.
- Can be included in Mölnlycke® ProcedurePak® trays to support operating room efficiency.

When to use BARRIER aperture drape:

You can use aperture drapes for a variety of procedures such as ophthalmic surgery and minor surgical interventions.

### Images

---




Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.





## Delivered items

63623-00

**Country of origin:** China**Shelf life:** 5 years**Sterilization method:** EtO**Production Responsibility:** Xuchang Zhengde Environstar Medical Products Co. Ltd., Area B, Baihua Road East, South Side of Erji river, Industry Cluster District, Yanling County, 461200 Xuchang, Henan, China**Packing information:** First packaging layer is a peel-open sterile barrier, papaer/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.**Is suitable for Tray:** No

Packing level	Quantity	GTIN Code	UDI-DI	Width x Length x Height	Volume	Weight gross / net
Piece	1	7323190241583			186.667 cm <sup>3</sup>	29.944 g / -
Consumer pack	1	6971004780868	6971004780868	290x150x0 mm		
Dispenser box	90	6971004780875	6971004780875	140x300x330 mm		
Transport box	180	6971004780882	6971004780882	300x320x350 mm	33.6 dm <sup>3</sup>	5.39 / 4.89 kg
Pallet	7200	7323190241590		800x1200x1880 mm	1804.8 dm <sup>3</sup>	240.6 kg / 64.8

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.  
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.



## Material

<b>Animal tissues:</b>	No
<b>Human blood derivatives:</b>	No
<b>Natural rubber latex:</b>	No
<b>Medicinal substances:</b>	No
<b>Phthalates:</b>	No
<b>Polyvinyl chloride:</b>	No

## Technical

## Dimension

Dimension text	Dimension value
Outer dimension	75 cm x 90 cm
Aperture	6 cm x 8 cm

## Instructions

## How product works:

The surgical patient drapes are intended to be used as protective patient coverings to create a sterile field for a surgical procedure.

## Classifications

Regulation type(s)	MDR Class IS
Intended Purpose:	Used as a protective covering or a barrier to isolate a site of surgical incision or a surgical field from contamination found in operating theatre or catheterization laboratory e.g. on instruments or tables.
Notified body medical devices/PPE:	Notified Body (see CE certificate)

**Applied standards: The standards presented below is a selection of the most essential standards that are adhered to.**

EN ISO 15223-1, EN ISO 14971, EN ISO 13485, EN 13795-1

## Removable label

Yes

## GMDN Code (Global Medical Device Nomenclature)

47783 Patient surgical drape, single-use

## EMDN Code (European Medical Device Nomenclature)

T02010101 INCISION DRAPES, WITHOUT ANTIBACTERIAL AGENT

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.  
©2025 Mölnlycke Health Care AB. All rights reserved.



## UNSPSC

42131701 Surgical drapes

## Commodity Code

3005100000 Wadding, Gauze, Dressings, drapes singlepacked - adhesive articles, Sets mainly consisting thereof

## CE Responsibility / Legal Manufacturer

Xuchang Zhengde Environstar Medical Products Co., Ltd. Area B, Baihua Road East, South Side of Erji river, Industry Cluster District, Yanling County, 461200 Xuchang, Henan, China

## Basic UDIDI

69317478SD000001GC

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.  
©2025 Mölnlycke Health Care AB. All rights reserved.

