



Actico® UlcerSys medical compression stocking system

REF 32510 – 32519, 39710 – 39719

February 2025

1. Composition of the product

The Actico® UlcerSys medical compression stocking system consists of:

- Two liners (white / sand-coloured), consisting of:
 - 70% Polyamide (nylon)
 - 30% Elastane (LYCRA®)
- One top stocking (sand-coloured / black), consisting of:
 - 65% Polyamide (nylon)
 - 35% Elastane (LYCRA®)

This product data sheet covers the following items:

REF	Product designation	Colour	Size
32510	Actico® UlcerSys medical compression stocking system standard	sand/white	S
32511	Actico® UlcerSys medical compression stocking system long	sand/white	S
32512	Actico® UlcerSys medical compression stocking system standard	sand/white	M
32513	Actico® UlcerSys medical compression stocking system long	sand/white	M
32514	Actico® UlcerSys medical compression stocking system standard	sand/white	L
32515	Actico® UlcerSys medical compression stocking system long	sand/white	L
32516	Actico® UlcerSys medical compression stocking system standard	sand/white	XL
32517	Actico® UlcerSys medical compression stocking system long	sand/white	XL
32518	Actico® UlcerSys medical compression stocking system standard	sand/white	XXL
32519	Actico® UlcerSys medical compression stocking system long	sand/white	XXL
39710	Actico® UlcerSys medical compression stocking system standard	black/sand	S

REF	Product designation	Colour	Size
39711	Actico® UlcerSys medical compression stocking system standard	black/sand	S
39712	Actico® UlcerSys medical compression stocking system long	black/sand	M
39713	Actico® UlcerSys medical compression stocking system standard	black/sand	M
39714	Actico® UlcerSys medical compression stocking system long	black/sand	L
39715	Actico® UlcerSys medical compression stocking system standard	black/sand	L
39716	Actico® UlcerSys medical compression stocking system long	black/sand	XL
39717	Actico® UlcerSys medical compression stocking system standard	black/sand	XL
39718	Actico® UlcerSys medical compression stocking system long	black/sand	XXL
39719	Actico® UlcerSys medical compression stocking system standard	black/sand	XXL
	Actico® UlcerSys medical compression stocking system long	black/sand	XXL

2. Packaging, structure and composition

2.1 Unit Container

- see 2.2

2.2 Shelf Container

1 package = 1 folding box (cellulose)

- instructions for use (cellulose)

2.3 Transit Container

16 packages = 16 folding boxes = 1 corrugated cardboard box (cellulose)

3. Manufacture

The Actico® UlcerSys medical compression stocking system is manufactured to specification under hygienic conditions and packaged as described in their respective packaging specifications.

4. Description

The Actico® UlcerSys medical compression stocking system is a kit containing two white / sand-coloured liners and one sand-coloured/ black top stocking.

5. Properties

Two-layered ready-for-use medical compression stocking with a pressure at ankle of approximately 40 mmHg. The foot of the liner is closed, the foot of the top-stocking is open. The stockings can be washed according to the marking or according to the instructions for care given in the IFU.

Natural rubber latex is not an intentional ingredient in the production of the above mentioned products or its packaging. However, to the best of our knowledge we can't completely exclude any potential contaminations with latex traces from the production or the environment.

6. Intended purpose

The Actico® UlcerSys medical compression stocking system is intended to be used for lower leg ulcer therapy over the treated wound.

7. Medical device classification

The Actico® UlcerSys medical compression stocking system is a medical device of Class I in terms of Rule 1 (according to the Regulation (EU) 2017/745 concerning medical devices (MDR)).

8. Biological evaluation and biocompatibility

The biological risk evaluation performed according to ISO 10993-1 confirms the safety and acceptable biocompatibility of the device Actico® UlcerSys medical compression stocking system within its intended purpose.

9. Stability

Stored appropriately (keep dry, at room temperature, protected from dust, moisture and direct sunlight), Actico® UlcerSys medical compression stocking system have a shelf life of 5 years from the date of manufacture.

10. Disposal

The user is advised to observe current national legislation, norms, and guidelines regulating the disposal of medical waste.
Packaging materials must also be disposed of in compliance with applicable national requirements.

Lohmann & Rauscher International GmbH & Co. KG
D-56579 Rengsdorf
signed by

Michael Hamedinger
Head of Global Marketing Bandages & Supports

Loay Karam
Global Brand Manager Bandages & Supports