

Product Specifications

Brand : Klinion®
Product : Blunt Fill Needles (with & without filter)

Product description

Klinion® Blunt Fill Needles are disposable needles characterized by their role in reducing needle stick injuries while preparing medicine. They are designed with a special needle tip to prevent rubber fragments. Available with and without a 5-micron filter to prevent particles contamination. Compatible with Luer connector.

Intended use

Klinion® Blunt Fill Needles are intended for drawing up medicine from a vial or ampoule before injection, by healthcare professionals.

Classification

Class Is under the Medical Devices Regulations (SI 2002 No. 618, as amended) and
Class Is under Regulation (EU) 2017/745 Medical Device Regulation (MDR).

Packaging & Product Reference

Product Reference	Description	Colour Code	Needle Gauge	Needle Length	Piece (s) per Peel Pouch	Pieces per Inner Carton	Pieces per Outer Carton
C0880440	KLINION BLUNT FILL NEEDLE WITH 5µ FILTER	Purple	18G	40mm	1	100	10 000 (100 x 100)
C0880441	KLINION BLUNT FILL NEEDLE WITHOUT FILTER	Pink	18G	40mm	1	100	10 000 (100 x 100)

Product Reference	Description	Colour Code	Needle Gauge	Needle length	EAN Peel Pouch	EAN Inner Carton	EAN Outer Carton
C0880440	KLINION BLUNT FILL NEEDLE WITH 5µ FILTER	Purple	18G	40mm	06971564462627	16971564462621	26971564462628
C0880441	KLINION BLUNT FILL NEEDLE WITHOUT FILTER	Pink	18G	40mm	06971564462631	16971564462638	26971564462635

Product Specifications

Features

- Special needle tip designed to prevent rubber fragment contamination.
- Compatible with ISO 80369-7 luer connector.

Contra-Indications

Do not use Klinion® Blunt Fill Needles for skin injection. Only for medicine preparation.

Precautions

- Sterile if pouch and seal are intact.
- Do not use if packaging is damaged or not intact.
- Do not use after expire date.
- Do not re-sterilize.
- If the needle is blocked or damaged, please do not use it.
- Do not re-use, for single use-only. Re-use may lead to contamination.
- For use under the guidance of a healthcare professional.

Presence of specific components

Component	Present
Latex	No
Phthalates	No
Substances of animal origin	No
Substances of human origin	No
Blood or blood derivates	No
BPA	No
Heavy metals	No

Sterilization

EO Sterilized.

Product Specifications

Storage conditions

Klinion® Blunt Fill Needles should be stored in a cool, dry environment away from direct sunlight, high heat, and high humidity.

Shelf life

Klinion® Blunt Fill Needles have a shelf-life of 3 years from production date, when stored under specified conditions.

Picture



Product ref.: C0880441



Product ref.: C0880440

Approval	Name	Date	Signature
Associate Category Management Director	Susanne van den Berg		
RA Specialist	Alexandra Maria da Silva Duarte		