

Product Specifications

Brand : Klinion®

Product : Klinifix® Cohesive



Ref: 132387

Product description

Klinifix® Cohesive is an elastic, comfortable, breathable, self-adhesive bandage. Due to the micro-structured coating, the bandage adheres only to itself—not to skin, hair, or clothing—eliminating the need for additional tapes, pins, or clips. The breathable, air-permeable structure promotes patient comfort and helps maintain skin integrity.

Intended use

Secure fixation of wound dressings on mobile or round body parts. Remains in place by adhering to itself. Its high elasticity and conformability ensure a smooth, crease-free application that adapts seamlessly to mobile or rounded body parts such as joints, providing lasting fixation.

Classification

Class I under the Medical Devices Regulations (SI 2002 No. 618, as amended) and Class I under Regulation (EU) 2017/745 Medical Device Regulation (MDR).

Packaging & Product Reference

| Product Reference | Description | Size (*) | Color | Piece(s) | Inner Carton (pcs) | Outer Carton (pcs) |
|-------------------|-------------------|----------|-------|----------|--------------------|--------------------|
| 132377 | KLINIFIX COHESIVE | 4MX2,5CM | White | 1 | 10 | 650 (65 x 10) |
| 132378 | KLINIFIX COHESIVE | 4MX4CM | White | 1 | 5 | 280 (56 x 5) |
| 132379 | KLINIFIX COHESIVE | 4MX6CM | White | 1 | 5 | 280 (56 x 5) |
| 132380 | KLINIFIX COHESIVE | 4MX8CM | White | 1 | 5 | 140 (28 x 5) |
| 132381 | KLINIFIX COHESIVE | 4MX10CM | White | 1 | 5 | 140 (28 x 5) |
| 132382 | KLINIFIX COHESIVE | 4MX12CM | White | 1 | 5 | 140 (28 x 5) |
| 132384 | KLINIFIX COHESIVE | 20MX4CM | White | 1 | 1 | 24 (24 x 1) |
| 132386 | KLINIFIX COHESIVE | 20MX6CM | White | 1 | 1 | 24 (24 x 1) |
| 132376 | KLINIFIX COHESIVE | 20MX6CM | Blue | 1 | 1 | 24 (24 x 1) |

Product Specifications

| | | | | | | |
|--------|-------------------|----------|-------|---|---|-------------|
| 132383 | KLINIFIX COHESIVE | 20MX6CM | Red | 1 | 1 | 24 (24 x 1) |
| 132387 | KLINIFIX COHESIVE | 20MX8CM | White | 1 | 1 | 24 (24 x 1) |
| 132388 | KLINIFIX COHESIVE | 20MX10CM | White | 1 | 1 | 24 (24 x 1) |
| 132389 | KLINIFIX COHESIVE | 20MX12CM | White | 1 | 1 | 24 (24 x 1) |

| Product Reference | Description | Size (*) | Color | Piece(s) | EAN Inner Carton | EAN Outer Carton |
|-------------------|-------------------|----------|-------|----------|------------------|------------------|
| 132377 | KLINIFIX COHESIVE | 4MX2,5CM | White | 1 | 8715343029671 | 8715343029688 |
| 132378 | KLINIFIX COHESIVE | 4MX4CM | White | 1 | 8715343060018 | 8715343008270 |
| 132379 | KLINIFIX COHESIVE | 4MX6CM | White | 1 | 8715343060087 | 8715343051597 |
| 132380 | KLINIFIX COHESIVE | 4MX8CM | White | 1 | 8715343060155 | 8715343051610 |
| 132381 | KLINIFIX COHESIVE | 4MX10CM | White | 1 | 8715343060223 | 8715343051641 |
| 132382 | KLINIFIX COHESIVE | 4MX12CM | White | 1 | 8715343060292 | 8715343051665 |
| 132384 | KLINIFIX COHESIVE | 20MX4CM | White | 1 | 8715343745229 | 8715343745212 |
| 132386 | KLINIFIX COHESIVE | 20MX6CM | White | 1 | 8715343060346 | 8715343008355 |
| 132376 | KLINIFIX COHESIVE | 20MX6CM | Blue | 1 | 8715343745274 | 8715343745250 |
| 132383 | KLINIFIX COHESIVE | 20MX6CM | Red | 1 | 8715343745243 | 8715343745236 |
| 132387 | KLINIFIX COHESIVE | 20MX8CM | White | 1 | 8715343060414 | 8715343008263 |
| 132388 | KLINIFIX COHESIVE | 20MX10CM | White | 1 | 8715343060483 | 8715343008362 |
| 132389 | KLINIFIX COHESIVE | 20MX12CM | White | 1 | 8715343745205 | 8715343745199 |

*Sizes are approximate

Features

- Elastic
- Breathable
- Self-adhesive
- Single use
- Latex free

Product Specifications

Precautions

- Non-sterile, do not reuse
- Do not apply with too much tension
- Do not use directly on open wounds

Specifications

| Ref | Size | Material | Extensibility of the material | Stretched length | Width | Weight |
|--------|----------|--|-------------------------------|------------------|-------------|----------------------|
| 132377 | 4mx2,5cm | Cotton 32%, viscose 25%, polyamide 23%, synthetic rubber 20% | 100% ± 20% | 4m±0.2m | 2.5cm±0.2cm | 39±5g/m ² |
| 132378 | 4mx4cm | Cotton 32%, viscose 25%, polyamide 23%, synthetic rubber 20% | 100% ± 20% | 4m±0.2m | 4cm±0.2cm | 39±5g/m ² |
| 132379 | 4mx6cm | Cotton 32%, viscose 25%, polyamide 23%, synthetic rubber 20% | 100% ± 20% | 4m±0.2m | 6cm±0.2cm | 39±5g/m ² |
| 132380 | 4mx8cm | Cotton 32%, viscose 25%, polyamide 23%, synthetic rubber 20% | 100% ± 20% | 4m±0.2m | 8cm±0.2cm | 39±5g/m ² |
| 132381 | 4mx10cm | Cotton 32%, viscose 25%, polyamide 23%, synthetic rubber 20% | 100% ± 20% | 4m±0.2m | 10cm±0.2cm | 39±5g/m ² |
| 132382 | 4mx12cm | Cotton 32%, viscose 25%, polyamide 23%, synthetic rubber 20% | 100% ± 20% | 4m±0.2m | 12cm±0.2cm | 39±5g/m ² |
| 132384 | 20mx4cm | Cotton 32%, viscose 25%, polyamide 23%, synthetic rubber 20% | 100% ± 20% | 20.0m±1.0m | 4cm±0.2cm | 39±5g/m ² |
| 132386 | 20mx6cm | Cotton 32%, viscose 25%, polyamide 23%, synthetic rubber 20% | 100% ± 20% | 20.0m±1.0m | 6cm±0.2cm | 39±5g/m ² |
| 132376 | 20mx6cm | Cotton 32%, viscose 25%, polyamide 23%, synthetic rubber 20% | 100% ± 20% | 20.0m±1.0m | 6cm±0.2cm | 39±5g/m ² |
| 132383 | 20mx6cm | Cotton 32%, viscose 25%, polyamide 23%, synthetic rubber 20% | 100% ± 20% | 20.0m±1.0m | 6cm±0.2cm | 39±5g/m ² |
| 132387 | 20mx8cm | Cotton 32%, viscose 25%, polyamide 23%, synthetic rubber 20% | 100% ± 20% | 20.0m±1.0m | 8cm±0.2cm | 39±5g/m ² |
| 132388 | 20mx10cm | Cotton 32%, viscose 25%, polyamide 23%, synthetic rubber 20% | 100% ± 20% | 20.0m±1.0m | 10cm±0.2cm | 39±5g/m ² |
| 132389 | 20mx12cm | Cotton 32%, viscose 25%, polyamide 23%, synthetic rubber 20% | 100% ± 20% | 20.0m±1.0m | 12cm±0.2cm | 39±5g/m ² |

All sizes have a folded edge, except for the 4mx2,5cm which has a paper-inlay as end-finder.

Product Specifications

Presence of specific components

| Component | Present |
|-----------------------------|---------|
| Latex | No |
| Phthalates | No |
| Substances of animal origin | No |
| Substances of human origin | No |
| Blood or blood derivatives | No |
| Medicinal substances | No |
| Heavy metals | No |

Sterilization

Non-sterile

Storage conditions

The products should be stored at dry conditions and protected from direct sunlight.

Shelf life

Klinifix® Cohesive has a shelf life of 3 years from production date, when stored under specified conditions.

| Approval | Name | Date | Signature |
|-------------------------------|--------------|------|-----------|
| Product Manager Wound Care | Florin Lazar | | |
| RA Specialist | Emil Laske | | |