


## 1) Product Content and Technical Drawing



1x Omslagsduk - Sterilization Wrap

40x50 cm

1x Skål Plast 3-delt- Tray 3 comp.

15x7,5x2 cm

10x Tupfer Non-woven - Nonwoven Tampon


ø30 mm

1x Pinsett m/Styrespor Plast - Plastic Forcep Ref: 489-365

13 cm

**Packagin information :**

**Packaging Details: Soft Blister**



**Photos:**




## 2) Product Labeling

### Paket Etiketi / Package Label


140 mm

# submed Safe Easy KIT

## Vaske-/skiftesett Sår Lite / Dressing Kit



- 1x Omslagsduk / Sterilization Wrap  
40x50 cm
- 1x Skål Plast 3-delt- Tray 3 comp.  
7,5x15x2
- 10x Tupfer Non-woven - Nonwoven Tampon  
ø30 mm
- 1x Pinsett m/Styrespor Plast - Plastic Forcep  
13 cm



YYYY/MM


LOT XXXXXXX

YYYY/MM

REF SKT-007-1



(01)08682780972726(17)YYMM00(30)1(10)XXXXXX



STERILE EO

MD



1984



Please scan the QR code for instructions for usage. If you need paper version please contact with the manufacturer.  
Scan QR-koden for brugsanvisning. Hvis du har brug for papirversion, kontakt venligst producenten.

Su Biyomedikal Sistemler ve Saglik Hiz. San. Tic. A.Ş. Orhangazi Mah. 1673. Sok. No:20  
Esenyurt/Istanbul/Turkey • Tel: +90 212 320 37 53 • info@submed.com.tr • www.submed.com.tr

115 mm

Barcode : 8682780972726

## 4) Instructions for Use



SUBMED – Su Biyomedikal Sistemler ve Saglik Sistemleri Sanayi Ticaret A. Ş.  
Orhangazi Mah. 1673 Sok. No:20 Esenyurt /Istanbul/ Türkiye  
Tel: (+ 90 212) 320 37 53 Fax: (+ 90 212) 320 53 51 e-mail: [info@submed.com.tr](mailto:info@submed.com.tr) [www.submed.com.tr](http://www.submed.com.tr)

### INSTRUCTIONS FOR USE

**Product Name:** Dressing Kit

**Reference Number:** SKT-007-1

#### Usage Instruction:

Sterilization wrap is used for bundling the materials together. Nonwoven tampon is used to clean and close wounds. The plastic 3 compartments tray is used to prepare the solution. Plastic tweezers used to hold swabs and tampons.

#### Symbols

	Date of manufacture
	Expiry date
	Reference number
	Lot number
	Sterile, Ethylene Oxide sterilization
	Do not re-use
	Do not re-sterile
	Do not use if package is damaged
	Keep dry
	Keep away from sunlight
	Manufacturer information
	CE mark and Notified Body number
	Consult instruction for use
	Medical Device

#### Intended Use

Surgical sets are products for ensure that the materials used by nurses and doctors in their standard procedures during the operation are sterile and combined in hospitals and clinics where surgical operations are performed. They are medical devices that prevent diseases by reducing the risk of cross contamination.

#### User and Patient Information

There is no age and gender restriction on the patient to be used. The products are used by physicians and healthcare professionals.

#### Storage

Store the product in a cool and dry place protected from light under normal conditions. The use of medical paper and film guarantees sterility during shelf life under normal transport and storage conditions.

#### Destruction

Product contains fabric, plastic, metallic components. These are considered as medical waste and should be treated according to national and regional laws.

#### Warnings

Product should only be used by qualified personnel. Check the integrity of the package before use. Do not use if the package is opened or damaged. Product is only for single use. Store under dustless conditions.

#### Indications



Medical devices that prevent diseases by reducing the risk of cross contamination.

#### Contraindications


There are no identified disadvantageous.

### 3)Boxing and Palleting Details

#### 3.1) Desing



Inner Carton Label



Outer Carton Label

#### 3.2) Details

Pallet/ Box Qty	32 boxes of 100 pcs = 3200 pcs
Box Qty	100 pcs.
Inner Box Qty	25 pcs
Box Dimensions	600x400x250 mm - 190x290x250 mm
Pallet Dimensions	800x1200 mm