

SPECIFICATION

servoprax Tourniquet



- Made of strong material, with a plastic fastener
- Single-handed use: Can be applied and removed using only one hand
- Once in place, it cannot open unintentionally
- Strap can be washed at 60 °C
- Packed singly in a display package Usage is simple!

- 1 Place the strap around the arm and insert the clip into the fastener until it "clicks" into place, indicating that it is closed.
- 2 To tighten, pull the end of the strap in the direction of the arrow. The fastener is securely locked.
- 3 Once the needle is inserted, press the tension button 1 to gradually loosen the tourniquet.
- 4 To remove the tourniquet, press the release button 2.

SPECIFICATION

servoprax Tourniquet

Type	Length	PZN	REF	Unit
Green	40 cm		G2 1001G	1 piece
Children's motif	40 cm		G2 1001KM	1 piece
Red	40 cm		G2 1001R	1 piece
Blue	40 cm		G2 1001	1 piece

