


Sol-M™ 3-piece Syringes

Luer Lock & Luer Slip without Needle

Sol-M™ 3-Piece Syringes accurately measure the quantity of liquid medications (according to ISO 7886-1) and are available in a wide range of sizes and connection combinations to address a wide variety of clinical needs.

Features:


- High-transparency barrel facilitates seeing syringe content
- Silicone lubrication ensures smooth plunger movement
- Syringe plunger backstop – helps to prevent accidental withdrawal of plunger
- Black graduation scale – designed to enhance contrast and readability
- Available in Luer Lock and Luer Slip tip
- Large finger grips support clinicians' better grip control during injections
- Designed according to ISO 7886-1
- Not manufactured with Latex, PVC and DEHP



CE 2797

Sol-M™ 3-piece Syringes

Luer Lock without Needle



Sol-M™ Luer Lock Syringe without Needle

Ref. Number	Volume	Box Qty	Case Qty
P180001PP*	1 ml	100	800
180001**	1 ml	100	800
P180001***	1 ml	100	800
P180003	3 ml	100	800
P180005	5 ml	100	800
P180010	10 ml	100	1200
P180020	20 ml	100	800
P180030	30 ml	50	400
P180060	60 ml	30	240

*Low Dead Space feature (polypropylene barrel), **Polycarbonate barrel, ***Polypropylene barrel

Sol-M™ 3-piece Syringes

Luer Slip without Needle



Sol-M™ Slip Tip Syringe without Needle - Luer Slip

Ref. Number	Volume	Box Qty	Case Qty
P180011	1 ml	100	800
P180011LDS*	1 ml	100	800
P180003ST	3 ml	100	800
P180005ST	5 ml	100	800
P180010ST	10 ml	100	1200
P180020ST	20 ml	100	800

* Low Dead Space feature

Sol-M™ Slip Tip Syringe without Needle - Eccentric Tip

Ref. Number	Volume	Box Qty	Case Qty
P180010ET	10 ml	100	1200
P180020ET	20 ml	100	800
P180060ET	60 ml	30	240

Sol - M™ 3-piece Syringes



Capacity	Luer Lock tip	Luer Slip tip	Eccentric tip
1 ml	P18001PP* 180001** P180001***	P180011 P180011LDS*	
3 ml	P180003	P180003ST	
5 ml	P180005	P180005ST	
10 ml	P180010	P180010ST	P180010ET
20 ml	P180020	P180020ST	P180020ET
30 ml	P180030		
60 ml	P180060		P180060ET



*Syringe with Low Dead Space feature (polypropylene barrel), **Polycarbonate barrel, ***Polypropylene barrel

Intended use: The Standard Syringe without Needle is used to inject fluids into, or withdraw fluids from, the body.

Intended user: Licensed healthcare professionals (HCP).


Cautions:

- Single use device. Re-use or use if the package is damaged may lead to infection or other illness/ injury.

The content of this material is directed at and restricted to a professional audience only (Healthcare Professionals). This material is not medical advice, guideline or a recommendation. Before using any particular product always read the package insert, product label and/or instructions for use. This material is distributed by Sol-Millennium Europe Sp z o.o. with a registered address at Twarda 18, 00-105 Warsaw, Poland, Registry Court: District Court for Warsaw, XIIIth Economic Department KRS 0000434783, EU VAT: PL5252535720, Share capital: 500.000 PLN. Sol-Millennium Europe Sp. z o.o. All rights reserved.

EMA-MD-451-0423

For more information, please contact:

 Sol-Millennium Europe Sp. z o.o.
Warsaw HUB, Rondo Daszyńskiego 2B
00-843 Warsaw, Poland

 info.eu@sol-m.com

 **Sol M** Building a Healthier Tomorrow