



GLORIOUS ENTERPRISES

PRODUCT SHEET

Products Name: KIDNEY TRAYS 170X100X35MM



Glorious Art nr.	60.543.00
Net weight:	210 gram (approximate)
Material:	Stainless steel AISI 202
Finish:	BRIGHT
Sterility:	Non-Sterile (Re-Usable)
Exposure Temperature:	134°C (274°F)
Exposure Time:	3 Minutes
CE Classification:	Class-I
Storage:	16°C to 30°C (Store the instruments in a closed cabinet, protected from dust, moisture and temperature fluctuations)
Shelf life:	3 years
life-span:	20 times of use:- after each use, product need to sterilize for ready to use in next operation or next surgery.
Disposal:	No special decomposition or disposal required for Surgical Instruments as they don't contain any toxic or hazardous material



P.O. Box 63, Naverland 11, DK 2600 Glostrup, Denmark
 Tel: +45 43 62 26 93, Fax: +45 43 62 26 11 , Email: glorious@glorious-enterprises.dk