


TECHNICAL SHEET

23/11/2017

Reference	IVS 01, IVS03, IVS05, IVS10, IVS20, IVS30, IVS60
Trade Name:	Luer Slip IV syringes
Product Description	Luer Slip IV syringes
Intended Use	Sterile syringes can be used together with hypodermic needles, and are intended to be used intravenously, to deliver drugs or take blood from the patient
Raw Materials:	<ol style="list-style-type: none"> 1. Plunger in polypropylene 2. Barrel in polypropylene 3. Piston in polyisoprene <p style="text-align: center;">Lubricant : silicone oil</p>
	
Packaging	Paper/PP-PE film Blister
Shelf Life:	5 years
Usage	Single use device
Sterilisation	EO Sterilisation
	Device should not be re-sterilised
Storage	Normal storage procedures and conditions
Disposal	The user must follow the legal regulations regarding disposal of hospital waste.
Manufacturer	Medicina Ltd Units 1-4, Rivington View Business Park Station Road, Blackrod Bolton BL6 5BN United Kingdom
Classification	Class IIa
Conformity Assessment Route:	MDD 93/42/EEC, Annex V & VII
Notified Body	SGS Ltd.,UK. NB Designation: 0120
Certificate Ref	CE GB12/84913