

Cornea-Care™ Eye Protection Dressing

ANAESTHESIA

Reduces the risk of corneal drying and corneal abrasion during general anaesthesia or sedation

'No touch' application technique avoids cross-contamination of fluids

Easy to apply and easy to remove

Latex free

One size

Reduces the loss of eyebrow and eyelash hair upon removal

Transparent windows allow the eyes to be observed



Mild adhesive on the transparent window ensures proper closure of eyelid margins

CORNEA-CARE™ EYE PROTECTION DRESSING

Reference	Unit	Units/Box
CCA-01(500)	Pair	500



BLUE BOX
MEDICAL LTD

Unit 29, New Forest Enterprise Centre,
Chapel Lane, Totton, Southampton SO40 9LA



Tel: 023 8066 9000 Fax: 023 8066 9002
enquiries@blueboxmedical.co.uk
www.blueboxmedical.co.uk



PS027

06/14 - A