


699054 Universal Set Standard

Product details

Sterile: Sterile

Images



Products included in set

Product ID	Quantity	Product name	Size	Descriptive feature
996106	1,0	Mayo Stand Cover	79cm x 145cm	Absorbent
910200	2,0	Adhesive OP-Towel	90cm x 75cm	3-ply, Adhesive, Patched
26153-20	4,0	Cellulose Towel	18cm x 25cm	Absorbent
911500	1,0	OP-Tape	9cm x 49cm	Adhesive
967772	1,0	Adhesive OP-Sheet	175cm x 175cm	3-ply, Adhesive, Patched
967773	1,0	Adhesive OP-Sheet	240cm x 150cm	3-ply, Adhesive, Patched
90190	1,0	Table Cover	150cm x 190cm	Absorbent

Delivered items

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2019 Mölnlycke Health Care AB. All rights reserved.

699054-07

Country of origin: Belgium**Shelf life:** 5 years**Sterilization method:** Irradiation**Packing information:** First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.**Is suitable for Tray:** No

Packing level	Quantity	GS1 code
Consumer pack	1	7332430983794
Dispenser box	6	7332430983787
Transport box	12	7332430983770
Pallet	144	7332430983763

Material

Animal tissues:	No
Human blood derivatives:	No
Natural rubber latex:	No
Medicinal substances:	No
Phthalates:	No
Polyvinyl chloride:	No

Technical

Disposal instructions

Non-hazardous waste used BARRIER products and sterility barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. BARRIER products do not contain any hazardous substances that can leach out if the products are land filled. Transport boxes are designed to fit existing recovery systems. The new BARRIER packaging system complies with the Packaging Waste Directive of the European Union.

Storage instructions

Mölnlycke Health Care recommends that BARRIER products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterility barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the package is opened and the product is used.

Classifications

Regulation type(s)	MDD Class IS	Locally Regulated	Unregulated
--------------------	--------------	-------------------	-------------

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2019 Mölnlycke Health Care AB. All rights reserved.



Applied standards: The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 13795, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Removable label

Yes

GMDN Code (Global Medical Device Nomenclature)

33961 General/plastic surgical procedure kit, non-medicated, single-use

Find out more at www.molnlycke.com


Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2019 Mölnlycke Health Care AB. All rights reserved.



Material information for the ingoing products

699054 Universal Set Standard, 996106 Mayo Stand Cover

- : 79cm x 145cm
- : Absorbent



- : No
- : No
- : No
- : No
- : No
- : No

Outer dimension	79 cm x 145 cm
Absorbent	65 cm x 85 cm

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2019 Mölnlycke Health Care AB. All rights reserved.



699054 Universal Set Standard, 910200 Adhesive OP-Towel

- : 90cm x 75cm, Patch 50 cm x 20 cm
- : 3-ply, Adhesive, Patched



- : No
- : No
- : No
- : No
- : No
- : No

Outer dimension	90 cm x 75 cm
Patch	50 cm x 20 cm

Find out more at www.molnlycke.com


Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2019 Mölnlycke Health Care AB. All rights reserved.



699054 Universal Set Standard, 911500 OP-Tape

: 9cm x 49cm

: Adhesive



: No

: No

: No

: No

: No

Outer dimension

9 cm x 49 cm

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2019 Mölnlycke Health Care AB. All rights reserved.

699054 Universal Set Standard, 967772 Adhesive OP-Sheet

- : 175cm x 175cm, Patch 20 cm x 55 cm
- : 3-ply, Adhesive, Patched



- : No
- : No
- : No
- : No
- : No
- : No

Outer dimension	175 cm x 175 cm
Patch	20 cm x 55 cm

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2019 Mölnlycke Health Care AB. All rights reserved.



699054 Universal Set Standard, 967773 Adhesive OP-Sheet

- : 240cm x 150cm, Patch 55 cm x 20 cm
- : 3-ply, Adhesive, Patched



- : No
- : No
- : No
- : No
- : No
- : No

Outer dimension	240 cm x 150 cm
Patch	55 cm x 20 cm

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2019 Mölnlycke Health Care AB. All rights reserved.



699054 Universal Set Standard, 90190 Table Cover

- : 150cm x 190cm
- : Absorbent



- : No
- : No
- : No
- : No
- : No

Outer dimension	150 cm x 190 cm
Absorbent	75 cm x 190 cm

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2019 Mölnlycke Health Care AB. All rights reserved.

