



# Biosyn™ suture

## PRODUCT INFORMATION

Structure	Monofilament
Suture Type	Mid-Term Absorbable
Composition	Synthetic polyester composed of glycolide, dioxanone, and trimethylene carbonate
Coating	Uncoated
Color	Violet, Undyed
Tensile Strength	2 weeks: 75% USP; 3 weeks: 40% USP
Absorption Profile	90-110 days
Sizes	6-0 to 1
Indications	Biosyn™ synthetic absorbable suture is indicated for use in general soft tissue approximation and/ or ligation including use in ophthalmic surgery, but not for use in cardiovascular or neurological surgery.
Contraindications	This suture, being absorbable, should not be used where extended approximation of tissue is required.
Box Quantities	1, 2 or 3 dozen

## ORDER INFORMATION

Biosyn™ codes can be identified with codes beginning with:  
**M, SM, CM, XM, UM, GM, GMT, CMT**

For more information, please visit  
[medtronic.eu/product-catalog](https://www.medtronic.eu/product-catalog)

**IMPORTANT:** Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 17-weu-biosyn-tech-sheet-en-1855595

**Medtronic**  
Further, Together