


650104 Surgical Gown CLASSIC

Surgical Gown

Product details

Product group name :	CLASSIC, Standard performance
Size :	L-L
Descriptive feature :	Alcohol repellent, Antistatic treated, Standard performance, Two towels, Wrapped
Colour :	Ocean
Sterile :	Sterile
Brand :	BARRIER®

Images



Delivered items

650104-00

Country of origin: Thailand

Shelf life: 5 years

Sterilization method: Irradiation

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, Amata Nakorn (Bang Pakong), Industrial Estate, 700/461 Moo 7 Bangna-Trad Road, Km 57 Tambol Donhuaroh, Amphur Muang Chonburi 20000, Thailand

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.

The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2024 Mölnlycke Health Care AB. All rights reserved.



opened the sterile barrier cannot be closed again. Second layer is a plastic liner / bag. Third layer is a corrugated board dispenser box. Fourth layer is a corrugated board transport box.

Is suitable for Tray: Yes

Packing level	Quantity	GTIN Code	UDI-DI	Width x Length x Height	Volume	Weight gross / net
Piece	1	7332430791474			1893.9 cm ³	272.7 g / -
Consumer pack	1	7332430875631	7332430875631			
Dispenser box	20	7332430791504	7332430791504	279x310x384 mm		
Transport box	40	7332430791498	7332430791498	397x580x329 mm	75.8 dm ³	10.9 / 8.5 kg
Pallet	800	7332430791481		800x1200x1789 mm		

Technical

Dimension

Dimension reference	Dimension text	Dimension value
A	Shoulder	22 cm
B	Chest (most narrow part at front)	54 cm
C	Top sleeve	58 cm
D	Armhole	29 cm
E	Cuff length	7 cm
F	Neck	66 cm
G	Length, highest point at shoulder	131 cm
H	Circumference	154 cm

Disposal instructions

Non-hazardous waste used BARRIER products and sterility barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. BARRIER products do not contain any hazardous substances that can leach out if the products are land filled. Transport boxes are designed to fit existing recovery systems. The new BARRIER packaging system complies with the Packaging Waste Directive of the European Union.

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.

The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2024 Mölnlycke Health Care AB. All rights reserved.



Storage instructions

Mölnlycke Health Care recommends that BARRIER products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterility barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the package is opened and the product is used.

Classifications

Regulation type(s)	MDR Class IS	Locally Regulated	Locally Regulated	Locally Regulated
CE Certificate Number :	MDR 722028			
Intended Purpose :	The Surgical Gowns are single use, disposable, fluid repellent garments intended to be used as sterile by operating room personnel during surgical procedures to protect both the patient and the operating room personnel (users) from the transfer of microorganisms, body fluids and particulate material.	The Surgical Gowns are single use, disposable, fluid repellent garments intended to be used as sterile by operating room personnel during surgical procedures to protect both the patient and the operating room personnel (users) from the transfer of microorganisms, body fluids and particulate material.	The Surgical Gowns are single use, disposable, fluid repellent garments intended to be used as sterile by operating room personnel during surgical procedures to protect both the patient and the operating room personnel (users) from the transfer of microorganisms, body fluids and particulate material.	The Surgical Gowns are single use, disposable, fluid repellent garments intended to be used as sterile by operating room personnel during surgical procedures to protect both the patient and the operating room personnel (users) from the transfer of microorganisms, body fluids and particulate material.
MDR Classification Rule :	1			
Conformity Annexes :	IX			
Measuring Function :	No			
Notified body medical devices/PPE :	BSi 2797	Local Authority MEA	Local Authority EU Non EEC	Local Authority APAC

Find out more at www.molnlycke.com



Classifications

Regulation type(s)	Locally Regulated	Locally Regulated	Locally Regulated
CE Certificate Number :			
Intended Purpose :	The Surgical Gowns are single use, disposable, fluid repellent garments intended to be used as sterile by operating room personnel during surgical procedures to protect both the patient and the operating room personnel (users) from the transfer of microorganisms, body fluids and particulate material.	The Surgical Gowns are single use, disposable, fluid repellent garments intended to be used as sterile by operating room personnel during surgical procedures to protect both the patient and the operating room personnel (users) from the transfer of microorganisms, body fluids and particulate material.	The Surgical Gowns are single use, disposable, fluid repellent garments intended to be used as sterile by operating room personnel during surgical procedures to protect both the patient and the operating room personnel (users) from the transfer of microorganisms, body fluids and particulate material.
MDR Classification Rule :			
Conformity Annexes :			
Measuring Function :			
Notified body medical devices/PPE :	Local Authority Americas	Local Authority United Arab Emirates	Local Authority Kuwait

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

ANSI/AAMI PB70, EN 1041, EN 556-1, EN 13795, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 10993-5, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, EN ISO 10993-10, ISO 14001, ISO 11137-1, ISO 11137-2

Removable label

Yes

GMDN Code (Global Medical Device Nomenclature)

35091 Operating room gown, single-use

UNSPSC

42131702 Surgical gowns

Commodity Code

4818500000 Wearing apparel mainly made of tissue/pulp paper

Find out more at www.molnlycke.com



CE Responsibility / Legal Manufacturer

Mölnlycke Health Care AB, Gamlestadsvägen 3c, Box 13080, SE-402 52 Goteborg Sweden

Basic UDIDI

73324300000000186KH

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.

The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2024 Mölnlycke Health Care AB. All rights reserved.

