

708131 Tube Holder

Product details

Size : 2.5cm x 30cm
Descriptive feature : Adhesive, Velcro
Sterile : Sterile

Images



Delivered items

708131-30

Country of origin: Czechia

Shelf life: 5 years

Sterilization method: Gamma

Packing information: First packaging layer is a peel-open sterile barrier, paper/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 code	WxLxH (mm)	Vol (dm3)	Weight gross/net (kg)
Consumer pack	1	7323190251261			
Dispenser box	50	7323190250707	138x227x186	5.8	
Transport box	200	7323190250691	291x461x198	26.6	2.6 / 1.8
Pallet	9600	7323190250684	800x1200x1736		149.8 / n/a

708131-20

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2021 Mölnlycke Health Care AB. All rights reserved.



Country of origin: Czechia

Shelf life: 5 years

Sterilization method: Gamma

Packing information: First packaging layer is a peel-open sterile barrier, paper/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 code	WxLxH (mm)	Vol (dm3)	Weight gross/net (kg)
Consumer pack	1	7323190177431			
Dispenser box	50	7323190168811	142x227x186		
Transport box	200	7323190168804	291x461x198	26.6	2.4 / 1.8
Pallet	9600	7323190168798	800x1200x1736		

Material

Animal tissues : No
Human blood derivatives : No
Natural rubber latex : No
Medicinal substances : No
Phthalates : No
Polyvinyl chloride : No

Product Composition Drape

Product Component	Composition
Main Material	Adhesive coated polypropylene (hook), polyester (loop)

Technical

Dimension

Dimension text	Dimension value
Outer dimension	2.5 cm x 30 cm

Disposal instructions

Used Mölnlycke products and sterile barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. Mölnlycke products do not contain any hazardous substances that can leach out if the products are land filled. Transport- and dispenser boxes are designed to fit existing recovery systems. The Mölnlycke packaging system complies with the EU Directive on packaging and packaging waste.

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2021 Mölnlycke Health Care AB. All rights reserved.



Storage instructions

Mölnlycke Health Care recommends that Mölnlycke products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterile barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the sterile barrier is opened and the product is used.

Classifications

Regulation type(s)	MDD Class IS	Locally Regulated	Unregulated
--------------------	--------------	-------------------	-------------

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Removable label

No

GMDN Code (Global Medical Device Nomenclature)

62137 Surgical tubing holder/clip, sterile

UNSPSC

42131701 Surgical drapes

Commodity Code

3926909790 Plastics and articles thereof

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2021 Mölnlycke Health Care AB. All rights reserved.

