

705820 Camera Drape

Allows sterile handling of equipment that cannot be sterilised

Product details

Product group name :	Equipment drapes
Size :	14cm x 250cm
Descriptive feature :	Perforated tip
Colour :	Transparent
Sterile :	Sterile
Brand :	BARRIER®
Equivalences :	815332

Product descriptions

Short description : Allows sterile handling of equipment that cannot be sterilised

Long description : BARRIER® equipment drapes provide a transparent and sterile draping solution for equipment in the operating room.

- Allows sterile handling of equipment that cannot be sterilised.
- Impermeable materials.
- Easy to drape.

Images



Delivered items

705820-21

Country of origin: Thailand

Shelf life: 5 years

Sterilization method: Irradiation

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, Amata Nakorn (Bang Pakong), Industrial

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.



Estate, 700/461 Moo 7 Bangna-Trad Road, Km 57 Tambol Donhuaroh, Amphur Muang Chonburi 20000, Thailand

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GTIN Code	UDI-DI	Width x Length x Height	Volume	Weight gross / net
Piece	1	7332430092472			344.312 cm3	49.583 g / -
Consumer pack	1	7323190292189	7323190292189			
Dispenser box	20	7323190289806	7323190289806	142x222x184 mm		
Transport box	120	7323190289790	7323190289790	298x683x203 mm	41.317 dm3	5.95 / 4.68 kg
Pallet	3840	7323190289783		800x1200x1768 mm	1697.28 dm3	215.4 / 149.76 kg

705820-15

Country of origin: Thailand

Shelf life: 5 years

Sterilization method: Irradiation

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, Amata Nakorn (Bang Pakong), Industrial Estate, 700/461 Moo 7 Bangna-Trad Road, Km 57 Tambol Donhuaroh, Amphur Muang Chonburi 20000, Thailand

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GTIN Code	UDI-DI	Width x Length x Height	Volume	Weight gross / net
Piece	1	7332430092472			334.088 cm3	52.75 g / -
Consumer pack	1	7323190041688	7323190041688			
Dispenser box	20	7323190041664	7323190041664			
Transport box	80	7323190041657	7323190041657	295x453x200 mm	26.727 dm3	4.22 / 3.16 kg

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.



Packing level	Quantity	GTIN Code	UDI-DI	Width x Length x Height	Volume	Weight gross / net
Pallet	3840	7323190041640		800x1200x1744 mm		

Material

Animal tissues :	No
Human blood derivatives :	No
Natural rubber latex :	No
Medicinal substances :	No
Phthalates :	No
Polyvinyl chloride :	No

Packaging Composition

Product Component	Composition
RET box	Corrugated board, 90-95% recycled content
Primary Packaging	Top web: Polyethylene film, Base web:PE/PA Laminate
TRP box	Corrugated board, 90-95% recycled content, FSC certified - C137267

Product Composition Drapes

Product Component	Composition
Drape material	Polyethylene film 50 µm, transparent
Adhesive material	Acrylate based

Product Performance Surgical Drapes, EN 13795-1:2019 High Performance

Characteristics	Test Method	Internal Test Method	Unit	Requirement Critical Product Area	Requirement Less Critical Product Area	Product Performance Critical Product Area	Product Performance Less Critical Product Area
Resistance to microbial penetration - Dry	ISO 22612	T-1004	CFU	Not required	≤300	N/A - no critical product area	0

Find out more at www.molnlycke.com



Characteristics	Test Method	Internal Test Method	Unit	Requirement Critical Product Area	Requirement Less Critical Product Area	Product Performance Critical Product Area	Product Performance Less Critical Product Area
Resistance to microbial penetration - Wet	ISO 22610	T-1005	BI	6.0	Not required	N/A - no critical product area	-
Cleanliness - Microbial	EN ISO 11737-1	T-303	CFU/100 cm ²	≤300	≤300	N/A - no critical product area	0
Particle release	EN ISO 9073-10	T-1006	Log ₁₀ (lint count)	≤4.0	≤4.0	N/A - no critical product area	N/A, 1-ply PE film
Resistance to liquid penetration	EN ISO 811	T-280	cm H ₂ O	≥100	≥10	N/A - no critical product area	>150
Bursting strength - Dry	EN ISO 13938-1	T-233 or T-1179	kPa	≥40	≥40	N/A - no critical product area	68
Bursting strength - Wet	EN ISO 13938-1	T-233 or T-1179	kPa	≥40	Not required	N/A - no critical product area	-
Tensile strength - Dry	EN 29073-3	T-229	N	≥20	≥20	N/A - no critical product area	34
Tensile strength - Wet	EN 29073-3	T-229	N	≥20	Not required	N/A - no critical product area	-

Product Performance Drapes, Additional Tests

Characteristics	Test Method	Internal Test Method	Unit	Product Performance
Absorption	ISO 9073-12	T-1158	g/dm ²	N/A
Flammability	16 CFR 1610.4	N/A	s	Class 1, >3.5s

Technical

Dimension

Dimension text	Dimension value
Outer dimension	14 cm x 250 cm

Find out more at www.molnlycke.com



Instructions

How product works :

The surgical equipment drapes are intended to be used as a protective equipment coverings, and to create a sterile field during the surgical procedure.

Classifications

Regulation type(s)	MDR Class IS
CE Certificate Number :	MDR 722028
Intended Purpose :	The surgical equipment drapes are intended to be used as a protective equipment covering, and to create a sterile field during the surgical procedure.
MDR Classification Rule :	1
Conformity Annexes :	IX
Measuring Function :	No
Notified body medical devices/PPE :	BSi 2797

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 13795, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Removable label

No

GMDN Code (Global Medical Device Nomenclature)

47783 Patient surgical drape, single-use

EMDN Code (European Medical Device Nomenclature)

T030199 COVERS, INSTRUMENTS AND EQUIPMENT - OTHER

Other specific product information

Equipment Category: Tube/Camera Cover

UNSPSC

42295415 Surgical equipment covers

Commodity Code

3926909790 Plastics and articles thereof

CE Responsibility / Legal Manufacturer

Mölnlycke Health Care AB, Entreprenorsstraket 21, SE-431 53 Mölndal Sweden

Basic UDIDI

73324300000000164K7

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.



Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.

The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.

