


# 694000 Dental Set

## Product details

**Sterile :** Sterile

**Brand :** BARRIER

## Images



## Products included in set

Product ID	Quantity	Product name	Size	Product description
	2	OP-Tape	9cm x 24,5cm	
815330	2	Camera Drape	8cm x 120cm	Plain tip
968004	2	OP-Towel	45cm x 75cm	2Ply, Absorbent
910176	1	Aperture Drape	150cm x 150cm, Aperture 10x14cm	Aperture
	1	Table Cover	75cm x 100cm	

## Delivered items

694000-07

**Country of origin:** Thailand

**Shelf life:** 5 years

**Sterilization method:** Irradiation

**Production Responsibility:** Mölnlycke Health Care (Thailand) Ltd, 160 Moo 17 Bangplee Industrial Estate, Bangna-Trad Road, Km 23 Tambol Bangsaothong, Amphur Bangsaothong, Samutprakarn 10570, Thailand

**Packing information:** First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a plastic liner/bag. Third layer is a corrugated board dispenser box. Fourth layer is a corrugated board transport box.

**Is suitable for Tray:** No

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7332551041472			
Consumer pack	1	7332551048020			
Dispenser box	14	7332551014469			
Transport box	28	7332551014452	397x580x329 mm	75.8 dm3	10.7 / 8.9 kg
Pallet	560	7332551014445	800x1200x1789 mm		

## Technical

### Disposal instructions

Used Mölnlycke products and sterile barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. Mölnlycke products do not contain any hazardous substances that can leach out if the products are land filled. Transport- and dispenser boxes are designed to fit existing recovery systems. The Mölnlycke packaging system complies with the EU Directive on packaging and packaging waste.

### Storage instructions

Mölnlycke Health Care recommends that Mölnlycke products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterile barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the sterile barrier is opened and the product is used.

### Classifications

Regulation type(s)	MDD Class IS	Locally Regulated
CE Certificate Number :	CE 01966	
Intended use MDD :	Surgical drapes, when sterilised, are intended to minimise the spread of micro-organisms, in order to reduce the risk for post operative wound infection.	
Intended Purpose :		Surgical drapes, when sterilised, are intended to minimise the spread of micro-organisms, in order to reduce the risk for post operative wound infection.
Conformity Annexes :	V	
Measuring Function :	No	
Notified body medical devices/PPE :	BSi 2797	Local Authority MEA

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



**Applied standards :** The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 13795, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

### **Removable label**

Yes

### **GMDN Code (Global Medical Device Nomenclature)**

33961 General/plastic surgical procedure kit, non-medicated, single-use

### **UNSPSC**

42131701 Surgical drapes

### **Commodity Code**

3005905000 Articles of Wadding, Gauze, Dressings, drapes singlepacked and Sets - textile

### **CE Responsibility**

Mölnlycke Health Care AB, Gamlestadsvägen 3c, Box 13080, SE-402 52 Goteborg Sweden

### **Basic UDIDI**

733243000000000182K9

**Find out more at [www.molnlycke.com](http://www.molnlycke.com)**

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.




Material information for the ingoing products

## 694000 Dental Set, 815330 Camera Drape

### Product details

Size : 8cm x 120cm  
 Descriptive feature : Plain tip  
 Colour : Transparent  
 Brand : BARRIER  
 Equivalences : 705740

### Images



### Technical

#### Dimension

Dimension text	Dimension value
Outer dimension	8 cm x 120 cm

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

# 694000 Dental Set, 968004 OP-Towel


## Product details

---

**Size :** 45cm x 75cm  
**Descriptive feature :** 2Ply, Absorbent  
**Brand :** BARRIER  
**Equivalences :** 800430

## Images

---



## Technical

---

### Dimension

Dimension text	Dimension value
Outer dimension	45 cm x 75 cm

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 13080, Gamlestadvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.  
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.  
 ©2022 Mölnlycke Health Care AB. All rights reserved.



## 694000 Dental Set, 910176 Aperture Drape


### Product details

**Size :** 150cm x 150cm, Aperture 10x14cm

**Descriptive feature :** Aperture

**Brand :** BARRIER

### Images



### Technical

#### Dimension

Dimension text	Dimension value
Outer dimension	150 cm x 150 cm
Aperture	10 cm x 14 cm
Aperture placement from top	50 cm

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.