



## PRODUCT SPECIFICATION



**Product:** CEDERROTH EYE&WOUND CLEANSING  
SPRAY 150 ML

**Reference code:** 726000



### 1. General description

A portable spray bottle filled with 150 mL physiological (0.9%) unbuffered sodium chloride solution, preservative free. The nozzle is suitable for use on both eyes and wounds. The actuator is protected by a plastic cap and a tamper seal. The nitrogen propellant is not in contact with the solution. The solution remains sterile until last drop.

### 2. Intended use

The Cederroth Eye&Wound Cleansing spray is intended to be used for rinsing and cleansing out dust and dirt from eyes and wounds.

### 3. Technical description

#### *Product dimensions*

#### Spray bottle (can)

Width (mm):	53
Height (mm)	162
Fill volume (ml):	150

#### *Material*

Liquid	Sterile 0.9% Sodium Chloride (NaCl). Preservative free.
Can	Aluminum.
Bag-on-valve	4 layer PP-OPA-ALU-PET
Actuator	White POM
Cover CAP	Transparent PP

### 4. Sterilization

Sterilization method: Irradiation

### 5. Composition

#### *Consumer package*

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>
726000	Cederroth Eye&Wound Cleansing Spray	1

#### *Transport package*

Number of consumer packages in transport package:	12
Dimensions (l x w x h) mm:	221 x 168 x 191
Weight (kg):	2,5



## PRODUCT SPECIFICATION



### 6. Labeling

#### *On product = consumer package*

Product name, volume, Reference code, CE+NB#, SterileR, Expiry date, LOT, EAN-code, temperature limitations, Instructions for use  
Cederroth address, telephone- and fax-number, e-mail and web-address

#### *On transport package*

Cederroth address, telephone- and faxnumber, web-address  
Product name, Reference code  
Lot-number, Expiry date, EAN-code, Limited quantity symbol  
Quantity, CE+NB#, SterileR

### 7. Storage conditions

Store in temperatures between 0°C and +40°C.

Shelf life: 3 years from production date

**Authorization:**

**QA&RA Mgr**

Digitally signed by QA&RA Mgr  
DN: cn=QA&RA Mgr, o=Wound Care Cederroth  
AB, ou, email=firstaid@cederroth.com, c=SE  
Reason: I am approving this document  
Date: 2012.01.30 11:45:43 +01'00'