




COVIDIEN

positive results for life™

Valleylab™ Tungsten Loop Electrodes and Stainless Steel Ball Electrodes For LLETZ procedures

Valleylab™ tungsten loop electrodes offer greater control and increased durability compared to stainless steel.

Assorted loop sizes are available for all your surgical needs. These fine wire loops minimize thermal damage and yield good histopathology specimens. Ball electrodes are also available for fulguration and desiccation.



Variety of offerings for best results



Convenient and safe


Single use Valleylab™ tungsten loop electrodes eliminate reprocessing and cross-contamination concerns for office and clinic staff.

ORDER INFORMATION

CATALOG NUMBER	DESCRIPTION	ORDER QUANTITY
E1559	Tungsten Loop, 10 mm x 10 mm	10/case
E1560	Tungsten Loop, 15 mm x 12 mm	10/case
E1561	Tungsten Loop, 20 mm x 12 mm	10/case
E1562	Tungsten Loop, 20 mm x 15 mm	10/case
E1565	Tungsten Square Loop, 5 mm x 5 mm	10/case
E1567	Tungsten Square Loop, 10 mm x 8 mm	10/case
E1563	Stainless Steel Ball, 3 mm diameter	10/case
E1564	Stainless Steel Ball, 5 mm diameter	10/case

Precise control and strength

Tungsten wire resists breakage even at higher power levels while providing better control against flexing and bending during use.



Compatible with

All Valleylab™ pencils and the Force TriVerse™ electrosurgical device


TECHNICAL SPECIFICATIONS

Shaft Diameter:

2.4 mm (0.0945 in.)

Shaft Length:

13 cm (5.1181 in.)



COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. ©2013 Covidien.

US130021 07/13

5920 LONGBOW DRIVE
BOULDER, CO
80301

303-530-2300 [T]
800-722-8772 [T]

WWW.COVIDIEN.COM