

Prolapse Ring

PER00 - PER09



Prolapse Ring

The flexible silicone Prolapse ring is commonly used for a mild first or second degree uterine prolapse or procidentia the Donut used for a third uterine prolapse as well as cystocele and rectocele and often used by post-menopausal women.

- The Prolapse ring are made from 100% medical grade silicone.
- The Prolapse rings are designed for single patient use.
- Fitting usually requires a trial of various sizes to determine the proper Prolapse ring size.
- The Prolapse ring Fitting Set is a valuable aid in selecting the correct Prolapse ring.
- A vaginal lubricant may be required for the patient's comfort before insertion of the Prolapse ring.
- The physician may suggest K-Y jelly®, or equivalent.

The physician should determine if the patient requires oestrogen therapy before prescribing a pessary. Within 24 to 48 hours, ensure the patient is not allergic to the pessary. Examine the vagina and ask the patient if there has been any discomfort, irritation, pressure, sensitivity, or unusual vaginal discharge. Also determine if there has been any improvement in her personal symptoms.


Guangzhou Fame Medical Co., Ltd. is a company specializing in the R&D, production, and sales of medical silicone products. Founded in 2007, Fame Medical has been dedicated to the R&D and production of medical silicone products for over ten years, possessing a profound understanding of silicone products and their technical processes. Fame Medical can provide comprehensive production process training and product processing services, including product development, mold design, and medical silicone accessories.

Fame Medical has a 100K-level clean workshop and a physical and chemical testing laboratory, and it has passed the ISO13485 quality management system certification and obtained CE ISO certification.

Material Silicone



PER00	Prolapse Ring 43mm (9613-0)
PER01	Prolapse Ring 51mm (9613-1)
PER02	Prolapse Ring 55mm (9613-2)
PER03	Prolapse Ring 62mm (9613-3)
PER04	Prolapse Ring 67mm (9613-4)
PER05	Prolapse Ring 74mm (9613-5)
PER06	Prolapse Ring 80mm (9613-6)
PER07	Prolapse Ring 87mm (9613-7)
PER08	Prolapse Ring 95mm (9613-8)
PER09	Prolapse Ring 99mm (9613-9)



ARTICLE NO.	DEPT. PACKAGING
PER00 - PER09	1 pc