




Fabian ventilatorer










Partnumber	Description
P/N 7200	fabianHFO Ventilator (without HFO module)
	<p>Neonatal/Pediatric ventilator with battery backup (conventional mode 2.5 hours HFO mode 1 hour). Automatic O2 sensor calibration. User interface: colour touchscreen Time cycled, pressure controlled Ventilation modes (TCPL): IPPV/IMV, SIPPV, SIMV, SIMV+PSV, PSV, CPAP, Manual 100 – 240 V, 50/60Hz</p> <p>Including standard accessories:</p> <ul style="list-style-type: none"> 1 pc. P/N 1031 Flow sensor for neonates, autoclavable/reusable 1 pc. P/N 1016 Flow sensor cable 1 pc. P/N 1116 Test lung 1 pc. P/N 7052 Membrane holder 1 pc. P/N 7007 Membrane for exhalation valve 1 pc. Powercord (specify version at order!) 1 pc. User manual

Partnumber5219	Description
7208	NIV (non-invasive ventilations with generators) Infant Flow/Inspire/MediJet nCPAP driver software. Modes: nCPAP and duoPAP Incl. 2 pcs. P/N 7067 Infant Flow connecting tube 60 cm.
7223	High Frequency Oscillation Module (upgrade kit) including Hard- and Software
7224	High Flow Oxygen Therapy HFT/HFNC. Flow: 1 – 20 LPM, FiO ₂ : 21 – 100 Vol. % O ₂
7225	Patient Data Management VueLink, AcuLink
7226	Data Trending, 5 days
7251	Volume Guarantee function VG
7252	Volume Limit function VL
7228	VG + VL package
7227	etCO ₂ Option, sidestream microstream technology. Oridion™
7210	etCO ₂ Option with Capnostat 5 mainstream sensor
7256	Lung recruitment for HFO ventilation



Neonatal/Paediatric ventilator, model fabian+ NCPAP Evolution	
Partnumber	Description
P/N 7250	fabian Evolution
	<p>Neonatal/Pediatric ventilator with battery backup 3.0 hours Automatic O2 sensor calibration. User interface: colour touchscreen Time cycled, pressure controlled Ventilation (TCPL) modes: IPPV/ IMV, SIPPV, SIMV, SIMV + PSV, PSV, CPAP, nCPAP and duoPAP, Manual 100-240 V, 50/60 Hz;</p> <p>Including standard accessories: 1 pc. P/N 1031 Flow Sensor for neonates, autoclavable/reusable 1 pc. P/N 1016 Flow Sensor cable 1 pc. P/N 1116 Test lung 1 pc. P/N 7052 Membrane holder 1 pc. P/N 7007 Membrane for exhalation valve 1 pc. P/N 7208 NIV (non-invasive ventilations with generators) InfantFlow/Inspire/Infant Flow LP/MediJet nCPAP driver software Modes: nCPAP and duoPAP Incl. 2 pcs. P/N 7067 Infant Flow connecting tube 60 cm 1 pc. Powercord (specify version with order!)</p>
Partnumber	Description
7252	Volume Limit function VL
7251	Volume Guarantee function VG
7228	VG + VL package
7224	Flow: 1 – 20 LPM, FiO2: 21 – 100 Vol. % O2
7225	Patient Data Management VueLink, AcuLink
7226	Data Trending, 5 days

nCPAP/duoPAP Driver, model fabian Therapy Evolution	
Partnumber	Description
P/N 7250.N	fabian Therapy evolution
	<p>NIV (non-invasive ventilations with generators) Infant Flow/Inspire/Infant Flow LP/MediJet</p> <p>nCPAP/duoPAP driver and O2 Therapy unit with battery backup 3.0 hours. Automatic O2 sensor calibration. User interface: colour touchscreen 100 – 240 V, 50/60Hz; nCPAP/duoPAP, High Flow Oxygen Therapy HFT/HFNC</p> <p>Including standard accessories: 2 pcs. P/N 7067 Infant Flow connecting tube, 60 cm 1 pc. Powercord (specify version at order!) 1 pc. User manual</p>
7225	Patient Data Management VueLink, AcuLink
7226	Data Trending, 5 days




Accessories and consumables	
Partnumber	Description
P/N 1031	
	Flow sensor, reusable, for fabian ventilators Autoclave at 134°C. Individually packed
P/N 7058.1	P/N 7058.1
	Flow sensor, single use, for fabian Individually packed, box of 10pcs. Clean but not sterile
P/N 7359	
	Pediatric Flow Sensor for fabian HFO Reusable Flow range from 0 up to 60 L/min Autoclavable at 134°C Individually packed
P/N 1016	
	Flow Sensor cable for fabian ventilators Except fabian Therapy
P/N 7007	
	Membrane for exhalation valve for fabian ventilators Autoclavable at 134°C
P/N 7052	
	Membrane holder for fabian nCPAP and fabian HFO ventilators. With 22 mm male / 15 mm female connection Made of Stainless steel Autoclavable at 134°C.
P/N 7360	
	Membrane holder for fabian nCPAP and fabian HFO ventilators. With 22 mm male / 15 mm female connection Made of special plastic for rainout reduction Autoclavable at 134°C.

<p>P/N 7041</p> 	<p>Transport Handle to fabian Therapy and +nCPAP</p>
<p>P/N 7209</p> 	<p>Interface adapter for fabian HFO. Connection between HFO port and inspiratory port. Hytrel, 30 cm long, with 15 mm cuff on HFO port and 22 mm cuff on inspiratory port with T-piece Reusable, soakable in disinfecting solution Autoclavable at 134°C</p>
<p>P/N 7209.e</p> 	<p>Interface adapter for fabian HFO. Connection between HFO port and inspiratory port, machine outlet. Silicone tube, 30 cm, with 15mm cuff on HFO port and 22 mm cuff on inspiratory port with T-piece Reusable Autoclavable at 134°C</p>
<p>P/N 7099</p> 	<p>Connector, right angled, 90° 22 mm O.D / 15 mm I.D. to 22 mm I.D. Autoclavable at 134°C</p>
<p>P/N 2858</p> 	<p>Adapter, right angled, with 15mm male to 22mm male connection to fit 7057.1 or 7067.1 on reusable humidifier chamber inlet. (to be used with flow generator circuits in combination with reusable water chamber only)</p>
<p>P/N 1116</p> 	<p>Test lung for neonatal ventilators, with 2 bellows, to simulate one lung ventilation and ET tube leakage Not autoclavable</p>
<p>P/N 7204</p> 	<p>Test lung silicone with 1 bellow Autoclavable at 134° C</p>




5-Wheel Trolley for fabian HFO - Standard Trolley 1

P/N 7243	
	<p>ACUTRONIC Cart, fabian HFO and monsoon III</p> <ul style="list-style-type: none"> • Height 85 cm • With 6 kg weight installed on base • 2 wheels with brake • Holder for cable and hoses • Bracket for humidifier included <p>Add mounting plate!</p>
P/N 7310	
	<p>Mounting plate for fabian HFO on ACUTRONIC Cart, 5-wheel & on Smart Trolley, 4-wheel</p> <ul style="list-style-type: none"> • With two DIN rails • Made of stainless steel




5-Wheel Trolleys for fabian + nCPAP evolution and fabian Therapy evolution – Standard Trolley 2

Partnumber	
P/N 7249	
	<p>ACUTRONIC Light Cart, for fabian evolution series: 5-star twin caster trolley</p> <ul style="list-style-type: none"> • Height 85 cm • 2 wheels with brake • Holder for cable and hoses • Bracket for humidifier included <p>Add mounting plate!</p>
P/N 7346	
	<p>Mounting plate for fabian +nCPAP evolution and fabian Therapy evolution ACUTRONIC Light Cart, 5-wheel & on Smart Trolley, 4-wheel</p> <ul style="list-style-type: none"> • With two DIN rails • Made of stainless steel
P/N 2947	
	<p>Korg för tillbehör för Acutronic cart and lightcart, storlek 300x200x180 mm</p>
P/N 2921	
	<p>Hållare för 2 gasflaskor upp till 5 l / st att monteras på rullstativ</p>






4-Wheel Trolley for fabian HFO, fabian +nCPAP evolution



Partnumber	Description
P/N 7063	
	<p>ACUTRONIC Smart Trolley to fabian HFO ventilator, fabian +nCPAP evolution, fabian Therapy evolution 4-casters trolley (with 2 brakes) Bracket for humidifier included Note: Add mounting plate! Note: Only to be offered/sold in connection with air compressor P/N 7065 & P/N 7066</p>
P/N 7346	
	<p>Mounting plate for fabian +nCPAP evolution and fabian Therapy evolution on ACUTRONIC Light Cart, 5-wheel & on Smart Trolley, 4-wheel</p> <ul style="list-style-type: none"> • With two DIN rails • Made of stainless steel
P/N 7310	
	<p>Mounting plate for fabian HFO on ACUTRONIC Cart, 5-wheel & on Smart Trolley, 4-wheel</p> <ul style="list-style-type: none"> • With two DIN rails • Made of stainless steel





Flexible Patient Arm, IV-pole, Rail Mounting Bracket



<p>Partnumber</p>	
<p>P/N 3590</p>	
	<p>Flexible arm for patient circuit for Fabian ventilators. One knob for fixation and adjustment</p>
<p>P/N 4188.N</p>	
	<p>IV-pole with support for humidifier. Adjustable length 1.2 m fits on standard rail bar</p>
<p>P/N 5909</p>	
	<p>Rail mounting bracket for humidifier (to fix humidifier on standard rail bar)</p> <p>(Optional item only – already included in ACUTRONIC Carts + Smart Trolley as standard)</p>






Air/O2 supply accessories

Partnumber	Description
P/N 5131	
	Air inlet water-trap with female NIST adapter on ventilator side and NIST adapter male for inlet.
P/N 1124	
	NIST Connector O ₂ with hose connection OD 7mm
P/N 1123	
	NIST connector AIR with hose connection OD 7mm
3030	
	High pressure AIR inlet hose, black/white colour, length 3 m, ISO standard with NIST connector, Swedish wall connector.
3050	
	High pressure Oxygen inlet hose, white colour, length 3 m, ISO standard with NIST connector, with Swedish wall adapter.

Medical Grade Air Compressors	
Partnumber	Description
P/N 7065	
	<p>Medical compressor ACUTRONIC "Silent", 230V/50Hz can be used on Smart Trolley Including pressure connecting hose between compressor and fabian ventilator. Compressor side DISS connection Ventilator side with NIST connection Technical details:</p> <ul style="list-style-type: none"> • 40 LPM at 3 bar • Air tank: 2 litres • Peak flow at 3 bar: 200 LPM • Air filtration: 5 micro • Dew point: 5°C below ambient temperature • Noise level: 49 dB [A] • Dimensions: 440 x 345 x 440 mm (WxDxH) • Weight: 30 kg
P/N 7248	
	<p>Medical compressor ACUTRONIC "Silent", 115V/60Hz can be used on Smart Trolley Including pressure connecting hose between compressor and fabian ventilator. Compressor side DISS connection Ventilator side with NIST connection Technical details: Same as P/N 7065</p>

F&P Humidifiers and related accessories	
P/N 5004	
	MR 850 Respiratory Humidifier base, 240 Volt/50Hz, including operators manual and powercord. Without emperature/flow probe and heater wires!
P/N 5013	
	Dual Temperature / Flow probe, right angle airway probe for use with circuit, length 1.45m, 900MR868
P/N 5023	
	Heater wire adapter for MR850 for dual heated single use breathing circuits, 900MR805 (clover leave type)
P/N 5029	
	Heater wire adapter for MR850 for single heated single use breathing circuits, 900MR806 (clover leave type)

P/N 5025	
	Heater wire adapter for MR850 for single heated reusable breathing circuits, 900MR858
P/N 5026	
	Heater wire adapter for MR850 for single heated single use breathing circuits, 900MR859. Old style connector, MOLEX
P/N 5048	
	Heater Wire - Inspiratory - Neonatal Low Flow - 24 W - 1.05m, used with MR850 and MR700 series. 900MR755
P/N 5045	
	Draw Wire - 1.5m
WM105107-1	
	neoPOD transportbefuktare


Water Chambers for Humidifiers	
Reusable chambers	
P/N 5912	
	Water chamber for humidifier, neonatal, low volume, reusable for PMH5000 and F&P models. If Infant Flow or Inspire single use circuits are used, need for adapter P/N 2858 (right angled 15 mm male to 22 mm male) Does not fit AIRcon humidifiers!
P/N 2858	
	Adapter, right angled with 15 mm male to 22 mm male connection to fit P/N 7057.1 or P/N 7067.1 on reusable humidifier chamber inlet. (To be used with flow generator circuits in combination with reusable water chamber only)
Single use chamber	
WM105113	
	Befuktningsskammare för användning med neoPOD transportbefuktare. neoPOD Lavabed kit (10/box).
AC914990	
	ventilatorslangar för transport med neoPOD-anlutning för Lavabed. Kan användas som dubbelslang eller enkelslang. (10/box)
P/N 5295	
	Humidification chamber, autofeed. Single use. Box of 20 pcs.
Disposable viral and bacterial filter	
P/N 5296	
	Bidirectional bacteria/virus filter for ventilator disposable box of 100





Reusable Patient Circuits for Fabian Ventilator




Partnumber	Description
P/N 7212	
	<p>Heated patient circuit for Fabian and Fabian HFO, reusable, complete set fits F&P humidifiers: including heater wire 900MR755 Includes: With Humidifier chamber low volume, Connecting tube ventilator to humidifier chamber Inspiratory limb with temperature probe adapter Extension tube for incubator use Patient y-piece with proximal pressure port Expiratory limb with water trap Proximal pressure tube</p>








Partnumber	Description
P/N 7213	<p>Heated patient circuit for fabian and fabian HFO, reusable, complete set fits F&P humidifiers: excluding heater wire! Includes: Without Humidifier chamber. Connecting tube ventilator to humidifier chamber Inspiratory limb with temperature probe adapter Extension tube for incubator use Patient y-piece with proximal pressure port Expiratory limb with water trap Proximal pressure tube Select appropriate inspiratory tube heater and water chamber, based on humidifier used with the system!</p>







	Single Use Patient Circuits for fabian Ventilators
Compatible with F&P and Similar Humidifiers	Description
P/N 5293	
	<p>Neonatal breathing circuit (>4 L/min), double limb system, 180cm, 10mm ID, single heated, expiratory water trap. Nitric Oxide kit. Extension line for incubators. Autofeed humidity chamber included. Single use. For use with fabian HFO, fabian +nCPAP evolution and other conventional ventilators Box of 10 pcs.</p>
P/N 5294	
	<p>Neonatal breathing circuit (>4 L/min), double limb system, 180cm, 10mm ID, single heated, expiratory water trap. Nitric Oxide kit. Extension line for incubators. Single use. For use with fabian HFO, fabian +nCPAP evolution and other conventional ventilators Box of 10 pcs</p>
P/N 5297	
	<p>Pediatric breathing Circuit dual limbs, inspiratory limb heated, expiratory water trap. Single use. For use with fabian HFO and fabian +nCPAP evolution and other ventilators. Box of 10 pcs</p>
P/N 5298	
	<p>HFO interface adaptor (T-Piece) with angled connector ready for inhaled Nitric Oxide. Single use. Box of 10 pcs. Note: only for fabian HFO with HFO module</p>



P/N 5298	
	<p>Transport breathing circuit, neonatal. Length 120 cm, with proximal pressure line 4x6 mm, smooth bore. Single use, Box of 10 pcs</p>
Heating & Dual Heating Cable for Wilamed A	
P/N 5306	NEW
	<p>Heating cable for Wilamed AIRcon humidifier To be used in combination with new single heated circuits Reusable</p>
P/N 5307	NEW
	<p>Dual heating cable for Wilamed AIRcon humidifier To be used in combination with dual heated limbs circuits Reusable</p>


	Transport Patient Circuits for fabian Ventilators
P/N 5299	
	<p>Transport breathing circuit, neonatal. Length 120 cm, with proximal pressure line 4x6 mm, smooth bore.</p> <p>Single use, Box of 10 pcs</p>
P/N7055	
	<p>Transport breathing circuit for fabian. 90 cm length, with proximal airway pressure tube, including patient Y-piece (Conventional ventilation only)</p> <p>Reusable</p> <p>Autoclavable at 134°C</p>
Single Use Patient Circuits with Dual Heated Limbs for fabian Ventilators	
5300	NEW
	<p>Breathing circuit, neonatal, dual limbs, heated, with humidity chamber</p> <p>Single use Box of 10 pcs.</p>
5301	
	<p>Breathing circuit, for adults, dual limbs, heated, with humidity chamber</p> <p>Single Use Box of 10 pcs.</p>
5302	
	<p>Breathing circuit, pediatric, dual limbs, heated, with humidity chamber</p> <p>Single use Box of 10 pcs.</p>





INSPIRE nCPAP Circuits for NIV Ventilation





P/N 5291	
	<p>NIV Patient circuit, heated with Infant Flow connecting tube (P/N 7067.1) included. Nitric Oxide Kit. Extension line for incubator. Compatible with Medijet generator by Medin. Compatible with nasal and RAM cannulas. Autofeed humidity chamber, included Single use. Box of 10 pcs. Note: needs Inspire nCPAP generator (P/N 7186)</p>
P/N 5292	
	<p>NIV Patient circuit, heated with Infant Flow connecting tube (P/N 7067.1) included. Nitric Oxide Kit. Extension line for incubator. Compatible with Medijet generator by Medin. Compatible with nasal and RAM cannulas. Single use. Box of 10 pcs. Note: needs Inspire nCPAP generator (P/N 7186)</p>
P/N 7186	
	<p>INSPIRE nCPAP generator, Includes generator with Coanda effect and Nasal prongs (S,M,L) Box of 10 pcs.</p>
P/N 7197	new
	<p>INSPIRE nCPAP Generator + Inspire/InfantFlow Patient circuit, heated Set of P/N 7186 and P/N 5292</p>

Connecting Tube for nCPAP Generators Used on fabian	
P/N 7067.1	
	InfantFlow connecting tube. fabian to humidifier chamber inlet. Length 60 cm, box of 10 pcs. If reusable chamber is used, #2858 is needed to fit chamber inlet!
Nasal Cannula for O₂ Therapy	
5303	NEW
	Nasal cannula for O ₂ Therapy small (S) Preemies Single use Box of 50 pcs.
5304	NEW
	Nasal cannula for O ₂ Therapy medium (M) Neonates Single use Box of 50 pcs.
5305	NEW
	Nasal cannula for O ₂ Therapy large (L) Infant Single use Box of 50 pcs.

Inspire nCPAP prongs and masks	
Partnumber	Description
P/N 7167	
	INSPIRE nCPAP Nasal Prong large, box of 10 pcs.
P/N 7168	INSPIRE nCPAP Nasal Prong medium, box of 10 pcs.
P/N 7169	INSPIRE nCPAP Nasal Prong small, box of 10 pcs.
P/N 7189	
	Inspire nCPAP Mask, size XL Box of 10 pcs.
P/N 7188	Inspire nCPAP Mask, size L Box of 10 pcs.
P/N 7182	Inspire nCPAP Mask, size M Box of 10 pcs.
P/N 7183	Inspire nCPAP Mask, size S Box of 10 pcs.

INSPIRE nCPAP bonnets	
Partnumber	Description
P/N 7170	
	INSPIRE nCPAP Bonnet, Size 9, navy, 38 - 40 cm, box of 10 pcs.
P/N 7171	INSPIRE nCPAP Bonnet, Size 8, dark green, 36 - 38 cm, box of 10 pcs.
P/N 7172	INSPIRE nCPAP Bonnet, Size 7, orange, 34 - 36 cm, box of 10 pcs.
P/N 7173	INSPIRE nCPAP Bonnet, Size 6, burgundy, 32 - 34 cm, box of 10 pcs.
P/N 7174	INSPIRE nCPAP Bonnet, Size 5, green, 30 - 32 cm, box of 10 pcs.
P/N 7175	INSPIRE nCPAP Bonnet, Size 4, gold, 28 - 30 cm, box of 10 pcs.
P/N 7176	INSPIRE nCPAP Bonnet, Size 3, blue, 26 - 28 cm, box of 10 pcs.
P/N 7177	INSPIRE nCPAP Bonnet, Size 2, yellow, 24 - 26 cm, box of 10 pcs.
P/N 7178	INSPIRE nCPAP Bonnet, Size 1, brown, 22 - 24 cm, box of 10 pcs.
P/N 7179	INSPIRE nCPAP Bonnet, Size 0, pink, 20 - 22 cm, box of 10 pcs.
P/N 7180	INSPIRE nCPAP Bonnet, Size 00, grey, 18 - 20 cm, box of 10 pcs.
P/N 7181	INSPIRE nCPAP Bonnet, Size 000, white, 16 - 18 cm, box of 10 pcs.

INSPIRE nCPAP transport circuits	
Partnumber	Description
P/N 7190	
	Transport single use InfantFlow/Inspire connecting tube for fabian. Hose between the unit and flow generator, box of 10 pcs
P/N 7191	
	Transport single use InfantFlow/Inspire set for fabian (includes: 7190 and 7186) box of 10 pcs.
Aerob Pro	
Partnumber	Description
7120	
	MiniMe 2, Pediatric Nasal Mask, Size S/M
7121	
	MiniMe 2, Pediatric Nasal Mask, Size M/L

Spareparts for fabian ventilators	
Partnumber	Description
P/N 7003	
	Oxygen sensor for fabian and fabian HFO, dual cathode, approximately 12 month lifetime fabian HFO fabian + nCPAP evolution
P/N 7714	
	Sinter filter (candle) for monsoon III V2.0, fabian HFO and fabian evolution from SN 20-00208 Fits NIST connector item 7970 and 7971
P/N 7265	
	Service Set, 1 year, for fabian HFO: Includes: <ul style="list-style-type: none"> • Oxygen sensor • Air Intake filter • Reusable flow sensor
P/N 7264	
	Service Set, 4 years, for fabian HFO Includes: <ul style="list-style-type: none"> • Annual preventive maintenance kit • Pressure regulators air/oxygen • Lithium battery • Calibration valves air/oxygen