

Prolapse Donut

PED00 - PED07



Prolapse Donut

The flexible silicone Prolapse ring is commonly used for a mild first or second degree uterine prolapse or procidentia the Donut used for a third uterine prolapse as well as cystocele and rectocele and often used by post-menopausal women.

- The Prolapse ring are made from 100% medical grade silicone.
- The Prolapse rings are designed for single patient use.
- Fitting usually requires a trial of various sizes to determine the proper Prolapse ring size.
- The Prolapse ring Fitting Set is a valuable aid in selecting the correct Prolapse ring.
- A vaginal lubricant may be required for the patient's comfort before insertion of the Prolapse ring.
- The physician may suggest K-Y jelly®, or equivalent.

The physician should determine if the patient requires oestrogen therapy before prescribing a pessary. Within 24 to 48 hours, ensure the patient is not allergic to the pessary. Examine the vagina and ask the patient if there has been any discomfort, irritation, pressure, sensitivity, or unusual vaginal discharge. Also determine if there has been any improvement in her personal symptoms.


Guangzhou Fame Medical Co., Ltd. is a company specializing in the R&D, production, and sales of medical silicone products. Founded in 2007, Fame Medical has been dedicated to the R&D and production of medical silicone products for over ten years, possessing a profound understanding of silicone products and their technical processes. Fame Medical can provide comprehensive production process training and product processing services, including product development, mold design, and medical silicone accessories.

Fame Medical has a 100K-level clean workshop and a physical and chemical testing laboratory, and it has passed the ISO13485 quality management system certification and obtained CE ISO certification.

Material Silicone



PED00	Prolapse Donut 50mm (9615-1)
PED01	Prolapse Donut 56mm (9615-2)
PED02	Prolapse Donut 64mm (9615-3)
PED03	Prolapse Donut 69mm (9615-4)
PED04	Prolapse Donut 74mm (9615-5)
PED05	Prolapse Donut 81mm (9615-6)
PED06	Prolapse Donut 87mm (9615-7)
PED07	Prolapse Donut 93mm (9615-8)



ARTICLE NO.	DEPT. PACKAGING
PED00-PED07	1 pc