


VWR® INDICATOR PAPERS/STRIPS

Quick and easy to use

Accurate, non bleeding

Extensive choice of pH ranges/scales

Peracetic acid and peroxide for disinfectant checks



Get answers in an instant!

A choice of rapid testing, pH measuring papers/strips ideal for applications in the field, laboratory, production and teaching. Providing a level of convenience surpassing even traditional stick meters which can be used by non scientific personnel.

One specialist area is verifying disinfectant solutions for cleaning – balancing efficiency with leftover residue can be tricky but with VWR peracetic acid or peroxide strips, cleaners can take control of the process.

More parameters are being introduced all the time so please go to vwr.com or contact your local VWR sales office if you cannot see what you'd like.



Rota pH test papers

- Easy to use: Simply match the colour to the scale
- Robust: Supplied on plastic reels for protection and stability
- Accurate: Scale and paper colour distinct match

Supplied in 5 m lengths, 10 mm wide with full instructions

pH range	Graduation	Cat. No.
1 - 11	1/2/3/4/5/6/7/8/9/10/11	85409.600
1 - 12	1/2/3/4/5/6/7/8/9/10/11/12	85403.600
1 - 14	1/2/3/5/6/7/8/9/10/12/14	85402.600
1 - 4,3	1,0/1,3/1,6/1,9/2,2/2,5/2,8/3,1/3,4/4,0/4,3	85404.600
3,5 - 8	3,5/3,8/4,1/4,4/4,7/5,0/5,3/5,6/5,9/6,2/6,5/6,8	85405.600
5 - 8	5,0/5,3/5,6/5,9/6,2/6,5/6,8/7,1/7,4/7,7/8,0	85406.600
7 - 10	7,0/7,3/7,6/7,9/8,2/8,5/8,8/9,1/9,4/9,7/10,0	85407.600
9,5 - 14	9,5/10,0/10,5/11,0/11,5/12,0/12,5/13,0/13,5/14,0	85408.600



pH indicator strips - DOSATEST

- Non bleeding, colour-fixed reactive pads - pH indicator is chemically bound to the cellulose fibre preventing bleeding and protects sample against contamination
- Long handle makes strips easy to work with for a variety of containers and prevents fingers from getting into the sample
- Precise and rapid pH determination with four different colour blocks for each pH value to match against a bright distinct colour chart
- Measurements even in weakly buffered or strongly alkaline solutions

Test strips 6x85 mm – supplied in boxes of 100

pH range	Graduation	Cat. No.
0 - 14	0/1/2/3/4/5/6/7/8/9/10/11/12/13/14	85410.601
0 - 6	0/0,5/1/1,5/2/2,5/3/3,5/4/4,5/5/5,5/6,0	85411.601
0,3 - 2,3	0,3/0,7/1,0/1,3/1,6/1,9/2,3	85421.601
1,7 - 3,8	1,7/2,0/2,3/2,6/2,9/3,2/3,5/3,8	85422.601
2 - 9	2,0/2,5/3,0/3,5/4,0/4,5/5,0/5,5/6,0/6,5/7,0/7,5/8,0/8,5/9,0	85412.601
3,6 - 6,1	3,6/4,1/4,4/4,7/5,0/5,3/5,6/6,1	85416.601
4 - 7	4,0/4,4/4,7/5,0/5,3/5,5/5,8/6,1/6,5/7,0	85423.601
4,5 - 10	4,5/5,0/5,5/6,0/6,5/7,0/7,5/8,0/8,5/9,0/9,5/10,0	85413.601
5,1 - 7,2	5,1/5,4/5,7/6,0/6,3/6,6/6,9/7,2	85417.601
6 - 10	6,0/6,4/6,7/7,0/7,3/7,6/7,9/8,2/8,4/8,6/8,8/9,1/9,5/10,0	85414.601
6,0 - 7,7	6,0/6,4/6,7/7,0/7,3/7,7	85418.601
7 - 14	7,0/7,5/8,0/8,5/9,0/9,5/10,0/10,5/11,0/11,5/12,0/12,5/13,0/13,5/14,0	85415.601
7,5 - 9,5	7,5/7,9/8,2/8,4/8,6/8,8/9,1/9,5	85419.601
7,9 - 9,8	7,9/8,3/8,6/8,9/9,1/9,4/9,8	85420.601



DOSATEST scale on the strip

Offering all the benefits of the DOSATEST strips above, the DOSATEST scale on the strip is ideal for coloured solutions as both scale and test area are equally impacted on by the sample. This means that there is no interference and ensures accurate results.

Handle protected by invisible barrier to prevent capillary action of sample up the strip.

Test strips 6x85 mm, supplied in packs of 200



pH scale range	Graduation	Cat. No.
1 - 12	1/2/3/4/5/6/7/8/9/10/11/12	85424.602
0 - 1,8	0/0,3/0,6/0,8/1,0/1,2/1,5/1,8	85425.602
1,8 - 3,8	1,8/2,1/2,4/2,7/3,0/3,2/3,5/3,8	85426.602
3,8 - 5,5	3,8/4,0/4,2/4,4/4,6/4,9/5,2/5,5	85427.602
5,2 - 6,8	5,2/5,5/5,7/5,9/6,1/6,3/6,5/6,8	85428.602
6,0 - 8,1	6,0/6,3/6,6/6,9/7,2/7,5/7,8/8,1	85429.602
8,0 - 9,7	8,0/8,2/8,4/8,6/8,8/9,1/9,4/9,7	85430.602
9,5 - 12,0	9,5/10,0/10,5/11,0/11,5/12,0	85431.602
12,0 - 14,0	12,0/12,5/13,0/13,5/14,0	85432.602

Disinfectant test strips

These easy to use strips enable you to make sure that your disinfectant solutions are at effective concentrations while leaving minimal undesirable residues on plant after cleaning.


- Rapid and convenient results in seconds without anything else needed
- Precise: Colours adjusted and checked using standard solutions
- Pictograms with clear instructions
- Low cost per test



Test strip	Max. mg/l	Graduations (mg/l)	Cat. No.
Peroxide	25	0/0,5/2/5/10/25	85433.601
	100	0/1/3/10/30/100	85434.601
	1000	0/50/150/300/500/800/1000	85435.601
Peracetic acid	50	0/5/10/20/30/50	85436.601
	500	0/50/100/200/300/400/500	85437.601
	2000	0/500/1000/1500/2000	85438.601

Have you got your copy of the VWR Chemicals catalogue yet?

Go to vwr.com or scan this QR code



Austria

VWR International GmbH
Graumannsgasse 7
1150 Wien
Tel.: +43 1 97 002 0
Email: info.at@vwr.com

Belgium

VWR International bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
3001 Leuven
Tel.: +32 (0) 16 385 011
Email: vwr.be@vwr.com

Czech Republic

VWR International s. r. o.
Veetee Business Park
Pražská 442
CZ - 281 67 Stříbrná Skalice
Tel.: +420 321 570 321
Email: info.cz@vwr.com

Denmark

VWR International A/S
Tobaksvej 21
2860 Søborg
Tel.: +45 43 86 87 88
Email: info.dk@vwr.com

Finland

VWR International Oy
Valimotie 9
00380 Helsinki
Tel.: +358 (0) 9 80 45 51
Email: info.fi@vwr.com

France

VWR International S.A.S.
Le Périgares – Bâtiment B
201, rue Carnot
94126 Fontenay-sous-Bois cedex
Tel.: 0 825 02 30 30* (national)
Tel.: +33 (0) 1 45 14 85 00 (international)
Email: info.fr@vwr.com
* 0,18 € TTC/min

Germany

VWR International GmbH
Hilpertstraße 20a
D - 64295 Darmstadt
Tel.: 0800 702 00 07* (national)
Tel.: +49 (0) 6151 3972 0 (international)
Email: info.de@vwr.com
*Freecall

Hungary

VWR International Kft.
Simon László u. 4.
4034 Debrecen
Tel.: +36 52 521130
Email: info.hu@vwr.com

Ireland / Northern Ireland

VWR International Ltd / VWR International
(Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 (0) 1 88 22 222
Email: sales.ie@vwr.com

Italy

VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02 3320311
Email: info.it@vwr.com

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: +31 (0) 20 4808 400
Email: info.nl@vwr.com

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: +47 22 90 00 00
Email: info.no@vwr.com

Poland

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: +48 58 32 38 200
Email: info.pl@vwr.com

Portugal

VWR International - Material de
Laboratório, Lda
Centro Empresarial de Alfragide
Rua da Indústria, nº 6
2610-088 Alfragide
Tel.: +351 21 3600 770
Email: info.pt@vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: +34 902 222 897
Email: info.es@vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: +46 (0) 8 621 34 00
Email: kundservice.se@vwr.com

Switzerland

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: +41 (0) 44 745 13 13
Email: info.ch@vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: +44 (0) 800 22 33 44
Email: uksales@vwr.com

China

VWR International China Co., Ltd.
Shanghai Branch
Room 256, No. 3058 Pusan Road
Pudong New District
Shanghai 200123
Tel.: +86 21 5898 6888
Email: info_china@vwr.com

India

VWR Lab Products Private Limited
No.139, BDA Industrial Suburb,
6th Main, Tumkur Road, Peenya Post,
Bangalore, India – 560058
Tel.: +91 80 28078400
Email: vwr_india@vwr.com

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel: +65 6505 0760
Email: sales.sg@vwr.com

GO TO VWR.COM FOR THE LATEST
NEWS, SPECIAL OFFERS AND
DETAILS FROM YOUR LOCAL
VWR SUPPORT TEAM