

708840 Adhesive Aperture Drape

Impermeable and absorbing drapes with a variety of apertures

Product details

Product group name :	Adhesive aperture drape
Size :	120cm x 150cm, Aperture 5x7cm
Descriptive feature :	Adhesive, Centered ellipse aperture
Sterile :	Sterile
Brand :	BARRIER®
Equivalences :	968840

Product descriptions

Short description : Impermeable and absorbing drapes with a variety of apertures

Long description : BARRIER® adhesive aperture drapes are available in a range of sizes and aperture shapes designed for different surgical needs. You can choose centred or off-centre aperture positions. We also offer variable aperture drapes which can be adjusted to different needs.

- Variety of adhesive aperture shapes and sizes to suit different procedures.
- Absorbent top layer
- Impermeable materials help prevent microbial transfer.
- Folded and packaged for easy, aseptic application by one person.
- Can be included in Mölnlycke® ProcedurePak® trays to support operating room efficiency.

When to use BARRIER aperture drape:

You can use aperture drapes for a variety of procedures such as ophthalmic surgery and minor surgical interventions.

Images



Delivered items

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.



708840-09

Country of origin: Thailand**Shelf life:** 5 years**Sterilization method:** Irradiation**Production Responsibility:** Mölnlycke Health Care (Thailand) Ltd, 160 Moo 17 Bangplee Industrial Estate, Bangna-Trad Road, Km 23 Tambol Bangsaothong, Amphur Bangsaothong, Samutprakarn 10570, Thailand**Packing information:** First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.**Is suitable for Tray:** No

Packing level	Quantity	GTIN Code	UDI-DI	Width x Length x Height	Volume	Weight gross / net
Piece	1	7332551041915			1278.864 cm3	157.5 g / -
Consumer pack	1	7323190250318	7323190250318			
Dispenser box	16	7323190250189	7323190250189	206x321x267 mm		
Transport box	32	7323190250172	7323190250172	332x428x288 mm	40.924 dm3	5.04 / 3.82 kg
Pallet	1152	7323190250165		800x1200x1872 mm	1797.12 dm3	206.44 / 137.52 kg

Material

Animal tissues :	No
Human blood derivatives :	No
Natural rubber latex :	No
Medicinal substances :	No
Phthalates :	No
Polyvinyl chloride :	No

Packaging Composition

Product Component	Composition
Primary Packaging	Top web: Polyethylene film, Base web:PE/PA Laminate
TRP box	Corrugated board, 90-95% recycled content, FSC certified - C137267
RET box	Corrugated board, 90-95% recycled content

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.



Product Composition Drapes

Product Component	Composition
Drape material	Viscose nonwoven 23 g/m ²
Drape material	Polyethylene film 40 µm
Adhesive material	Synthetic rubber based

Product Performance Surgical Drapes, EN 13795-1:2019 High Performance

Characteristics	Test Method	Internal Test Method	Unit	Requirement Critical Product Area	Requirement Less Critical Product Area	Product Performance Critical Product Area	Product Performance Less Critical Product Area
Resistance to microbial penetration - Dry	ISO 22612	T-1004	CFU	Not required	≤300	-	0
Resistance to microbial penetration - Wet	ISO 22610	T-1005	BI	6.0	Not required	6.0	-
Cleanliness - Microbial	EN ISO 11737-1	T-303	CFU/100 cm ²	≤300	≤300	0	0
Particle release	EN ISO 9073-10	T-1006	Log ₁₀ (lint count)	≤4.0	≤4.0	2.5	2.5
Resistance to liquid penetration	EN ISO 811	T-280	cm H ₂ O	≥100	≥10	127	127
Bursting strength - Dry	EN ISO 13938-1	T-233 or T-1179	kPa	≥40	≥40	67	67
Bursting strength - Wet	EN ISO 13938-1	T-233 or T-1179	kPa	≥40	Not required	45	-
Tensile strength - Dry	EN 29073-3	T-229	N	≥20	≥20	25	25
Tensile strength - Wet	EN 29073-3	T-229	N	≥20	Not required	22	-

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.



Product Performance Drapes, Additional Tests

Characteristics	Test Method	Internal Test Method	Unit	Product Performance
Absorption	ISO 9073-12	T-1158	g/dm ²	1,10
Flammability	16 CFR 1610.4	N/A	s	Class I, >3.5s

Technical

Dimension

Dimension text	Dimension value
Outer dimension	120 cm x 150 cm
Aperture	5 cm x 7 cm

Instructions

How product works :

The surgical patient drapes are intended to be used as protective patient coverings to create a sterile field for a surgical procedure.

Classifications

Regulation type(s)	MDR Class IS
Intended Purpose :	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.
Notified body medical devices/PPE :	BSi 2797

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 13795, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Removable label

Yes

GMDN Code (Global Medical Device Nomenclature)

47783 Patient surgical drape, single-use

EMDN Code (European Medical Device Nomenclature)

T020102 SPECIALISTIC SURGERY DRAPES

UNSPSC

42131701 Surgical drapes

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2025 Mölnlycke Health Care AB. All rights reserved.



Commodity Code

3005100000 Wadding, Gauze, Dressings, drapes singlepacked - adhesive articles, Sets mainly consisting thereof

CE Responsibility / Legal Manufacturer

Mölnlycke Health Care AB, Entreprenorsstraket 21, SE-431 53 Mölndal Sweden

Basic UDIDI

73324300000000160JX

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2025 Mölnlycke Health Care AB. All rights reserved.

