

# 888222 TUR Set

Procedure-specific urology drapes with effective solutions for fluid management

## Product details

**Product group name :** TUR drapes and sets


**Sterile :** Sterile

**Brand :** BARRIER®

## Product descriptions

**Short description :** Procedure-specific urology drapes with effective solutions for fluid management

## Images



## Products included in set

Product ID	Quantity	Product name	Size
	1	Cellulose Towel	18cm x 25cm
904762	1	Tube Holder	9 x 18cm
965520	1	TUR Drape	175cm / 270cm x 180cm, Aperture 8cm, 5cm
90190	1	Table Cover	150cm x 190cm

## Delivered items

888222-78

**Country of origin:** Thailand

**Shelf life:** 5 years

**Sterilization method:** Irradiation

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.

**Production Responsibility:** Mölnlycke Health Care (Thailand) Ltd, 160 Moo 17 Bangplee Industrial Estate, Bangna-Trad Road, Km 23 Tambol Bangsaothong, Amphur Bangsaothong, Samutprakarn 10570, Thailand

**Packing information:** First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a plastic liner/bag. Third layer is a corrugated board dispenser box. Fourth layer is a corrugated board transport box.

**Is suitable for Tray:** No

Packing level	Quantity	GTIN Code	UDI-DI	Width x Length x Height	Volume	Weight gross / net
Piece	1	7332551041694				892.857 g / -
Consumer pack	1	7323190287789	7323190287789			
Dispenser box	7	7323190287772	7323190287772	291x356x421 mm		
Transport box	14	07333350294700	07333350294700	371x593x440 mm	96.801 dm3	12.5 / 9.68 kg
Pallet	224	7313661376371		800x1200x1904 mm	1827.84 dm3	225 kg / -

888222-22

**Country of origin:** Belgium

**Shelf life:** 5 years

**Sterilization method:** Irradiation

**Production Responsibility:** Mölnlycke Health Care S.A., Chaussée Romaine 176, 4300 Waremme, Belgium

**Packing information:** First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

**Is suitable for Tray:** No

Packing level	Quantity	GTIN Code	UDI-DI	Width x Length x Height	Volume	Weight gross / net
Piece	1	7332551041694				8933.047 cm3 879.429 g / -
Consumer pack	1	07333350301491	07333350301491			
Dispenser box	7	07333350711016	07333350711016	366x290x448 mm		
Transport box	14	07333350239367	07333350239367	391x596x460 mm	107.197 dm3	12.312 / 9.781 kg

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.  
©2025 Mölnlycke Health Care AB. All rights reserved.



Packing level	Quantity	GTIN Code	UDI-DI	Width x Length x Height	Volume	Weight gross / net
Pallet	168	7313661257854		800x1200x1532 mm	14.707 dm3	172.744 kg / -

## Material

<b>Animal tissues :</b>	No
<b>Human blood derivatives :</b>	No
<b>Natural rubber latex :</b>	No
<b>Medicinal substances :</b>	No
<b>Phthalates :</b>	No
<b>Polyvinyl chloride :</b>	No

## Technical

## Instructions

## How product works :

The surgical patient drapes are intended to be used as protective patient coverings to create a sterile field for a surgical procedure.

## Classifications

Regulation type(s)	MDR Class IS	MDR Class IS
CE Certificate Number :	MDR 722028	MDR 722028
Intended Purpose :	The surgical drape sets consists of both surgical drapes and supplementary products intended be used to create a sterile field for a surgical procedure.	The surgical drape sets consists of both surgical drapes and supplementary products intended be used to create a sterile field for a surgical procedure.
MDR Classification Rule :	5	1
Conformity Annexes :	IX	IX
Measuring Function :	No	No
Notified body medical devices/PPE :	BSi 2797	BSi 2797

**Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.**

EN 1041, EN 556-1, EN 13795, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.



## Removable label

Yes

## GMDN Code (Global Medical Device Nomenclature)

33961 General/plastic surgical procedure kit, non-medicated, single-use

## EMDN Code (European Medical Device Nomenclature)

T0202 SURGICAL PROCEDURAL KITS (EXCLUDING SURGICAL INSTRUMENT KITS)

## UNSPSC

42131701 Surgical drapes

## Commodity Code

3005905000 Articles of Wadding, Gauze, Dressings, drapes singlepacked and Sets - textile

## CE Responsibility / Legal Manufacturer

Mölnlycke Health Care AB, Entreprenorsstraket 21, SE-431 53 Mölndal Sweden

## Basic UDIDI

733243000000000182K9

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.  
©2025 Mölnlycke Health Care AB. All rights reserved.



Material information for the ingoing products

## 888222 TUR Set, 904762 Tube Holder

Secure fixation of tubes and cables to surgical drapes

### Product details

---

**Product group name :** OP-Tape & Tube holders


**Size :** 9 x 18cm

**Descriptive feature :** Adhesive

**Brand :** BARRIER®

### Images

---



### Material

---

**Animal tissues :** No

**Human blood derivatives :** No

**Natural rubber latex :** No

**Medicinal substances :** No

**Phthalates :** No

**Polyvinyl chloride :** No

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.

## 888222 TUR Set, 965520 TUR Drape

Procedure-specific urology drapes with effective solutions for fluid management

### Product details

---

<b>Product group name :</b>	TUR drapes and sets
<b>Size :</b>	175cm / 270cm x 180cm, Aperture 8cm, 5cm
<b>Descriptive feature :</b>	Aperture, Finger cot, Leggings, Pouch
<b>Brand :</b>	BARRIER®

### Images

---



### Material

---

<b>Animal tissues :</b>	No
<b>Human blood derivatives :</b>	No
<b>Natural rubber latex :</b>	No
<b>Medicinal substances :</b>	No
<b>Phthalates :</b>	No
<b>Polyvinyl chloride :</b>	No

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.

## 888222 TUR Set, 90190 Table Cover

Impermeable and absorbing covers that conforms to the shape of the table

### Product details

---

**Product group name :** Table covers  
**Size :** 150cm x 190cm  
**Descriptive feature :** Absorbent  
**Brand :** BARRIER®  
**Equivalences :** 80190

### Images

---



### Material

---

**Animal tissues :** No  
**Human blood derivatives :** No  
**Natural rubber latex :** No  
**Medicinal substances :** No  
**Phthalates :** No  
**Polyvinyl chloride :** No

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 6, SE-431 21 Mölndal, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2025 Mölnlycke Health Care AB. All rights reserved.