




# MD300M

## HANDHÅLLEN PULSOXIMETER

ChoiceM Med Handhållen Pulsoximeter MD300M kan ersätta patientmonitorn och användas som patientövervakning på sjukhus och kliniker. Denna serien pulsoximetrar erbjuder praktiska kostnadseffektiva lösningar för mindre kliniker och vårdavdelningar.



- ♥ Högupplösande TFT skärm visar SpO<sub>2</sub>, PR, Perfusions Index (PI)\*, pulse stapel, och kurvform
- ♥ Justerbara akustiska och visuella larm; svagt-batteri indikator
- ♥ Justerbar bakgrundsbelysning och ljud
- ♥ 127 ID setup; 72-timmar datalagring och återuppspelning
- ♥ MedView programvara för data analys
- ♥ 3 AA alkaliska batterier eller nätadapter (tillval)
- ♥ Lämplig för vuxen och pediatrikt bruk
- ♥ Flerspråkig (Meny): Engelska, Franska, Tyska, Spanska, Italienska, Japanska, Ryska, Kinesiska
- ♥ Certifikat: FDA, CE, CFDA



### TEKNISKA SPECIFIKATIONER

Application Range			
Adult and pediatric			
SpO <sub>2</sub>			
Measurement Range:	70%~100%	Resolution:	1%
Measurement Accuracy:	70%~100%±2%; <70% unspecified		
Pulse Rate			
Measurement Range:	30bpm~250bpm	Resolution:	1bpm
Measurement Accuracy:	±2 bpm or 2% (The larger)		
PI*			
Measurement Range:	0.1%~20%		
Display			
Type:	TFT	Parameters:	SpO <sub>2</sub> %, PR, PI*
Environment Requirements			
Operation Temperature:	5°C~40°C	Operation Humidity:	≤ 80%, no condensation
Storage Temperature:	-20°C ~ 70°C	Storage Humidity:	≤93%, no condensation
Dimensions and Weight			
Dimension:	125mmX60mmX30mm	Weight:	195g (excluding the batteries)
Power supply			
Alkaline Batteries	Power supply:	Three AA-size alkaline batteries	
	Working time:	More than 10 hours	
AC Adapter(Optional)	Input Voltage:	AC 100V~240V	
	Input Frequency:	50Hz~60Hz	
	Output Voltage:	DC 5V±5%	
	Output Current:	2A MAX	

Classifications	
Protection Class:	Type BF equipment applied part is SpO <sub>2</sub> probe
Protection Against harmful Ingress of Water:	IPX1
Alarm	
Alarm:	SpO <sub>2</sub> and pulse rate value, probe off, low battery
Alarm Mode:	Audio alarm, visual alarm and information

\*PI: Perfusion Index, endast för enheter sålda inom EU  
**Beijing Choice Electronic Tech Co., Ltd.**

Room 4104, No. A12 Yuquan Road, Haidian District,  
Beijing 100143, P.R.China  
Website: [www.choicemmed.com](http://www.choicemmed.com) Fax: 86 10-88204632  
Email: [market@choicemmed.com](mailto:market@choicemmed.com) Tel: 86 10-88204630/31  
NOTE The specifications are subject to changes without prior notice.  
Revised Date: Jan 11, 2019

**ChoiceMMed**  
*Providing a Better Life*

Svensk importör:  
ViTri Medical AB  
Stensättravägen 9a  
127 39 Skärholmen  
Tel: 08 7477660  
[www.vitri.se](http://www.vitri.se)

