

Enteral Data Sheet

Revision 3.1

Ref Doc: PSD0001		No.: DBISO-XX, DBISO-XXX, DBISO-XXXX & DBISON-500		Rev. 8	
Product:	Drainage Bag		Application:	Enteral	
Technical Qualification					
Raw Materials:	Drainage Bag - EVA Fim		Filtering Tube - EVA		
	Inlet and Outlet Tube - EVA		Female ENFit Connector - ABS		
	Male ENFit Connector - ABS		Tether - LDPE		
	Cap - HDPE		Air Venting Valve Cover - ABS		
	Air Venting Valve Housing - ABS		Ink – Mixture 701 Black		
	Air Venting Valve Filter Film - PTFE				
Packing:	Primary – PE and PET Secondary - Hard Corrugated Paper Tertiary - Double Corrugated Paper				
Regulatory:	Class I Sterile				
Others:	Shelf life – 3 years				
Products:					
Code:	Description:				
DBISO-50	Enteral Drainage Bag, 50ml, outlet tube with EN-Fit connector and cap (Sterile).				
DBISO-100	Enteral Drainage Bag, 100ml outlet tube with EN-Fit connector and cap (Sterile).				
DBISO-250	Enteral Drainage Bag, 250ml, outlet tube with EN-Fit connector and cap (Sterile).				
DBISO-500	Enteral Drainage Bag, 500ml, outlet tube with funnel end and cap (Sterile).				
DBISO-1000	Enteral Drainage Bag, 1000ml, outlet tube with funnel end and cap (Sterile).				
DBISON-500	Enteral Drainage Bag, 500ml outlet tube with EN-Fit connector and cap (Sterile).				
Intended Use:					
Enteral Drainage Bags are designed for the passive drainage of gastric fluid.					
Using Recommendations:					
See IFU					

Pictures:

Examples:



DBISO-500



DBISO-500

Packaging:

Code:	Packing:	Units:
DBISO-50	Primary – Bag 210mm(H) x 175mm(W) x 10mm(D)	1
	Secondary – Shelf Box 380mm(H) x 180mm(W) x 120mm(D)	20
	Tertiary – Shipper Carton 400mm(H) x 390mm(W) x 390mm(D)	120
DBISO-100	Primary – Bag 210mm(H) x 175mm(W) x 10mm(D)	1
	Secondary – Shelf Box 380mm(H) x 180mm(W) x 120mm(D)	20
	Tertiary – Shipper Carton 400mm(H) x 390mm(W) x 390mm(D)	120
DBISO-250	Primary – Bag 265mm(H) x 175mm(W) x 10mm(D)	1
	Secondary – Shelf Box 380mm(H) x 180mm(W) x 120mm(D)	20
	Tertiary – Shipper Carton 400mm(H) x 390mm(W) x 390mm(D)	120
DBISO-500	Primary – Bag 320mm(H) x 205mm(W) x 10mm(D)	1
	Secondary – Shelf Box 287mm(H) x 245mm(W) x 165mm(D)	20
	Tertiary – Shipper Carton 592mm(H) x 262m(W) x 360mm(D)	80

Enteral Data Sheet

Revision 3.1

DBISO-1000	Primary – Bag 350mm(H) x 252mm(W) x 10mm(D)	1
	Secondary – Shelf Box 287mm(H) x 245mm(W) x 165mm(D)	20
	Tertiary – Shipper Carton 592mm(H) x 262m(W) x 360mm(D)	80
DBISON-500	Primary – Bag 265mm(H) x 175mm(W) x 10mm(D)	1
	Secondary – Shelf Box 380mm(H) x 180mm(W) x 120mm(D)	20
	Tertiary – Shipper Carton 400mm(H) x 390mm(W) x 390mm(D)	120

Code:	GTIN Primary:	GTIN Secondary:	GTIN Tertiary:
DBISO-50	05060460843727	05060460843734	05060460843741
DBISO-100	05060460845707	05060460845714	05060460845721
DBISO-250	05060460843697	05060460843703	05060460843710
DBISO-500	05060460843758	05060460843765	05060460843772
DBISO-1000	05060460843666	05060460843673	05060460843680
DBISON-500	05060460848043	05060460848050	05060460848067


Code:	Quantity per Pallet:	Pallet Dimensions:
DBISO-50	2880	1200 x 800 mm
DBISO-100	2880	1200 x 800 mm
DBISO-250	2880	1200 x 800 mm
DBISO-500	2400	1200 x 800 mm
DBISO-1000	2400	1200 x 800 mm
DBISON-500	2880	1200 x 800 mm

Drawings:

THIS DRAWING IS CONFIDENTIAL AND THE PROPERTY OF THE COMPANY
“GBUK Group” IT CANNOT BE GIVEN, COMMUNICATED OR COPIED
WITHOUT OUR AUTHORIZATION


DIMENSIONS THE ON THE
DRAWING NOT TO SCALE

Example:




DBISO-250


Description: Product Print




DBISO-50




DBISO-100



DBISO-250



DBISO-500



DBISO-1000

DBISO-500

Description: Primary Artwork

Figure 1: DBISO-500

*** Variations: DBISO-50, DBISO-100, DBISO-250 ***

Figure 2: DBISO-1000

*** Variations: DBISO-500 ***



Drainage Bag

EVA Drainage Bag with short inlet tube, female ENFit® ISO 80369-3 connector and male ENFit® port & cap/tether

The device and contents must be discarded after use in accordance with local waste guidelines.

500ml

- FR** Poche de drainage entérale - Tube d'accès court - avec connecteur femelle ENFit® - port ENFit® mâle
Le dispositif et son contenu doivent être mis au rebut après usage conformément aux directives locales d'élimination des déchets.
- IT** Sacca di drenaggio enterale - Tubo di ingresso corto - con connettore femmina ENFit® - porta ENFit® maschio
Dopo l'utilizzo, il dispositivo e i suoi contenuti devono essere smaltiti in conformità con le normative locali sullo smaltimento dei rifiuti.
- ES** Bolsa de drenaje entérico - Tubo de entrada corto - con conector hembra ENFit® - puerto ENFit® macho
El dispositivo y todo el contenido deben eliminarse después de su uso según las normas locales sobre desechos.
- DE** Enteraler Drainagebeutel - Kurzer Einlassschlauch - mit weiblicher Verbindung ENFit® - männlicher ENFit®-Anschluss
Nach Gebrauch müssen der Drainagebeutel und der Inhalt nach den örtlich geltenden Abfallrichtlinien entsorgt werden.
- DA** Enteral drænagepose - Kort tilførselsrør - med hun-konnektor ENFit® - han-ENFit®-port
Udstyr og indhold skal bortskaffes efter brug i henhold til lokale retningslinjer.
- NO** Enteral drenasjepose - Kort innløpslange - med hunnkobling ENFit® - mannlig ENFit®-port
Enheten og innholdet må kastes etter bruk i henhold til lokale retningslinjer for avfall.
- SV** Enteral dränagepåse - Kort inloppslang - med honkoppling ENFit® - ENFit®-port för hane
Efter användning måste produkt och innehåll slängas i enlighet med lokala avfallsriktlinjer.
- NL** Enterale drainagezak - Korte invoerbuis - met vrouwelijke aansluiting ENFit® - mannelijke ENFit®-poort
Het apparaat en de inhoud moeten na gebruik worden weggegooid in overeenstemming met de plaatselijke afvalrichtlijnen.
- PL** Worek do drenażu dojelitowego - Krótka rurka wlotowa - z łącznikiem żeńskim ENFit® - męski port ENFit®
Urządzenie i zawartość należy wyrzucić po użyciu zgodnie z lokalnymi przepisami dotyczącymi odpadów.
- PT** Saco de drenagem entérica - Tubo de entrada curto - com conector fêmea ENFit® - porta ENFit® macho
O dispositivo e o conteúdo devem ser descartados após o uso de acordo com as diretrizes locais de descarte.
- EL** Εντερικός οσός παραχέτευσης - Κοντός σωλήνας εισόδου - με θηλυκό σύνδεσμο ENFit® - αρσενική θύρα ENFit®
Η συσκευή και το περιεχόμενο πρέπει να απορριφθούν μετά τη χρήση σύμφωνα με τις τοπικές οδηγίες για τα απόβλητα.
- TR** Enteral Drenaj Torbası - Kısa giriş tüpü - dişi konektörlü ENFit® - erkek ENFit® bağlantı noktası
Cihaz ve içeriği, kullanımdan sonra yerel atık yönetmelilerine uygun olarak atılmalıdır.
- CS** Enteralní drenážní vak - Krátká přívodní hadička - se samčím konektorem ENFit® - mužský port ENFit®
Zařízení a obsah musí být po použití zlikvidovány v souladu s místními předpisy o odpadu.
- ET** Enteraalne drenaažikott - Lühike sisselasketoru - sisekeermega liitmikuga ENFit® - meeste ENFit®-port
Pärast kasutamist tuleb seadme ja sisu ära visata vastavalt kohalikele jäätmekäitlussuhistele.
- HU** Enteralis levezető zacskó - Rövid bevezető cső - Csatlakozóhüvelyvel ENFit® - férfi ENFit® port
Az eszközt és tartalmát használat után a helyi hulladék irányleveleknek megfelelően meg kell semmisíteni.
- LV** Enterālais drenāžas maisis - Īsa ieviešanas caurule - ar līdzdas tipa savienotāju ENFit® - vīriešu ENFit® ports
Ierīce un saturs pēc lietošanas jāiznīcina saskaņā ar vietējām atkritumu vadlīnijām.
- MT** Borża għat-Tbattil tal-Intestini - Tubu qasir għad-dhul - b'konnnettur femminili ENFit® - port maskili ENFit®
L-apparat u l-kontenut għandhom jintremew wara l-użu skont il-linji gwida lokali dwar l-iskart.
- RO** Sac de drenaj enteral - Tub de admisie scurt - cu conector mamă ENFit® - port ENFit® masculin
Dispozitivul și conținutul trebuie aruncate după utilizare în conformitate cu instrucțiunile locale privind deșeurile.
- SL** Enteralna drenážna vrečka - Kratka dovodna cevka - z ženskim konektorjem ENFit® - moški ENFit® priključek
Napravo in vsebino po uporabi zavrzite v skladu z lokalnimi smernicami za ravnanje z odpadki.

GBUK Group Ltd.
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 282 945
US: 18665 289 422
www.gbukenteral.com
info@gbukenteral.com

GBUK GmbH
Wirtschaftsprüfungsgesellschaft
Im Klapperhof 7-23,
50670 Köln, Germany

EC REP **QualRep Services BV**
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP **Swiss AR Services AG**
Industriestrasse 47,
6300 Zug, Switzerland



Symbol glossary available at
www.gbukenteral.com
Not made with natural rubber latex.

ENFit® is a Registered Trademark of
GEDSA and used with permission.

Quantity:
1 units

REF DBISON-500

LOT 123456

YYYY-MM-DD

YYYY-MM-DD

Made in XXXX.



(01)05060460848043
(17)YYMMDD
(10)123456

DBISON-500-04 2410 21

Figure 1



Drainage Bag

EVA Drainage Bag with 50cm inlet tube, female ENFit® ISO 80369-3 connector and funnel connector to outlet

The device and contents must be discarded after use in accordance with local waste guidelines.

500ml

- FR** Poche de drainage entéral - Tube d'accès 50cm - avec connecteur femelle ENFit® - connecteur d'entonnoir
Le dispositif et son contenu doivent être mis au rebut après usage conformément aux directives locales d'élimination des déchets.
- IT** Sacca di drenaggio enterale - Tubo di ingresso 50cm - con connettore femmina ENFit® - connettore ad imbuto
Dopo l'utilizzo, il dispositivo e i suoi contenuti devono essere smaltiti in conformità con le normative locali sullo smaltimento dei rifiuti.
- ES** Bolsa de drenaje entérico - Tubo de entrada 50cm - con conector hembra ENFit® - conector de embudo
El dispositivo y todo el contenido deben eliminarse después de su uso según las normas locales sobre desechos.
- DE** Enteraler Drainagebeutel - 50cm Einlassschlauch - mit weiblicher Verbindung ENFit® - Trichteranschluss
Nach Gebrauch müssen der Drainagebeutel und der Inhalt nach den örtlich geltenden Abfallrichtlinien entsorgt werden.
- DA** Enteral drænagepose - 50cm tilførselsrør - med hun-konnektor ENFit® - tragstik
Udstyr og indhold skal bortskaffes efter brug i henhold til lokale retningslinjer.
- NO** Enteral drenasjepose - 50cm innløpslange - med hunnkobling ENFit® - traktkontakt
Enheter og innholdet må kastes etter bruk i henhold til lokale retningslinjer for avfall.
- SV** Enteral dränagepåse - 50cm inloppslang - med honkoppling ENFit® - trattanslutning
Efter användning måste produkt och innehåll slängas i enlighet med lokala avfallsriktlinjer.
- NL** Enterale drainagezak - 50cm invoerbuis - met vrouwelijke aansluiting ENFit® - trechter connector
Het apparaat en de inhoud moeten na gebruik worden weggegooid in overeenstemming met de plaatselijke afvalrichtlijnen.
- PL** Worek do drenażu dojelitowego - 50cm rurka wlotowa - z łącznikiem żeńskim ENFit® - złącze lejka
Urządzenie i zawartość należy wyrzucić po użyciu zgodnie z lokalnymi przepisami dotyczącymi odpadów.
- PT** Saco de drenagem entérica - Tubo de entrada 50cm - com conector fêmea ENFit® - conector de funil
O dispositivo e o conteúdo devem ser descartados após o uso de acordo com as diretrizes locais de descarte.
- EL** Εντερικός σάκος παροχέτευσης - 50cm σωλήνας εισόδου - με θηλυκό σύνδεσμο ENFit® - σύνδεσμος διοχέτευσης
Η συσκευή και το περιεχόμενό της πρέπει να απορριφθούν μετά τη χρήση σύμφωνα με τις τοπικές οδηγίες για τα απόβλητα.
- TR** Enteral Drenaj Torbası - 50cm giriş tüpü - dişi konektörlü ENFit® - huni konektörü
Çihaz ve içeriği, kullanımdan sonra yerel atık yönetmeliklerine uygun olarak atılmalıdır.
- CS** Enterální drenážní vak - 50cm přírodní hadička - se samičím konektorem ENFit® - trychtýřový konektor
Zařízení a obsah musí být po použití zlikvidovány v souladu s místními předpisy o odpadu.
- ET** Enteraalne drenaažikott - 50cm sisselasketoru - sisekeermeega liitmikuga ENFit® - lehtri pistik
Pärast kasutamist tuleb seade ja sisu ära visata vastavalt kohalikele jäätmekäitlusjuhistele.
- HU** Enterális levezető zacskó - 50cm bevezető cső - Csatlakozóhüvelyvel ENFit® - tölcésér csatlakozó
Az eszközt és tartalmát használat után a helyi hulladék irányelveknek megfelelően meg kell semmisíteni.
- LV** Enterālās drenāžas maisis - 50cm ieplūdes caurule - ar ligzdas tipa savienotāju ENFit® - piltuves savienotājs
Ierīce un saturs pēc lietošanas jāiznīcina saskaņā ar vietējām atkritumu vadlīnijām.
- MT** Borża għat-Tbattli tal-Intestini - Tubu 50cm għad-dhul - b'konnettur femminili ENFit® - konnettur tal-lembut
L-apparat u l-kontenut għandhom jintremew wara l-użu skont il-linji gwida lokali dwar l-iskart.
- RO** Sac de drenaj enteral - Tub de admisie 50cm - cu conector mamă ENFit® - conector pălnie
Dispozitivul și conținutul trebuie aruncate după utilizare în conformitate cu instrucțiunile locale privind deșeurile.
- SL** Enteralna drenažna vrečka - 50cm dovodna cevka - z ženskim konektorjem ENFit® - prikjuček lijaka
Napravo in vsebino po uporabi zavrzite v skladu z lokalnimi smericami za ravnanje z odpadki.

GBUK Group Ltd.
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 282 945
US: 18665 289 422
www.gbukenteral.com
info@gbukenteral.com

GBUK GmbH
Wirtschaftsprüfungsgesellschaft
Im Klapperhof 7-23,
50670 Köln, Germany

EC REP **QualRep Services BV**
Ulrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP **Swiss AR Services AG**
Industriestrasse 47,
6300 Zug, Switzerland



CE 1639 **UK** CA0120

Rx Only Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Symbol glossary available at
www.gbukenteral.com
Not made with natural rubber latex.

ENFit® is a Registered Trademark of
GEDSA and used with permission.

Quantity:
1 units

REF **DBISO-500**

LOT 123456

YYYY-MM-DD

YYYY-MM-DD

Made in XXXX.



(01)05060460843758
(17)YYMMDD
(10)123456

DBISO-500-04 2410 21

Figure 2


Description: Secondary Artwork

Figure 3: DBISON-500

*** Variations: DBISO-50, DBISO-100, DBISO-250 ***

Figure 4: DBISO-500

*** Variations: DBISO-1000 ***




Drainage Bag
EVA Drainage Bag with short inlet tube, female ENFit® ISO 80369-3 connector and male ENFit® port & cap/tether

500ml

FR Poche de drainage entéral - Tube d'accès court - avec connecteur femelle ENFit® - port ENFit® mâle
IT Sacca di drenaggio enterale - Tubo di ingresso corto - con connettore femmina ENFit® - porta ENFit® maschio
ES Bolsa de drenaje entérico - Tubo de entrada corto - con conector hembra ENFit® - puerto ENFit® macho
DE Enteraler Drainagebeutel - Kurzer Einlassschlauch - mit weiblicher Verbindung ENFit® - männlicher ENFit®-Anschluss
DA Enteral drænagepose - Kort tilførselsrør - med hun-konnektor ENFit® - han-ENFit®-port
NO Enteral drenasjepose - Kort innløpsslange - med hunnkobling ENFit® - mannlig ENFit®-port
SV Enteral dränagepåse - Kort inloppslang - med honkoppling ENFit® - ENFit®-port för hane
NL Enterale drainagezak - Korte invoerbuis - met vrouwelijke aansluiting ENFit® - mannelijke ENFit®-poort
PL Worek do drenażu dojelitowego - Krótka rurka wlotowa - z łącznikiem żeńskim ENFit® - męski port ENFit®
PT Saco de drenagem entérica - Tubo de entrada curto - com conector fêmea ENFit® - porta ENFit® macho
EL Εντερικός σάκος παροχέτευσης - Κοντός σωλήνας εισόδου - με θηλυκό σύνδεσμο ENFit® - αρσενική θύρα ENFit®
TR Enteral Drenaj Torbası - Kısa giriş tüpü - dişi konektörü ENFit® - erkek ENFit® bağlantı noktası
CS Enteralní drenážní vak - Krátká přírodní hadička - se samičím konektorem ENFit® - mužský port ENFit®
ET Enteraalne drenaažikott - Lühike sisselasketoru - sisekeermega liitmikuga ENFit® - meeste ENFit®-port
HU Enteralis levezető zacskó - Rövid bevezető cső - Csatlakozóhüvelyvel ENFit® - férfi ENFit® port
LV Enteralās drenāžas maiss - Īsa ievādes caurule - ar ligzdas tipa savienotāju ENFit® - vīriešu ENFit® ports
MT Borża għat-Tabattil tal-Intestini - Tubu qasir għad-dhul - b'konnettur femminili ENFit® - port maskili ENFit®
RO Sac de drenaj enteral - Tub de admisie scurt - cu conector mamă ENFit® - port ENFit® masculin
SL Enteralna drenažna vrečka - Kratka dovodna cevka - z ženskim konektorjem ENFit® - moški ENFit® priključek

Rx Only Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.




UK 1639 CA0120 **STERILE EO**

Symbol glossary available at www.gbukenteral.com
Not made with natural rubber latex.
ENFit® is a Registered Trademark of GEDSA and used with permission.

GBUK Group Ltd.
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 282 945
US: 18665 289 422
www.gbukenteral.com
info@gbukenteral.com

GBUK GmbH
Wirtschaftsprüfungsgesellschaft
Im Klapperhof 7-23,
50670 Köln, Germany



QualRep Services BV
Ulrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands



Swiss AR Services AG
Industriestrasse 47,
6300 Zug, Switzerland

Quantity:

20 units

REF DBISON-500

LOT 123456

YYYY-MM-DD

YYYY-MM-DD

Made in XXXX.



(01)05060460848050

(17)YYMMDD

(10)123456

DBISON-500-02 2410 22

Figure 3



Drainage Bag

EVA Drainage Bag with 50cm inlet tube, female ENFit® ISO 80369-3 connector and funnel connector to outlet

500ml

- FR Poche de drainage entérale - Tube d'accès 50cm - avec connecteur femelle ENFit® - connecteur d'entonnoir
- IT Sacca di drenaggio enterale - Tubo di ingresso 50cm - con connettore femmina ENFit® - connettore ad imbuto
- ES Bolsa de drenaje entérico - Tubo de entrada 50cm - con conector hembra ENFit® - conector de embudo
- DE Enteraler Drainagebeutel - 50cm Einlassschlauch - mit weiblicher Verbindung ENFit® - Trichteranschluss
- DA Enteral drænagepose - 50cm tilførselsrør - med hun-konnektor ENFit® - tragtstik
- NO Enteral drænagepose - 50cm innløpsslange - med hunnkobling ENFit® - traktkontakt
- SV Enteral dränagepåse - 50cm inloppslang - med honkoppling ENFit® - trattanslutning
- NL Enteral drainagezak - 50cm invoerbuis - met vrouwelijke aansluiting ENFit® - trechter connector
- PL Worek do drenażu dojelitowego - 50cm rurka wlotowa - z łącznikiem żeńskim ENFit® - złączce lejka
- PT Saco de drenagem entérica - Tubo de entrada 50cm - com conector fêmea ENFit® - conector de funil
- EL Εντερικός κύριος παροχέτευσης - 50cm σωλήνας εισόδου - με θηλυκό σύνδεσμο ENFit® - σύνδεσμος διαχέτευσης
- TR Enteral Drenaj Torbası - 50cm giriş tüpü - dişi konektörlü ENFit® - huni konektörü
- CS Enteralní drenážní vak - 50cm přírodní hadička - se samičím konektorem ENFit® - trychťřívový konektor
- ET Enteralne drenaažkott - 50cm sisselasketoru - sisekeermeega liitmikuga ENFit® - lehtri pistik
- HU Enteralis levezető zacskó - 50cm bevezető cső - Csatlakozóhüvellyel ENFit® - tölszér csatlakozó
- LV Enteralās drenāžas maisis - 50cm ieplūdes caurule - ar ligzdas tipa savienotāju ENFit® - piltuves savienotājs
- MT Borża għat-Tabattil tal-Intestini - Tubu 50cm għad-dhul - b'konnettur femminili ENFit® - konnettur tal-fembut
- RO Sac de drenaj enteral - Tub de admisie 50cm - cu conector mamă ENFit® - conector pălnie
- SL Enteralna drenažna vrečka - 50cm dovodna cevka - z ženskim konektorjem ENFit® - priključek lijaka

Rx Only Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.



CE 1639 UK CA0120 STERILE EO

Symbol glossary available at www.gbukenteral.com
Not made with natural rubber latex.
ENFit® is a Registered Trademark of GEDSA and used with permission.

GBUK Group Ltd.
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 282 945
US: 18665 289 422
www.gbukenteral.com
info@gbukenteral.com

GBUK GmbH
Wirtschaftsprüfungsgesellschaft
Im Klappertof 7-23,
50670 Köln, Germany

QualRep Services BV
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

Swiss AR Services AG
Industriestrasse 47,
6300 Zug, Switzerland

EC REP

CH REP

Quantity:
20 units

REF **DBISO-500**

LOT **123456**

YYYY-MM-DD

YYYY-MM-DD

Made in XXXX.



(01)05060460843765
(17)YYMMDD
(10)123456

DBISO-500-02 2410 21

Figure 4

Description: Tertiary Artwork

Figure 5: DBISON-500
*** Variations: DBISO-50, DBISO-100, DBISO-250 ***

Figure 6: DBISO-500
*** Variations: DBISO-1000 ***



GBUK Group Ltd.
Woodland House,
Blackwood Hill Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 282 945
US: 85665 289 422
www.gbukenteral.com
info@gbukenteral.com

GBUK GmbH
Wirtschaftsprüfungsgesellschaft
im Kappelerhof 7, 23,
50670 Köln, Germany

EC REP: **QualRep Services BV**
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP: **Swiss AR Services AG**
Industriestrasse 47,
6300 Zug, Switzerland

REF DBISON-500

Enteral ISOSAF™ Drainage Bag, EVA Drainage Bag with short inlet tube, female ENFit® ISO 80369-3 connector and male ENFit® port & cap/tether

500ml

LOT 123456

YYYY-MM-DD

YYYY-MM-DD

Made in XXXXX

Quantity:
120 units

Weight:
XX kg
XX lb XX oz



Rx Only Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician. Symbol glossary available at www.gbukenteral.com. Not made with natural rubber latex. ENFit® is a Registered Trademark of GEDSA and used with permission.




Figure 5



GBUK Group Ltd.
Woodland House,
Blackwood Hill Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 282 945
US: 85665 289 422
www.gbukenteral.com
info@gbukenteral.com

GBUK GmbH
Wirtschaftsprüfungsgesellschaft
im Kappelerhof 7, 23,
50670 Köln, Germany

EC REP: **QualRep Services BV**
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP: **Swiss AR Services AG**
Industriestrasse 47,
6300 Zug, Switzerland

REF DBISO-500

Enteral ISOSAF™ Drainage Bag, EVA Drainage Bag with 50cm inlet tube, female ENFit® ISO 80369-3 connector and funnel connector to outlet

500ml

LOT 123456

YYYY-MM-DD

YYYY-MM-DD

Made in XXXXX

Quantity:
80 units

Weight:
XX kg
XX lb XX oz



Rx Only Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician. Symbol glossary available at www.gbukenteral.com. Not made with natural rubber latex. ENFit® is a Registered Trademark of GEDSA and used with permission.



Figure 6



GBUK Group Ltd.
Woodland House,
Blackwood Hill Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 282 945
US: 85665 289 422
www.gbukenteral.com
info@gbukenteral.com

GBUK GmbH
Wirtschaftsprüfungsgesellschaft
im Kappelerhof 7, 23,
50670 Köln, Germany

EC REP: **QualRep Services BV**
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP: **Swiss AR Services AG**
Industriestrasse 47,
6300 Zug, Switzerland

REF DBISON-500

Enteral ISOSAF™ Drainage Bag, EVA Drainage Bag with short inlet tube, female ENFit® ISO 80369-3 connector and male ENFit® port & cap/tether

500ml

LOT 123456

YYYY-MM-DD

YYYY-MM-DD

Made in XXXXX

Quantity:
120 units

Weight:
XX kg
XX lb XX oz



Rx Only Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician. Symbol glossary available at www.gbukenteral.com. Not made with natural rubber latex. ENFit® is a Registered Trademark of GEDSA and used with permission.



GBUK Group Ltd.
Woodland House,
Blackwood Hill Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 282 945
US: 85665 289 422
www.gbukenteral.com
info@gbukenteral.com

GBUK GmbH
Wirtschaftsprüfungsgesellschaft
im Kappelerhof 7, 23,
50670 Köln, Germany

EC REP: **QualRep Services BV**
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP: **Swiss AR Services AG**
Industriestrasse 47,
6300 Zug, Switzerland

REF DBISO-500

Enteral ISOSAF™ Drainage Bag, EVA Drainage Bag with 50cm inlet tube, female ENFit® ISO 80369-3 connector and funnel connector to outlet

500ml

LOT 123456

YYYY-MM-DD

YYYY-MM-DD

Made in XXXXX




Quantity:
80 units

Weight:
XX kg
XX lb XX oz



Rx Only Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician. Symbol glossary available at www.gbukenteral.com. Not made with natural rubber latex. ENFit® is a Registered Trademark of GEDSA and used with permission.



Modified by:	Name: Aleksandra Malek	Signature:  Signed by Aleksandra Malek Aleksandra Malek	Date: I have reviewed this document 09-Dec-24 9:18:56 AM GMT 09-Dec-24
Reviewed by:	Name: Nicola Neal	Signature:  Signed by Nicola Neal Nicola Neal	Date: I approve this document 09-Dec-24 15:10 PM GMT 09-Dec-24
 <p>GBUK Enteral is a trading name of GBUK Enteral Ltd. GBUK Enteral Ltd. Is part of GBUK Group Ltd.</p>		<p>GBUK Enteral Woodland House Blackwood Hall Business Park North Duffield, North Yorkshire YO8 5DD, UK</p> <p>TEL: +44 (0) 1757 282 945 https://www.gbukgroup.com/ https://www.gbukgroup.com/brands/gbuk-ental/</p>	