

Rosidal® 1C

REF 162154 - 162157

April 2025

1. Composition of the product

The Rosidal® 1C consists of:

- 48% polyamide
- 42% cotton
- 6% polyester
- 4% elastane

This product data sheet covers the following item:

REF	Produkt description	feature
162154	Rosidal® 1C	8 cm / 6 m stretched, with fixed ends and transport lock
162155	Rosidal® 1C	8 cm / 6 m stretched, with fixed ends and transport lock
162156	Rosidal® 1C	10 cm / 6 m stretched, with fixed ends and transport lock
162157	Rosidal® 1C	10 cm / 6 m stretched, with fixed ends and transport lock

2. Packaging, structure and composition

2.1 Unit container

see 2.2

2.2 Shelf container

1 cellulose folding box
1 instructions for use (cellulose)

REF	pieces
162154, 162156	1
162155, 162157	6

2.3 Transit container

1 corrugated cardboard box (cellulose):

REF	pieces	Folding boxes
162154	32	32

REF	pieces	Folding boxes
162155, 162157	24	4
162156	60	60

3. Manufacture

Rosidal® 1C is produced according to specification in hygienic conditions and packed as described in its relevant packaging specification.

4. Description

Rosidal® 1C is a bi-elastic bandage made by warp knitting technology. Due to its structure, the short stretch bandage combines three layers in one product: skin layer, padding layer and compression layer.

5. Properties

- short-stretch characteristics
- high working and low resting pressure
- textile terry structure
- the slim design allows it to be worn in most types of shoes in comparison to multi component systems
- air- and water vapour permeable
- may be worn directly on intact skin
- Oekotex certified yarn
- does not contain latex*

**Natural rubber latex is not an intended ingredient in the production of the above mentioned product or their packaging. We cannot completely exclude any potential contamination with traces of latex from production of the environment.*

6. Intended purpose (see valid instructions for use)

The Rosidal® 1C is used for compression, retention of dressings, support and relief and for immobilization of limbs.

7. Medical device classification

The Rosidal® 1C is a medical device of Class I according to Rule 1 (according to the Regulation (EU) 2017/745 on medical devices (MDR)).

8. Biological evaluation and biocompatibility

The biological risk evaluation performed according to ISO 10933-1 confirms the safety and acceptable biocompatibility of the Rosidal® 1C within its intended purpose.

9. Stability

When appropriate stored (keep dry, protected from dust and sun) the Rosidal® 1C has a shelf life of 5 years from the date of manufacture.

10. Disposal

The user is advised to observe current national legislation, norms and guidelines, regulating the disposal of medical refuse.

Packaging materials must also be disposed of in compliance with applicable national requirements.

Lohmann & Rauscher International GmbH & Co. KG
D-56579 Rengsdorf
signed by

Michael Hamedinger
Head of Global Marketing Bandages & Supports

Tanja Geiermann
Global Brand Manager Bandages & Supports