


Ref Doc: PSD0014A		No.: SBAISO		Rev. 1	
Product:	Stepped Bottle Adaptor		Application:	Enteral	
Technical Qualification					
Raw Materials:	Body PE, Core PUR				
Packing:	Single Blistered. Shelf Box 50				
Regulatory:	Class I Non Sterile				
Others:	Shelf Life 5 Years Phthalates and DEHP Free				
Products:					
Code:	Description:				
SBAISO	Stepped Bottle Adaptor for Oral and Enteral Syringes				
Intended Use:					
Used with ISO 80369-3 ENFit® Enteral syringes and/or GBUK Oral Syringes for putting into bottles and drawing enteral feed or medication from a bottle					
Using Recommendations:					
Can be attached to both GBUK Oral Syringes and Enteral ENFit Syringes via 2 different connectors. Single Use.					
Associated Devices:					
ENFit Enteral Syringes GBUK Oral Syringes					
Production Environment:					
Cleanroom					
Pictures:					
					

Packaging:			
Code:	Packing:	Units:	
SBAISO	Primary Blister 98mm x 127mm	1	
	Secondary - Shelf Box 175mm x 245mm x 490mm	50	
	Tertiary – Shipper box 375mm x 265mm x 510mm	100	
Code:	GTIN Primary:	GTIN Secondary:	GTIN Tertiary:
SBAISO	05060590883945	05060590883952	05060590886618
Code:	Quantity per Pallet:	Pallet Dimensions:	
SBAISO	TBC	1200 x 800mm	
Drawings:			
THIS DRAWING IS CONFIDENTIAL AND THE PROPERTY OF THE COMPANY "GBUK Group" IT CANNOT BE GIVEN, COMMUNICATED OR COPIED WITHOUT OUR AUTHORIZATION			DIMENSIONS THE ON THE DRAWING NOT TO SCALE
Description:	Primary Artwork		


ISOSAF™

Stepped Bottle Adaptor

ISO 80369-3 ENFit®

Used with ISO 80369-3 ENFit® Enteral syringes for putting into bottles and drawing enteral feed or medication from a bottle.

Bottle adaptors are not interchangeable between medicine bottles once in use; use one bottle adaptor per medicine bottle, dispose once detached.




Part of the GBUK Group

ENFit® is a Registered Trademark of GEDSA and used with permission.

GBUK Group Ltd.
Woodland House
Blackwood Hall Business Park,
North Duffield, Selby, North Yorkshire
YO8 5DD, UK

EU: +44 (0)1757 282 945
US: 18665 289 422
info@gbukenteral.com
www.gbukenteral.com



Symbol glossary available at www.gbukenteral.com

REF **SBAISO**

LOT **123456**

YYYY-MM-DD


Made in XXXXX

EXAMPLE BARCODE

(01)05060590883945

(11)YYMMDD

(10)123456

SBAISO-04 

Description: Secondary Artwork


ISOSAF™

50 units

Stepped Bottle Adaptor

ISO 80369-3 ENFit®


Used with ISO 80369-3 ENFit® Enteral syringes for putting into bottles and drawing enteral feed or medication from a bottle.



Part of the GBUK Group

GBUK Group Ltd.
Woodland House
Blackwood Hall Business Park,
North Duffield, Selby, North Yorkshire
YO8 5DD, UK

EU: +44 (0)1757 282 945
US: 18665 289 422
info@gbukenteral.com
www.gbukenteral.com



Symbol glossary available at www.gbukenteral.com

REF **SBAISO**

LOT **123456**

YYYY-MM-DD


Made in XXXXX

EXAMPLE BARCODE

(01)05060590883945

(11)YYMMDD

(10)123456

SBAISO-04 

Description: Tertiary Artwork

3 | Page

QF191

Issue 2

	<p>GBUK Group Ltd. Woodland House, Blackwood Hall Business Park, North Duffield, Selby, North Yorkshire, YO8 5DD, UK EU: +44 (0)1757 282 945 US: 18665 289 422 info@gbukenteral.com www.gbukenteral.com</p>
<p>REF SBAISO</p>	<p>Quantity: 100 units</p>
<p>ISOSAF™ Stepped Bottle Adaptor, ISO 80369-3 ENFit®</p>	<p>Weight: XX kg XX lb XX oz</p>
<p>LOT 123456 YYYY-MM-DD</p>	
<p>Made in XXXXX. </p>	
<p>Symbol glossary available at www.gbukenteral.com ENFit® is a Registered Trademark of GEDSA and used with permission.</p>	
<div style="border: 1px solid black; border-radius: 15px; padding: 10px; text-align: center;">  <p>EXAMPLE BARCODE</p> <p>(01)05060590886618 (11)YYMMDD (10)123456</p> <p style="font-size: small;">SBAISO-08 2007</p> </div>	

Modified by: Nicola Robertson	Date: 30.10.20
Reviewed by: Nicola Robertson	Date: 30.10.20

<p> Enteral UK is a trading name of GBUK Enteral Ltd. ™ Trademark of GBUK Enteral Ltd GBUK Enteral Ltd. Is part of GBUK Group Ltd.</p>	<p>GBUK Enteral Ltd. Blackwood Hall Business Park North Duffield, Selby, North Yorkshire Selby, YO8 5DD, UK TEL: +44 (0) 1757 282 945 FAX: +44 (0) 1757 600 545 www.gbukenteral.com</p>
---	--