




CENAMAN

Single use ENFit™ syringes

Cenaman® ENFit syringe - single use are syringes used for delivering liquid medication and enteral nutrition orally or via a feeding tube into the gastro-intestinal tract, hydration and for flushing of enteral feeding devices.

Classification

The Cenaman ENFit™ syringes reduce the risk of an unintended misconnection. All syringes are equipped with the safe enteral connector ENFit™, according to international safety standard ISO 80369-3, thus offering the highest safety in the clinical practice.

Class Is under the Medical Devices Regulations (SI 2002 No. 618, as amended) and Class Is under Regulation (EU) 2017/745 Medical Device Regulation (MDR).

Packaging and product reference

REF code	Size	description	pcs per inner box	pcs per outer box
C062001	1 ml	Cenaman Single-use EnFit syringe 1 ml	30 (30 x 1)	900 (30x30)
C062002	2,5 ml	Cenaman Single-use EnFit syringe 2,5 ml	30 (30 x 1)	900 (30x30)
C062005	5 ml	Cenaman Single-use EnFit syringe 5 ml	30 (30 x 1)	900 (30x30)
C062010	10 ml	Cenaman Single-use EnFit syringe 10 ml	30 (30 x 1)	900 (30x30)
C062020	20 ml	Cenaman Single-use EnFit syringe 20 ml	30 (30 x 1)	480 (16x30)
C062060	60 ml	Cenaman Single-use EnFit syringe CENTRIC TIP 60 ml	30 (30 x 1)	120 (4x30)
C062061	60 ml	Cenaman Single-use EnFit syringe ECCENTRIC TIP 60 ml	30 (30 x 1)	120 (4x30)
C062100	100 ml	Cenaman Single-use EnFit syringe 100 ml	30 (30 x 1)	120 (4x30)

REF code	Size	description	pcs per inner box	pcs per outer box
C062001	1 ml	Cenaman Single-use EnFit syringe 1 ml	16971564466414	26971564466411
C062002	2,5 ml	Cenaman Single-use EnFit syringe 2,5 ml	16971564466421	26971564466428
C062005	5 ml	Cenaman Single-use EnFit syringe 5 ml	16971564466438	26971564466435
C062010	10 ml	Cenaman Single-use EnFit syringe 10 ml	16971564466445	26971564466442
C062020	20 ml	Cenaman Single-use EnFit syringe 20 ml	16971564466452	26971564466459
C062060	60 ml	Cenaman Single-use EnFit syringe CENTRIC TIP 60 ml	16971564466469	26971564466466
C062061	60 ml	Cenaman Single-use EnFit syringe ECCENTRIC TIP 60 ml	16971564466544	26971564466541
C062100	100 ml	Cenaman Single-use EnFit syringe 100 ml	16971564466476	26971564466473

Features

- Enfit connector (ISO 80369-3)
- Clear number markings
- High visibility through barrel
- Purple coloured plunger
- Low dose tip in 1ml, 2.5ml and 5ml
- Smooth plunger movement
- Sterile (EO sterilized)

Instructions for use



Intended for the delivery of liquid medication, enteral feed and hydration.
Do not autoclave or sterilize syringe as this may cause syringe to malfunction.

Absence of specific components

Free of latex, phthalates, substances of animal and/or human origin, blood or blood derivatives, medicinal substances, heavy metals and PVC.

Shelflife

Cenaman® ENFit syringe – single patient use have a shelf-life of 5 years from production date, when stored under specified conditions.