



ACBS ☒ GMS ☒ PBS ☒

## Description

### Food for Special Medical Purposes

K.Quik is a ready to use, liquid emulsion of medium chain triglycerides (MCT).  
With sweetener.

## Indications

For the dietary management of conditions requiring a source of medium chain triglycerides (MCT) including use in the ketogenic diet.

## Recommended Intake and Administration

To be determined by the clinician or dietitian and is dependent on the age, body weight and medical condition of the patient.

## Preparation Guidelines

Shake well before use.

Best served chilled.

Further serving suggestions and recipes for K.Quik are available on request.

### IMPORTANT NOTICE

Use under medical supervision.

Not suitable for use as a sole source of nutrition.

Suitable from 3 years of age onwards.

Not suitable for people with medium chain acyl-CoA dehydrogenase deficiency (MCADD).

For enteral use only.

## Nutritional Information

		per 100ml	per 225ml (1 bottle)			per 100ml	per 225ml (1 bottle)
Energy	kJ	777	1750	Protein	g	0	0
	kcal	189	426	Salt	g	0.13	0.29
Fat	g	21	47.3	<b>Minerals</b>			
of which				Sodium	mg	52	117
saturates	g	19.8	44.6		mmol	2.2	5.0
MCT *	g	20	45	Chloride	mg	60	135
LCT **	g	1.0	2.3		mmol	1.7	3.8
Carbohydrate	g	0	0	<b>Osmolality</b>			
of which sugars	g	0	0	50 mOsm/Kg			
				<b>Osmolarity</b>			
				40 mOsm/L			

\*Medium Chain Triglycerides

\*\*Long Chain Triglycerides

### Storage

Unopened: Store in a cool dry place.

Once opened: Recap unused portion, refrigerate and use within 48 hours.

Packaged in a protective atmosphere.

### Flavour

Unflavoured

### Sample Box

2 x 225ml = 450ml

1 measuring cup

### Pack Size/Weight

15 x 225ml e bottles = 3.37 L

Mediq varenr. 7404047

