

glycosade[®]

Disorders of Carbohydrate Metabolism

Glycosade Lemon

ACBS ☒ GMS ☐ PBS ☐

Description

Food for special medical purposes, with sweetener.

Glycosade[®] lemon is a lemon flavoured, hydrothermally processed high amylopectin waxy maize starch.

Indications

Glycosade lemon is for the dietary management of hepatic glycogen storage disease (GSD) where the use of long acting starch is indicated, from 3 years of age.

Dosage and Administration

To be determined by the clinician or dietitian and is dependent on the age, bodyweight and metabolic response of the patient. A fasting Glycosade load is recommended.

Preparation Guidelines

To be determined by the clinician or dietitian. Glycosade can be either taken as a drink or added to cold food.

DO NOT HEAT OR WARM GLYCOSADE AS THIS DESTROYS THE PROPERTIES OF THE STARCH.

Recommended preparation



1. Measure the required amount of liquid into a beaker with a lid.



2. Add the required amount of Glycosade lemon to the beaker. Typically a sachet of Glycosade can be added to 100ml of cold water or other suitable fluids.



3. Secure lid and shake until the powder is dissolved.



4. Glycosade lemon is now ready to use. The product should be freshly prepared and consumed immediately.

Glycosade may also be given through a feeding tube. Administer as a bolus to prevent possible tube blockage.

IMPORTANT NOTICE

Must only be used in individuals with proven hepatic GSD where the use of long acting starch is indicated under strict medical supervision

Not for use as a sole source of nutrition.

Suitable from 3 years of age.

For enteral use only.

Nutritional Information

		per 100g	per 60g sachet			per 100g	per 60g sachet
Energy	kJ	1496	901	Minerals			
	kcal	352	212		Sodium	mg	<200
Fat	g	0	0		mmol	<0.90	<0.50
of which saturates	g	0	0				
Carbohydrate	g	88	53	Osmolality: 40 mOsm/kg (1 x 60g sachet + 100ml water)			
of which sugars	g	0	0				
Protein	g	0	0				
Salt	g	<0.50	<0.30				

Storage

Store in a cool, dry place. Sachets are designed for single use and should be used immediately once opened.

If your dietitian/clinician prescribes less than one full sachet, any unused powder in a sachet must be tightly sealed, and used within 24 hours.

Pack size/Weight

30 x 60g sachets = 1.8kg

Flavour

Lemon

