

Product Specifications

Brand : Klinion®
Product : Kliniplast® Fix Stretch

Product description

Kliniplast® Fix is a non-sterile, elastic, adhesive nonwoven fixing plaster on a roll. Nonwoven backing and solvent pressure sensitive adhesive layer.

Intended use

Elastic fixation of dressings, catheters and drains.

Classification

Class I under the Medical Devices Regulations (SI 2002 No. 618, as amended) and
Class I under Regulation (EU) 2017/745 Medical Device Regulation (MDR).

Packaging & Product Reference

Product Reference	Description	Size (*)	Piece(s)	Inner Carton (pcs)	Outer Carton (pcs)
40294340	Kliniplast® Fix Stretch 10 m x 5 cm	5CM X 10M	N/A	1	48
40294341	Kliniplast® Fix Stretch 10 m x 10 cm	10CM X 10M	N/A	1	24
40294342	Kliniplast® Fix Stretch 10 m x 15 cm	15CM X 10M	N/A	1	24
40294343	Kliniplast® Fix Stretch 10 m x 20 cm	20CM X 10M	N/A	1	16

*Sizes are approximate

Product Reference	Description	Size (*)	Piece(s)	EAN Inner Carton	EAN Outer Carton
40294340	Kliniplast® Fix Stretch 10 m x 5 cm	5CM X 10M	N/A	8715343747582	8715343747575
40294341	Kliniplast® Fix Stretch 10 m x 10 cm	10CM X 10M	N/A	8715343747568	8715343747551
40294342	Kliniplast® Fix Stretch 10 m x 15 cm	15CM X 10M	N/A	8715343747544	8715343747537
40294343	Kliniplast® Fix Stretch 10 m x 20 cm	20CM X 10M	N/A	8715343747520	8715343747513

*Sizes are approximate

Product Specifications

Features

- Elastic, soft and comfortable
- Easy to cut to size

Precautions

- Single use
- Non-sterile

Specifications

Component	Material
Backing	White water repellent nonwoven 100% polyester - 50g/m ²
Adhesive	Transparent solvent pressure sensitive adhesive - Ethyl acetate + 2-ethyl hexyl acrylate
Release liner	White glassine paper - 87% wood pulp + 8% water + 5% organosilicon

Performance Characteristic	Specification	Test method
180° Peel Strength	≥2.5N/2,5 cm	ISO 29862-2018
Holding power	≥5min	ISO 29863-2018
Tensile strength machine direction	≥35%	EN 29073-3
Tensile strength cross direction	≥90%	EN 29073-3

Product Specifications

Presence of specific components

Component	Present
Latex	No
Phthalates	No
Substances of animal origin	No
Substances of human origin	No
Blood or blood derivatives	No
Medicinal substances	No

Sterilization

Non-sterile.

Storage conditions

The products should be stored in the original packaging, in cool, well-ventilated and dry place, protected from direct sunlight.

Shelf life

Kliniplast® Fix Stretch has a shelf-life of 5 years from production date, when stored under specified conditions.

Product Specifications

Picture



Product ref.: 40294341

Approval	Name	Date	Signature
Product Manager Wound Care	Florin Lazar	23.05.2025	
RA Specialist	Emil Laske	23.05.2025	