



SMOKE EVACUATION PENCILS



510(k)NO.: registered with the U.S. FDA



EN ISO 9001:2008
EN ISO 13485:2003+AC:20

Notice: Single use to avoid cross-infection

SURGICAL SMOKE EVACUATION PENCILS


Telescopic



2B-D013-SG0035

telescopic, 70mm coated blade
3.5m hose, 3m cable
with holster
sterile, single use

Standard



2A-D009-SG0072

70mm coated needle electrodes
3m hose, 3m cable,
with holster, sterile, single use

2A-D009-SG0073

70mm coated blade
3m hose, 3m cable,
with holster, sterile, single use

2A-D009-SG0074

70mm coated blade
3.5m hose, 3.5m cable,
with holster, sterile, single use

Hose



3A-D000-SG0042

extraction hose
2 meters
sterile, single use

- **Large flow rate, which can be as high as 77 LPM. The effect of smoke evacuation is perfect, which protects the health of surgeons to the greatest extent**
- **Ultra-long telescopic distance, which is suitable for different surgical sites during the operation**
- **Non-slip design on the hand-held part, which makes the hand-held more stable**
- **The 360°rotating design at the tail, which makes the operation comfortable and flexible**

 **SHUYOU
SURGICAL**

West Lake Trade Building Floor12
West Lake Avenue 1
Shangcheng District
Hangzhou, China

+86 / 571 / 87802003 (T)
+86 / 571 / 87177609 (F)
sales@shuyou-china.com

www.shuyou-china.com