



Medline Industries Ltd
 3rd Floor, Quayside Wilderspool Business Park
 Greenalls Avenue
 Warrington WA4 6HL
 United Kingdom

Tel.: +44 844 334 5237
 Fax: +44 844 334 5238
www.medline.eu/uk



Nipple and Nut Adaptor

Ref: 2555CE

Tubing Adaptor

Ref: 1423CE

Technical Sheet

Oxygen Adaptors Non-Sterile

Reference(s)

Item	Description	Use	Symbol	Units/case
1423CE	Tubing Adaptor, 22mm ID x 5-7mm tapered nipple, can be used with tracheostomy mask (ref 1075)	Single-use only, do not reuse		50
2555CE	Nipple and Nut hose adaptor for standard DISS oxygen outlets, swivel connector	Reusable up to 30 uses		50

Description

Connect standard oxygen tubing to a medical gas source and/or the patient.

Main materials

1423CE: Polypropylene
 2555CE: Polypropylene

Pack factor

Refer to the table

Packaging

Polyethylene bag (PE)

Legal Manufacturer

Medline Industries LP
 Three Lakes Drive, Northfield, IL 60093, USA
 EN ISO 13485

Regulation

In compliance with Medical Devices Regulation (EU) 2017/745 and UK MDR
 MD class: IIa non-sterile

•Please refer to complete labeling before using the product.

Other Information

NA

Revision: 9
 Document Title: FINT.441 - Technical sheet
 Document Number: EU1-00004-F-00025

© 2023 Medline Industries LP or one of its subsidiaries. All rights reserved
 TDS_Oxygen_Tubing_Reduction_Adaptors_Non-Sterile_EN01

Retention period: End of life and support period

