


PRODUCT SHEET



Products Name:	Micro Ear Forceps 8mm, 0.5 x 4 mm, straight, serrated
Glorious Art nr.	49.0940.00
Net weight:	24 gram
Material:	Stainless steel AISI 410
Sterility:	Non-Sterile (Re-Usable)
Exposure Temperature:	134°C (274°F)
Exposure Time:	10 Minutes
CE Classification:	Class-I
Shelf life:	3 years
Storage:	16°C to 30°C (Store the instruments in a closed cabinet, protected from dust, moisture and temperature fluctuations)
Life-span:	20 times of use: after each use, product needs to sterilize for ready to use in the next operation or next surgery.
Disposal:	No special decomposition or disposal required for Surgical Instruments as they don't contain any toxic or hazardous material

Glorious Danmark ApS