
ETHILON™

(Polymide 6) Suture

SUTURE

Description

Synthetic, Non-absorbable, Monofilament

Colour

Black, green or clear

Range

11-0 to 2 (USP). Supplied as needled sutures armed with, for example, MULTIPASS* and PRIME* needles. ETHILON* is also available as ligatures.

Indications












ETHILON* sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurosurgical procedures.

Please refer always to the Instructions for use/package insert that follows with the device for the most current and complete instructions.












ETHILON* Monofilament Polyamide 6

3/8 Circle Reverse Cutting Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 M-2	13 mm	6-0	BLACK	45	250 microns	36	667H
		5-0	BLACK	45	250 microns	36	668H
 C-3	13 mm	6-0	BLACK	75	350 microns	36	F2403
		5-0	BLACK	75	350 microns	36	F2412H
		4-0	BLACK	75	450 microns	36	F3206BH
 C-2	13 mm	5-0	BLACK	45	1/4 CIRCLE, 360 microns	36	644H
 FS-3	16 mm	6-0	BLACK	45		36	660H
		6-0	BLACK	45	VISI-BLACK*	36	EH7177H
		5-0	BLACK	45		36	1628H
		4-0	GREEN	75		36	G2435H
		4-0	BLACK	45		36	1162H
		3-0	BLACK	75		36	F3211BH
 FS-2	19 mm	5-0	BLACK	45		36	661H
		4-0	GREEN	45		12	G662G
		4-0	GREEN	75		36	G7143H
		4-0	BLACK	75		36	EH7143BH
		3-0	BLACK	45		36	653H
		3-0	BLACK	75		36	EH7797BH
		2-0	BLACK	90		36	F2448H
 FS-2S	19 mm	4-0	BLACK	45	Slimmed	36	662SLH
 FS-1	24 mm	4-0	BLACK	45		36	1629H
		4-0	BLACK	75		36	FJ5272BH
		3-0	BLACK	45		36	663H
		3-0	BLACK	75		36	669H
		3-0	BLACK	90		36	F2417H
		2-0	BLACK	45		36	EH7826BH
		2-0	BLACK	75		36	EH7827BH
		2-0	BLACK	90		36	F3223BH
 FS-1	26 mm	4-0	BLACK	45		12	W345BG
		3-0	BLACK	45		12	W320BG
		2-0	BLACK	45		36	664H
 FSL	30 mm	4-0	BLACK	75		36	1670BH
		3-0	BLACK	45		36	1171H
		3-0	BLACK	75		36	1671H
		2-0	BLACK	75		36	EH7798BH
 FSLX	36 mm	3-0	BLACK	75		36	1673BH
		3-0	BLACK	90		36	F2421H
		2-0	BLACK	75		36	1674BH
		0	BLACK	75		36	P7085H
 FSLX	45 mm	2-0	BLACK	100		12	W736G

ETHILON* Monofilament Polyamide 6

3/8 Circle Reverse Cutting PRIME* Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 P-6	8 mm	7-0	BLACK	45		36 1647H
 P-1	11 mm	7-0	BLACK	45	MULTIPASS* ETHALLOY*	36 1696G
		6-0	BLACK	45		36 697H
		5-0	BLACK	45		12 W1611G
 P-3	13 mm	6-0	BLACK	45	MULTIPASS* ETHALLOY*	36 1698H
		5-0	BLACK	45		36 698H
		4-0	UNDYED	45		12 691G
		4-0	BLACK	45		36 699H
 PS-3	16 mm	6-0	BLACK	45	MULTIPASS* ETHALLOY*	12 1665G
		5-0	BLACK	45		36 1668H
		5-0	GREEN	45		36 G1668H
		4-0	GREEN	45		36 GEH7761H
		4-0	BLACK	45		36 EH7761H
 PS-2	19 mm	5-0	BLACK	45	MULTIPASS* ETHALLOY*	36 1666H
		5-0	GREEN	45		12 G666G
		4-0	GREEN	45		12 G667G
		4-0	UNDYED	45		12 1611G
		3-0	BLACK	45		36 1669H
 PS-2 S	19 mm	4-0	BLACK	45	MULTIPASS* ETHALLOY*	36 1667SLH
		4-0	GREEN	45		36 G1667SLH
 PS-1	24 mm	4-0	BLACK	45	MULTIPASS* ETHALLOY*	36 EH7957BH
		3-0	BLACK	45		36 1663H
		3-0	BLACK	75		36 EH7563H
		2-0	BLACK	45		36 EH7553H
 PS	26 mm	3-0	BLACK	45	MULTIPASS* ETHALLOY*	36 W1625BH
		3-0	BLACK	75		36 W1685BH
		2-0	BLACK	75	ETHALLOY*	24 W1626BH
 PSL	30 mm	3-0	BLACK	75	ETHALLOY*	36 1691H

ETHILON* Monofilament Polyamide 6

3/8 Circle Conventional Cutting PRIME* Needles




Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PC-1	13 mm	6-0	BLACK	45	MULTIPASS* ETHALLOY*	12 1856G
 PC-3	16 mm	6-0	BLACK	45	MULTIPASS* ETHALLOY*	12 1866G
		5-0	BLACK	45		36 1865H
 PC-5	19 mm	4-0	BLACK	45	MULTIPASS* ETHALLOY*	12 1894G

ETHILON* Monofilament Polyamide 6


3/8 Circle Conventional Cutting PRIME* Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 PC-1	13 mm	6-0	BLACK	45	MULTIPASS* ETHALLOY*	12	1856G
 PC-3	16 mm	6-0	BLACK	45	MULTIPASS* ETHALLOY*	12	1866G
		5-0	BLACK	45		36	1865H
 PC-5	19 mm	4-0	BLACK	45	MULTIPASS* ETHALLOY*	12	1894G
 PC-25	26 mm	5-0	BLACK	45		12	W1858BG
 LR	90 mm	1	BLACK	2x50	With retention tubing	12	W2797

1/2 Circle Reverse Cutting Needles



Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 M-1	13 mm	5-0	BLACK	45		36	648H
 PS-4	16 mm	4-0	UNDYED	45	PRIME* MULTIPASS* ETHALLOY*	12	1603G
 J-1	17 mm	5-0	BLACK	75		36	F2423H
 X-1	22 mm	3-0	BLACK	45		12	642G
 CPX	48 mm	1	BLACK	100		12	W738G

Straight Cutting Needles













Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 KS	60 mm	3-0	BLACK	75		36	627H
		2-0	BLACK	75		36	628H

ETHILON* Monofilament Polyamide 6



1/2 Circle Taper Point Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CT 40 mm	0	BLACK	150	Loop	12	W740G
 CTX 48 mm	1	BLACK	150	Loop	12	W748




3/8 Circle Taper Point Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 BV130-3 3,7 mm, 150 microns	10-0	BLACK	13	ETHALLOY*	12	2820G
	9-0	BLACK	13	ETHALLOY*	12	2819G
 BV75-3 3,8 mm, 75 microns	10-0	BLACK	13		12	W2870
	9-0	BLACK	13		12	W2871
 BV50-3 3,8 mm, 50	11-0	BLACK	13		12	2881G
 BV130-4 5.0 mm, 150 microns	9-0	BLACK	13		12	W2813
	8-0	BLACK	13		12	W2812
 BVH100-3 5 mm, 100 microns	10-0	BLACK	13	1/2 circle	12	W2860
 BV75-4 5,1 mm, 75 microns	10-0	BLACK	13		12	W2850
 BV100-4 5,1 mm, 100 microns	10-0	BLACK	13		12	W2830
	10-0	BLACK	15		12	2814G
	9-0	BLACK	13		12	W2829
	9-0	BLACK	15		12	EH7448G
 BV-2 6,5 mm, 150 microns	10-0	BLACK	15		12	2810G
	9-0	BLACK	15		12	2809G
 BV130-5 6,5 mm, 150 microns	10-0	BLACK	13		12	W2810
	8-0	BLACK	15	MULTIPASS*	12	2808G
	8-0	BLACK	13		12	W2808
 BV130-5 6,5 mm, 150 microns	10-0	BLACK	30		12	W2801
 BV-1 9,3 mm, 250 microns	8-0	BLACK	45		12	W8170
 BV-1 9,3 mm, 250 microns	7-0	BLACK	45		36	EH7446H


3/8 Circle CS ULTIMA* Spatula

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-140-6 	6,5 mm, 150 microns	10-0	BLACK	15		12	7060G
		10-0	BLACK	30		12	U7061
CS-140-6 	6,5 mm, 150 microns	10-0	BLACK	30		12	U7003


3/8 Circle MICROPOINT* Spatula

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-35-C 	6 mm, 150 microns	10-0	BLACK	15		12	EH7980G
GS-9 GS-9 	6,5 mm, 200 microns	10-0	BLACK	30	Non-stocked	12	1718G
		9-0	BLACK	30		12	1717G
TG140-6  TG140-6	6,5 mm, 150 microns	10-0	BLACK	30	Reverse spatula	12	W1770
		9-0	BLACK	30		12	W1769


1/2 Circle MICROPOINT* Reverse Cutting Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
G-7 	8 mm	8-0	BLACK	45		12	W1765


1/2 Circle MICROPOINT* Spatula

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
TG160-6 	5,5 mm, 150 microns	10-0	BLACK	30		12	W1756


1/2 Circle CS ULTIMA* Spatula

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-160-6 	5,5 mm, 150 microns	10-0	BLACK	30		12	U7000



1/4 Circle MICROPOINT* Spatula

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
TG100-8 	6,5 mm, 220 microns	8-0	BLACK	30		12	B7714

1/4 Circle Spatula Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
S-14 S-14 	8 mm	5-0	BLACK	30	200 microns	12	B7712

1/4 Circle TAPERCUT* Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
V100-3 	3,5 mm, 100 microns	10-0	BLACK	13		12	W2899
		9-0	BLACK	13		12	W2898
V150-3 	4,5 mm, 150 microns	10-0	BLACK	15			W1701G