


YOUR CLASSIC CVC RANGE




REDUCE
PATIENT **RISKS**
WITH ALTIUS



A CLASSIC CATHETER WITH AN INSERTION KIT

Altius® Classic is manufactured using our proprietary Multi-Tube™ technology, which creates perfectly round lumens.

Perfectly round lumens deliver laminar flow, which allows for higher flows than traditionally extruded catheters.




THE RANGE

2 to 5 lumens in a variety of lengths.
See the Ordering Guide for offering & flow rates.

INSERTION KIT

Altius Classic insertion kit has been designed to reduce the risks during insertion.



Echogenic Needle:
Easy visualisation under ultrasound to prevent posterior wall puncture.




5ml Guidewire Introducer Luer Slip Syringe:
Facilitates guidewire advancement.




Altius Guidewire Dispenser: Ergonomic design to facilitate guidewire advancement even with wet gloves.




Kink-Resistant Nitinol Guidewire: To prevent guidewire kinking even in difficult anatomy.



No. 11 Mini Scalpel:
To make a small incision prior to dilation.



Hydrophilic Dilator:
For smooth dilation to prevent vessel trauma.



Secondary Fixation Wing

For more information, pricing or samples, speak to your local Kimal Representative, or contact:

UK

Tel: +44 (0)845 437 9542
Email: customercare@kimal.com

INTERNATIONAL

Tel: +44 (0)845 437 9543
Email: internationalenquiries@kimal.com

MIDDLE EAST, ASIA & AFRICA

Tel: +97 16 552 8146
Email: mea@kimal.com

© 2025 Kimal PLC. All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from Kimal PLC.

Kimal PLC reserve the right to change the specification of its products at any time.
Kimal PLC, Clayfield Rd, Worcester WR4 0AE.

BRO-155-(3) | AUGUST 2025

KIMAL
delivering healthcare innovation