

PRODUCT SHEET



Products Name:	42.0394.14LQ SYGEPLEJESAKS, 14 CM, RUSTFRIT STÅL, MED KNOP, LQ
Producer name:	M.A. Arain & Brothers (PVT)LTD.
Finish:	Satin
Net weight:	32 gram (approximate)
Material:	Stainless steel AISI 410
Sterility:	Non-Sterile (Re-Usable)
Exposure Temperature:	134°C (274°F)
Exposure Time:	10 Minutes
CE Classification:	Class-I
Storage:	16°C to 30°C (Store the instruments in a closed cabinet, protected from dust, moisture and temperature fluctuations)
Shelf life:	5 years
Life-span:	20 times of use: after each use, product needs to sterilize for ready to use in the next operation or next surgery.
Disposal:	No special decomposition or disposal required for Surgical Instruments as they don't contain any toxic or hazardous material

Producent



M. A. ARAIN C BROTHERS (PVT)LTD.
Adresse: Wazirabad Troad, Sialkot – 51310 Pakistan
Kontakt: Tel. (+92) 52 4264571, 52 4264926
Email: info@maarain.com

EU Importør/Distributør

Glorious Danmark ApS Adresse: Naverland 11, DK2600 Glostrup
Email: Glorious@glorious-enterprises.dk Web: www.glorious-danmark.dk