

Atpro Whey 3 kg - English Translation

Temperature: MIN 18°C, MAX 25°C

Brand: Toft Care®

Product Name: Atpro Whey 3 kg

Net Content: 3 Kilograms

Functional Name: Whey Protein Powder

Regulated Name: Whey Protein Concentrate

GTIN: 05701477501377

Supplier No.: 3501370

Country of Origin: EU

Packaging: Bucket

Guaranteed Shelf Life from Delivery: 405 days

Total Shelf Life from Production: 540 days

Shelf Life After Opening: -

Ingredients: Whey protein concentrate (milk).

Directions for Use: The powder is mixed with the desired amount of liquid. Add ice cubes if desired to chill the drink and make it even more palatable. Suggestion: 5-10 g per 100 ml of liquid and stir.

Storage: Store in a dry place, not too warm, and in tightly sealed packaging.

Allergens: Cereals containing gluten, crustaceans, eggs, fish, peanuts, soy, nuts, celery, mustard, sesame seeds, sulphur dioxide and sulphites, lupin and molluscs: Free From. Milk: Contains.

Nutritional Information per 100 g

Energy: 1855 kJ / 442 kcal

Fat: 17 g

of which saturated fat: 10 g

Carbohydrates: 0.5 g

of which sugars: 0.5 g

Protein: 70.2 g

Salt: 0.23 g

Preparation Instructions: Measure 20 g of protein powder and add 100-150 ml of liquid. Stir well.

Contains: Milk

Company Information: Toft Care A/S (5790002259166)

Date of Last Change: 29-06-2026 13:38:27