

Phase in Announcement

Date : 16-07-2026

Brand : Klinion®

Product : Kliniderm® Super Fiber + Antimicrobial Polymer Technology (APT™)

Excited to announce the phase in of new items in the Klinion® Kliniderm® range!

PHASE IN ITEM		
Product Reference	Medeco article #	Article description
40511813	3074333	KLINION® Kliniderm® Super Fiber + APT™ 5x5cm
40511814	3074370	KLINION® Kliniderm® Super Fiber + APT™ 10x10cm
40511815	3074371	KLINION® Kliniderm® Super Fiber + APT™ 15x15cm
40511816	3074372	KLINION® Kliniderm® Super Fiber + APT™ 2x45cm

Product Description

Kliniderm Super Fiber + APT™ are 200 gsm fibers, they are soft and conformable non-woven dressings made from sodium, carboxymethylcellulose (CMC) fibers and strengthening cellulose fibers with antimicrobial polymer, quaternary polyethyleneimine (QPEI). The QPEI polymer when incorporated into the wound dressing, and in contact with wound exudate, prevents bacterial colonization in the dressing due to its antimicrobial activity. The absorption of exudate into the dressing forms a gel that assists in maintaining a moist wound environment, supports autolytic debridement and improves the wound edge and surrounding skin from maceration. Through the gel formation the structure of the dressing remains intact. When dry, Kliniderm Super Fiber + APT™ can easily be cut to the size of the wound. Even when the dressing is moist, as gelled fiber, its structure remains intact. The high vertical absorption of exudate into the fiber dressing improve the wound environment and the wound edge. The product can also be used under compression. Kliniderm Super Fiber + APT™ are provided in different sizes, in rope or pad, with a basic weight of 200gsm and composed of 80% w/w CMC and 20% Tencel/Lyocell fibers.

Phase in Announcement

Indications

The product is for single use only. Under the supervision of a healthcare professional, Kliniderm Super Fiber + APT™ wound dressings may be used for the management of:

- Wounds with moderate to heavy exudate
- Acute or chronic wounds
- Leg ulcers, pressure ulcers (stage ii to iv) and diabetic ulcers
- Surgical wounds (e.g. post-operative, wounds left to heal by secondary intent and donor sites)
- Traumatic wounds (e.g. abrasions and lacerations)
- Exuding oncology wounds
- Superficial or cavity wounds
- Improves the wound edge and the surrounding skin from maceration¹.

¹ The properties mentioned above refer to the current in vitro data.

Changes

Please find details on the new product in the tables below.

Product Reference	Medeco article #	Description	Pieces per peel pouch	Pieces per inner carton	Pieces per outer carton
40511813	3074333	KLINION® Kliniderm® Super Fiber + APT™ 5x5cm	1	10	200
40511814	3074370	KLINION® Kliniderm® Super Fiber + APT™ 10x10cm	1	10	100
40511815	3074371	KLINION® Kliniderm® Super Fiber + APT™ 15x15cm	1	5	50
40511816	3074372	KLINION® Kliniderm® Super Fiber + APT™ 2x45cm	1	5	100

Product Reference	Medeco article #	Description	EAN code peel pouch	EAN code inner carton	EAN code outer carton
40511813	3074333	KLINION® Kliniderm® Super Fiber + APT™ 5x5cm	7640328070458	7640328070373	7640328070410
40511814	3074370	KLINION® Kliniderm® Super Fiber + APT™ 10x10cm	7640328070342	7640328070380	7640328070427
40511815	3074371	KLINION® Kliniderm® Super Fiber + APT™ 15x15cm	7640328070359	7640328070397	7640328070434
40511816	3074372	KLINION® Kliniderm® Super Fiber + APT™ 2x45cm	7640328070366	7640328070403	7640328070441

Phase in Announcement

Product Reference	Medeco article #	Description	Dimension peel pouch	Dimension inner carton	Dimension outer carton
40511813	3074333	KLINION® Kliniderm® Super Fiber + APT™ 5x5cm	90x142 mm	97x40x147 mm	406x205x152 mm
40511814	3074370	KLINION® Kliniderm® Super Fiber + APT™ 10x10cm	144x184 mm	146x40x193 mm	298x199x206 mm
40511815	3074371	KLINION® Kliniderm® Super Fiber + APT™ 15x15cm	210x240 mm	213x26x245 mm	215x265x250 mm
40511816	3074372	KLINION® Kliniderm® Super Fiber + APT™ 2x45cm	95x260 mm	97x40x263 mm	297x269x280 mm

Product Reference	Medeco article #	Description	Gross weight peel pouch	Gross weight inner carton	Gross weight outer carton
40511813	3074333	KLINION® Kliniderm® Super Fiber + APT™ 5x5cm	2.43 grams	48.9 grams	1150 grams
40511814	3074370	KLINION® Kliniderm® Super Fiber + APT™ 10x10cm	6 grams	95 grams	1100 grams
40511815	3074371	KLINION® Kliniderm® Super Fiber + APT™ 15x15cm	12.1 grams	108.5 grams	1230 grams
40511816	3074372	KLINION® Kliniderm® Super Fiber + APT™ 2x45cm	5.5 grams	61.3 grams	1460 grams

Timeline

This item is available as from today. If you wish to purchase this item, please make sure to submit your official forecast through Slimstock. In case you do not use the Slimstock system, add it to your monthly forecast file which is shared with Medeco.

Phase in Announcement

Picture

