


Art. No.	1030128
NAME OF PRODUCT	Wooden applicator
REVISION NO	A/1
REPLACED REVISION	A/0
VALID FROM	From 21NS-S193
CHANGE DESCRIPTION	printing
DESCRIPTION OF PRODUCT	Wooden applicator single small head
	Size:150×2.2mm
	Remark: Specifications pls. refer to attachments I, II
INTENDED USE OF PRODUCT	To be a testing tool for examining inthe hospital and laboratory


RAW MATERIAL		
CODE FOR THE ENCLOSED	RAW MATERIAL	ADDITIONAL INFORMATION
RAW MATERIAL SPECIFICATION	Wooden applicator	Birch
	Cotton head	Nature cotton ,white
	Closure tape	PP transparent tape
	Inner zip bag	PE, transparent
	Inner box	Single Wall Corrugated Fibreboard, white Weight: 450g/m ²
	Transport carton	Double Wall Corrugated Fibreboard Weight: 600-700g/m ²

QUALITY SPECIFICATION OF THE PRODUCT			
CHARACTERISTIC	TARGET VALUE	QUALITY TOLERANCES	TEST METHOD
FINISHED SIZE	Stick=148(length)×2.2(Dia.)mm	+/- 1mm(length), +/- 0.1mm(Dia.)	Measuring
	Head=12(length)×4.5(Dia.)mm	+/-2mm(length,), +/- 0.5mm(Dia.)	Measuring
	Total length=150mm	149~152mm	Measuring

DESCRIPTION OF THE PACKAGING SYSTEM			
PACKING			
100pcs/bag, 20bags/box, 10boxes/ctn			
ZIP PE BAG SIZE (cm)			
L	W	T	
18	11	50mic	
BOX SIZE (cm)			
L	W	H	
21	18.5	12	
CARTON SIZE (cm)			
L	W	H	
62	23	38.5	

Wooden applicator

Attachment I



Length 1= 149~152mm

Length 2= 10~14mm

Dia. 1 = 2.1~2.3mm

Dia. 2 = 4~5mm


Printing/label

Attachment II

Bag printing



Box printing



Carton label printing

