

Cenaman® Multi-use ENFit™ syringes

Description

The Enteral Syringes are intended to be used for delivering liquid medication and / or formulated enteral nutrition / feed either orally or via a feeding tube into the gastrointestinal tract.

The Cenaman ENFit™ syringes reduce the risk of an unintended misconnection.

By allowing multiple uses on the same patient, these syringes reduce medical waste and contribute to a more sustainable healthcare practice.



Article overview

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3069013	Cenaman Multi-use Enfit Syringe 1ml 409330 Sterile	100 pcs.	3600 pcs.			409330
	3069014	Cenaman Multi-use Enfit Syringe 2,5ml 409340 Sterile	100 pcs.	3200 pcs.			409340
	3069015	Cenaman Multi-use Enfit Syringe 5ml 409350 Sterile	100 pcs.	2400 pcs.			409350
	3069016	Cenaman Multi-use Enfit Syringe 10ml 409360 Sterile	100 pcs.	1500 pcs.			409360
	3069017	Cenaman Multi-use Enfit Syringe 20ml 409370 Sterile	100 pcs.	900 pcs.			409370
	3069268	Cenaman Multi-use Enfit Syringe 20ml X30 409370-30 Sterile	30 pcs.	240 pcs.			409370-30
	3069018	Cenaman Multi-use Enfit Syringe 60ml Decentral 409380 Sterile	50 pcs.	200 pcs.			409380
	3069269	Cenaman Multi-use Enfit Syringe 60ml X30 409380-30 Sterile	30 pcs.	90 pcs.			409380-30