

ETHILON*

(Polyamide 6 or 6,6)
Suture



ETHICON

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

ETHILON*

Monofilament Polyamide 6 / 6,6

SUTURE

Description: **Monofilament** For smooth passage through tissue.
 Synthetic For minimal tissue reaction.
 Non-absorbable Provides prolonged tensile strength retention in tissue.

Color: Green, black or clear.

Range: 11-0 to 2 (USP). Supplied as needled sutures armed with, for example, MULTIPASS* and PRIME* needles. ETHILON* is also available as ligatures.

Indications: ETHILON* (black, green and clear) sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurosurgical procedures. Blue ETHILON* sutures are intended for use in skin closure.

Please refer always to the Instructions for use / package insert that follows with the device for the most current and complete instructions.

