

KLINION - Nitrile Sensitive Long Gloves

Description

The Klinion Protection examination gloves are intended to be worn on the hand for medical purposes and to help prevent contamination between patient and examiner. This information does not reflect the actual duration of protection in the workspace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400mm – where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. This product provides protection against bacteria, fungi and virus. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. Before usage inspect the gloves for any defect or imperfections. For single use only. Re-use can cause contamination or infection. Dispose according to local procedures and guidelines. Not made with natural rubber latex. Medical device.

Product information

In compliance with: Regulation (EU) 2017/745 Medical Device Regulation (MDR) Class I & PPE Regulation (EU) 2016/425 Category III










Klinion Nitrile gloves are available in:

- Ultra Comfort (Thin glove in a white or blue color)
- Ultra Save (Better protection than Ultra Comfort)
- Sensitive (Accelerator-free gloves to reduce the risk on allergies)

The total length of our standard gloves is minimal 240 mm and our long cuff gloves 290 mm.



Article overview

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3051843	KLINION PERSONAL PROTECTION EXAMINATION GLOVE SENSITIVE LONG CUFF NITRILE XS INDIGO 102488	150 pcs.	1500 pcs.			102488
	3050618	KLINION PERSONAL PROTECTION EXAMINATION GLOVE SENSITIVE LONG CUFF NITRILE S INDIGO 102484	150 pcs.	1500 pcs.			102484
	3050655	KLINION PERSONAL PROTECTION EXAMINATION GLOVE SENSITIVE LONG CUFF NITRILE M INDIGO 102485	150 pcs.	1500 pcs.			102485
	3050656	KLINION PERSONAL PROTECTION EXAMINATION GLOVE SENSITIVE LONG CUFF NITRILE L INDIGO 102486	150 pcs.	1500 pcs.			102486
	3050657	KLINION PERSONAL PROTECTION EXAMINATION GLOVE SENSITIVE LONG CUFF NITRILE XL INDIGO 102487	150 pcs.	1500 pcs.			102487
	3050697	KLINION PERSONAL PROTECTION EXAMINATION GLOVE SENSITIVE LONG CUFF NITRILE XS BLUE 102534	150 pcs.	1500 pcs.			102534
	3050698	KLINION PERSONAL PROTECTION EXAMINATION GLOVE SENSITIVE LONG CUFF NITRILE S BLUE 102535	150 pcs.	1500 pcs.			102535
	3050699	KLINION PERSONAL PROTECTION EXAMINATION GLOVE SENSITIVE LONG CUFF NITRILE M BLUE 102536	150 pcs.	1500 pcs.			102536
	3050700	KLINION PERSONAL PROTECTION EXAMINATION GLOVE SENSITIVE LONG CUFF NITRILE L BLUE 102537	150 pcs.	1500 pcs.			102537
	3050701	KLINION PERSONAL PROTECTION EXAMINATION GLOVE SENSITIVE LONG CUFF NITRILE XL BLUE 102538	150 pcs.	1500 pcs.			102538